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Quick Setup Guide for InSUS Head for Optical Single-Use O₂ Sensors

1 Introduction

Thank you for buying this optical oxygen sensor from METTLER TOLEDO. The InSUS™ Head is intended solely for in-line measurement of the oxygen partial pressure in combination with METTLER TOLEDO optical single-use O₂ sensors, and the InSUS fiber optical cable with integrated Pt1000 temperature probe. The O₂ single-use sensors are typically integrated in single-use devices such as single-use bags and other similar applications. Please also refer to the documentation provided by your single-use device supplier.

2 Safety instructions

 **Note:** Before every start-up, the sensor head must be checked for:

- Damage to the connections, fastenings, etc.
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

 **Note:** The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

 **Caution:** Before the sensor head is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

 **Caution:** A defective sensor head may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

 **Caution:** Some components inside the sensor head are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero potential before any operation on the wiring terminal.


 **Caution:** If a sensor head is broken it must be replaced immediately.

 **Caution:** It is essential that personnel involved in the operation or maintenance of the sensor head and any of the associated equipment (e.g. single-use bags, InSUS fiber optical cable, transmitters, etc.) is properly trained in the use and handling of the equipment.

 **Warning:** Before dismantling the sensor head or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (e.g. power supplies switched off).

3 Product description

The sensor heads are delivered fully assembled, factory tested and calibrated for correct function. A certificate of calibration is provided together with the sensor head:

- 1: 5 Pin connector
- 2: InSUS head
- 3: InSUS head connector
- 4: InSUS fiber optical cable
- 5: Sensor connection (M16 thread)
- 6: Single-use sensor (integrated in single-use device)

The sensor head is connected to a single-use O₂ sensor with the InSUS fiber optical cable.

 **Note:** Cable assignment can be found in the METTLER TOLEDO cable instruction manual.

3.1 Digital connection of the InSUS Head to a transmitter

Transmitter M400 or M800

The sensor head is connected to the transmitter via 5 pin data cable which is available in different lengths.

Wiring information

		M400	M800 1-channel	M800 2/4-channel
Color	Function	TB3	TB3	TB2
brown	24 V	18	7	9
black	GND24V	17	8	10
gray	shield	13	2	12
yellow	shield	13	5	15
blue	RS485B	14	3	13
white	RS485A	15	4	14

3.2 Analog connection of the InSUS Head to a transmitter or biocontroller

A METTLER TOLEDO M50 signal converter RS 485 / nano-Ampere mimics a polarographic Clark sensor. The M50 converter is connected to an analog signal input of a transmitter, e.g. M300, or biocontroller.

4 Operation

4.1 Calibration

Calibration must be performed every time a new single-use sensor is connected to the InSUS Head via the InSUS optical fiber cable. For most applications, a single point calibration should be sufficient. The calibration medium can be either air or a calibration gas with known O₂ concentration or liquid medium with a known O₂ concentration. Before starting the calibration in gas, the correct pressure and the correct humidity have to be set in the transmitter. For detailed information please refer to the transmitter or biocontroller manual.

5 Maintenance

5.1 Visual inspection

To check your sensor head, we recommend the following procedure:

- The contacts of the connector must be dry: moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Check the contacts of the InSUS fiber optical cable.

6 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

7 Warranty

On manufacturing defects, 12 months after delivery.

8 Storage conditions

Please inspect the instrument immediately after receiving for any transport damage. The sensor Head should be stored in a dry, clean and protected area until the time of installation. If the sensor Head has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

9 EC Declaration of Conformity

The complete declaration is available in the declaration of conformity certificate.

InSUS is a trademark of the METTLER TOLEDO Group.

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 本表依据SJ/T 11364的规定编制，本产品符合以下标志规范：
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 03/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 30 313 265 C

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 3 0 3 1 3 2 6 5 C *