Quick Setup Guide InPro 6000/InTap Optical O₂ Sensors

cs	Průvodce rychlým nastavením	2
da	Lynvejledning	8
de	Kurzanleitung	14
en	Quick Setup Guide	20
es	Guía de configuración rápida	26
fi	Pika-asetusopas	32
fr	Guide de paramétrage rapide	38
hu	Gyorsbeállítási útmutató	44
it	Guida alla configurazione rapida	50
ja	クイック セットアップ ガイド	56
ko	사용 매뉴얼	62
nl	Beknopte handleiding	68
pl	Podręcznik szybkiej konfiguracji	74
pt	Guia de Configuração Rápida	80
ru	Руководство по быстрой настройке	86
sv	Snabbguide	92
th	คู่มือการตั้งค่าอย่างรวดเร็ว	98
zh	快速设置指南	104

InPro 6000/InTap Optické senzory O₂

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Obsah

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Quick Setup Guide

Contents

1	Introduction	21
2	Safety instructions	21
3	Product description	22
3.1	Analog integration	22
4	Operation	23
5	Maintenance	23
5.1	Replacing the OptoCap	24
6	Environmental protection.....	25
7	Warranty	25
8	Storage conditions	25
9	EC Declaration of Conformity	25

Appendix

Technical Drawings	112
--------------------------	-----

InPro, InTap and ISM are trademarks of METTLER TOLEDO Group.

1 Introduction

Thank you for buying this optical oxygen sensor from METTLER TOLEDO. The InPro™ 6000 optical O₂ sensors are intended solely for in-line measurement of the oxygen partial pressure. The optical oxygen sensors with integrated temperature probes are used for measurement of oxygen in hygienic applications such as bio-tech fermentation and food & beverage applications.

The sensors are sterilizable and compatible with CIP (cleaning-in-place). For fermentation applications the InPro 6860i is applicable and fully sterilizable/autoclavable. The InPro 6960i and 6970i are intended for use in the food & beverage industry for measuring high and low ppb levels of oxygen respectively. The InTap™ sensor is intended to be used in the portable oDO analyzer (30 425 550/30 457 912).

2 Safety instructions

Notice: Before every start-up, the sensor must be checked for:

- Damage to the connections, fastenings, etc.
- Leakage
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

Notice: The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

Caution: Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

Caution: A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

Caution: The sensor is not equipped with heat protection. During steam-sterilization procedures in pipes, the surface of the sensor can reach high temperatures and cause burns.

Caution: Some components inside the sensor are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero potential before any operation on the wiring terminal.

Warning: Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (pressurize, no explosion risk, empty, rinse, vent, etc.).

3 Product description

The sensors are delivered fully assembled, factory tested and calibrated for correct function. A quality certificate and an inspection certificate 3.1 acc. EN10204.3/1B is provided.

Drawings of the sensors see Appendix A (InTap sensor not shown).

The sensors are installed by using a VP-Cable (InPro 6860i) or a 5 pin cable (InPro 6970i/InPro 6960i/InTap sensor). Cables in different lengths 2–50 m are available from METTLER TOLEDO.

Notice: Cable assignment can be found in the METTLER TOLEDO cable installation manual.

3.1 Analog integration

Integration without transmitter (InPro 6860i nA/InPro 6860i mA):

For analog integration two versions of the InPro 6860i exist.

- The simulated nA signal is used for legacy biocontroller integration
- The 4/20 mA-output incl. HART signal is used to integrate the sensor directly into a process logistic controller (PLC).

RS485 cable for InPro 6870i/InPro 6960i/InPro 6870i

Color	Function	M400	M800	M800	M400
		TB4	1/2 channel TB2	4 channel TB2 or TB4	2w/FF/PA TB2
brown	24 VDC+	1	9	9	n.c.
black	24 VDC–	2	10	10	n.c.
gray	shield	6	12	12	M
yellow	shield	6	15	15	n.c.
blue	RS485–	7	13	13	N
white	RS485+	8	14	14	O

VP8 cable InPro 6860i/nA

Color	Function	M400	M800	M800	M400
		TB4	1/2 channel TB2	4 channel TB2 or TB4	2w/FF/PA TB2
gray	24 VDC+	1	9	9	n.c.
blue	24 VDC–	2	10	10	n.c.
pink	RS485+	8	14	14	O
brown	RS485–	7	13	13	N
green/ yellow	shield	4	12	12	M

For M400 2-w, M400FF and M400PA, the optical oxygen sensor needs to be powered separately with an appropriate power supply.

Warning: Never connect an InPro 6860i with mA output to a nA input. This could damage the electronics of the sensor and/or the instrument the sensor is connected to.

The RS-485 signal also includes a MODBUS RTU signal, for a MODBUS installation and additional material including configuration tools can be downloaded at **www.mt.com/pro-ODO-documentation**

4 Operation

The sensor is supplied ready for use. Before using the sensor, remove the protecting cap. For continuous applications we recommend periodic recalibration in line with your requirements on accuracy.

In case of any failure or malfunction, the equipment in which the sensor is installed, needs to be in a safe condition (depressurized, empty, rinse, vent, etc.).

5 Maintenance

The incorporated Intelligent Sensor Management™ (ISM) monitors and displays the remaining lifetime of the sensing element as the Dynamic Lifetime Indicator.


Depending on the aggressiveness (CIP/SIP) of the application conditions the sensor's sensing element, the OptoCap experiences a degree of stress. Typically, after 6 to 12 months this part needs to be exchanged followed by a zero point and air calibration.

This calibration can be performed by using the transmitter's calibration menu and/or using METTLER TOLEDO's calibration and maintenance software tool, iSense.


www.mt.com/iSense

www.mt.com/pro-transmitters

5.1 Replacing the OptoCap

Unscrew and remove the spent OptoCap and the old O-ring and replace with new ones. After this procedure the sensor needs to be calibrated with an air and zero-point calibration using the correct humidity and ambient pressure settings.

Type of calibration
Air and Zero-Point Calibration

Currently used salinity 0.0 g/kg

Humidity in % 40.0 %

Pressure in mbar 980 mbar

Use USB

6 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

7 Warranty

On manufacturing defects, 12 months after delivery.

8 Storage conditions

Inspect the sensor immediately after receiving for any transport damage. The sensor should be stored in a dry, clean and protected area until the time of installation. If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

9 EC Declaration of Conformity

The complete declaration is available in the declaration of conformity certificate.

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
本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

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
For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 10/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 30 313 394 C

Management System
certified according to
ISO 9001 / ISO 14001

**UK
CA**₂₅₀₃

CE₁₂₅₈

EAC

www.mt.com/pro

* 3 0 3 1 3 3 9 4 C *