## Quick Setup Guide InPro 6000 Amperometric O<sub>2</sub> Sensors

cs	Průvodce rychlým nastavením	2
da	Lynvejledning	8
de	Kurzanleitung	14
en	Quick Setup Guide	20
es	Guía de configuración rápida	26
fi	Pika-asetusopas	32
fr	Guide de paramétrage rapide	38
hu	Gyorsbeállítási útmutató	44
it	Guida alla configurazione rapida	50
ja	クイック セットアップ ガイド	56
ko	사용 매뉴얼	62
nl	Beknopte handleiding	68
pl	Podręcznik szybkiej konfiguracji	74
pt	Guia de Configuração Rápida	80
ru	Руководство по быстрой настройке	86
sv	Snabbguide	92
th	คู่มือการตั้งค่าอย่างรวดเร็ว	98
zh	快速设置指南	104

# InPro 6000

## Ampérometrické senzory O<sub>2</sub>

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---

#### Obsah

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---

### Contents

1	Introduction .....	21
2	Safety instructions .....	21
3	Product description .....	21
4	Installation .....	22
4.1	Connection .....	22
4.2	Mounting the sensor .....	22
5	Operation .....	23
5.1	Start-up and polarizing .....	23
5.2	Calibration .....	23
6	Maintenance .....	23
6.1	Refilling electrolyte; exchanging membrane body .....	24
7	Disposal .....	25
8	Warranty .....	25
9	Storage conditions .....	25
10	EC Declaration of Conformity .....	25

### Appendix

Technical Drawings .....	112
--------------------------	-----

InPro and ISM are trademarks of the METTLER TOLEDO Group.

## 1 Introduction

Thank you for buying this amperometric oxygen sensor from METTLER TOLEDO. The InPro™ 6000(i) series amperometric O<sub>2</sub> sensors are intended solely for in-line measurement of the oxygen partial pressure e.g., oxygen measurement in liquids.

Application examples are (but not limited to), fermentation in biotechnology, oxygen measurement in brewing and chemical applications.

The following sensor types are covered within this quick set-up guide.

InPro 6800/InPro 6850 i/InPro 6900/InPro 6900 i/  
InPro 6950/InPro 6950 i/easySense.

Sensors without 'i' are analog sensors and with 'i' are digital sensors incorporating Intelligent Sensor Management (ISM™). EasySense is digital sensor without ISM. ISM functionality offer Plug and Measure as well as enhanced diagnostics features.

## 2 Safety instructions

**Notice:** Before every start-up, the sensor must be checked for:

- Damage to the connections, fastenings, etc.
- Leakage
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

**Notice:** The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

**Caution:** Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

**Caution:** A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

**Caution:** The sensor is not equipped with heat protection. During steam-sterilization procedures in pipes, the surface of the sensor can reach high temperatures and cause burns.

**Caution:** Some components inside the sensor are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero potential before any operation on the wiring terminal.

**Warning:** Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (pressurize, no explosion risk, empty, rinse, vent, etc.).

## 3 Product description

**Notice:** Digital sensors (ISM): All sensors ending with an 'i' are delivered dry and must be filled with electrolyte before start-up. **Electrolyte has to be ordered separately.**

Analog sensors are delivered filled with electrolyte, however, it is advised to replace the electrolyte reservoir with fresh solution.

**Article Nr.: Electrolytes**

59 907 065 – O<sub>2</sub> Electrolyte 25 mL (InPro 6050; InPro 6800; InPro 6850 i; easySense)

52 206 107 – O<sub>2</sub> Electrolyte 25 mL High Salinity [InPro 6900 (i); InPro 6950 (i)]

30 298 424 – O<sub>2</sub> electrolyte pack 3 × 25 mL (InPro 6050; InPro 6800; InPro 6850 i; easySense)

30 298 425 – O<sub>2</sub> electrolyte pack – high salinity 3 × 25 mL [InPro 6900 (i); InPro 6950 (i)]

**Warning:** The O<sub>2</sub> electrolyte has an alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.** If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.



## 4 Installation

### 4.1 Connection

The sensors are installed to a transmitter by using a VP (analog sensors) or AK9 (ISM sensors) connector cable. Cables in different lengths 2–50 m are available from METTLER TOLEDO.

**Notice:** Cable and terminal assignment can be found in the METTLER TOLEDO cable and/or transmitter instruction manual.

### 4.2 Mounting the sensor

Before mounting the sensor into a housing remove the protection cap on the sensor. The sensor can then be mounted into an InFit/InTrac Series housing from METTLER TOLEDO. Alternatively, the 12 mm sensors can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.

**Notice:** For connecting the cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual ([www.mt.com/pro-transmitters](http://www.mt.com/pro-transmitters))

## 5 Operation

### 5.1 Start-up and polarizing

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O<sub>2</sub> sensor master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O<sub>2</sub> transmitter or to a sensor master. After 6 hours, the sensor will be fully polarized and ready for operation.

For polarizing and storage the cleaning and conditioning solution (order Nr. 52 200 255) is recommended. This solution contains an oxygen scavenger which ensures that the electrolyte of the sensor does not come into contact with oxygen when not in use.

**Attention:** The InPro 6950 i should never be polarized at air!

It is recommended to store the sensors in the according electrolyte or cleaning and conditioning solution during the polarization time. The InPro 6900 i and 6950 i use a special electrolyte which contains an oxygen scavenger.

**Warning: Setting of the polarization voltage on the transmitter for correct measurements:**

Standard applications: –675 mV.

Measurements of continuously lowest oxygen concentrations (<500 ppb in liquids or <10,000 ppm [vol.] in gases) and in presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries): –500 mV.

### 5.2 Calibration

With digital ISM sensors all calibration data are stored in the sensor. Once filled with electrolyte and polarized they can be connected to an ISM transmitter and immediately used for measurement. Analog sensors are delivered wet and after polarization can be slope 1-pt, 2pt- or process calibrated using a transmitter.

With ISM sensors METTLER TOLEDO's calibration and maintenance software tools iSense/iSense light are recommended for sensor checking and calibration. ([www.mt.com/iSense](http://www.mt.com/iSense)).

Polarization voltage for calibration purposes are always put to –675 mV.

For measurements in carbon dioxide containing solutions the transmitter setting has to be switched to –500 mV accordingly after the calibration.

## 6 Maintenance


If the membrane and/or the interior body exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body and/or the interior body has to be replaced.

The incorporated Intelligent Sensor Management (ISM) technology monitors and reports the remaining lifetime of the inner bodies as the Dynamic Lifetime Indicator.

Depending on the aggressiveness (CIP/SIP) of the application conditions a replacement of electrolyte, the membrane body or the interior body is recommended. This maintenance can be performed by using the transmitter's calibration menu and/or using METTLER TOLEDO's calibration and maintenance tool iSense ([www.mt.com/iSense](http://www.mt.com/iSense)).

For detailed maintenance description and criteria please refer to iSense or the main sensor manual.

## 6.1 Refilling electrolyte; exchanging membrane body

**Procedure**

Cleanse with demin water; turn cap sleeve twice, do not unmount yet. From the bottom push the membrane body out by a slight upwards push; unscrew completely. Remove cap sleeve and membrane body. Clean inner body. If necessary use new membrane body and fill with fresh electrolyte solution. Make sure no air bubbles are trapped. Replace O-ring each time when replacing the membrane body and screw the cap sleeve on carefully.

**7 Environmental protection**

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

**8 Warranty**

On manufacturing defects, 12 months after delivery.

**9 Storage conditions**

Inspect the sensor immediately after receiving for any transport damage. The sensor should be stored in a dry, clean and protected area until the time of installation. If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

**10 EC Declaration of Conformity**

The complete declaration is available in the declaration of conformity certificate.

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 本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：  
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下  
 ○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。  
 x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

产品中有害物质的名称及含量 Toxic and hazardous substance name and containment in product						
部件名称 Part Name	有害物质或元素 Toxic and hazardous substances					
	铅 Lead (Pb)	汞 Mercury (Hg)	镉 Cadmium (Cd)	六价铬 Hexavalent chromium (Cr6+)	多溴联苯 Polybrominated biphenyls (PBB)	多溴二苯醚 Polybrominated diphenyl ethers (PBDE)
插头 Elect. Connector	x	○	○	○	○	○
电路板 PCBA	x	○	○	○	○	○


本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：  
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下  
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。  
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

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For addresses of METTLER TOLEDO  
Market Organizations please go to:  
**[www.mt.com/contacts](http://www.mt.com/contacts)**

**METTLER TOLEDO Group**  
Process Analytics  
Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

Subject to technical changes  
© 10/2022 METTLER TOLEDO  
All rights reserved  
Printed in Switzerland. 30 313 395 C

Management System  
certified according to  
ISO 9001 / ISO 14001

**UK  
CA**<sub>2503</sub>

**CE**<sub>1258</sub>

**EAC**

**[www.mt.com/pro](http://www.mt.com/pro)**

\* 3 0 3 1 3 3 9 5 C \*