METTLER TOLEDO

InPro 3200 (SG) pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction manual	14
es	Manual de instrucciones	18
fi	Käyttöohje	22
fr	Instructions d'utilisation	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
jp	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

Document Color Code:

magenta = check translation

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
Instruction Manual

Contents

- 1 Introduction 15
- 2 Safety instructions 15
- 3 Product description 15
- 4 Installation and commissioning 16
- 5 Operation 16
 - 5.1 Calibration of the sensor and pH transmitter 16
 - 5.2 Sterilization of sensor 16
 - 5.3 Temperature sensor 17
- 6 Maintenance 17
- 7 Trouble-shooting 17
- 8 Disposal 17
- 9 Guarantee 17

InPro is a registered trademark of **METTLER TOLEDO Group** in Switzerland, the USA, the European Union and a further five countries.

1 Introduction

METTLER TOLEDO pH sensors type InPro™ 3200 are low-maintenance electrodes with pressurized gel-electrolyte, silver-ion trap and built-in temperature sensor.

The InPro 3200 sensors are suitable for the same applications as the corresponding DPA(S) sensors without RTD, and provide an additional temperature signal for automatic temperature compensation of the pH signal during calibration and operation.

InPro 3200 sensors are autoclavable and in-situ sterilizable.

2 Safety instructions

The sensor consists partly of glass and should be handled carefully. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. Local safety regulations are applicable to the handling of the sensor. Please follow the guidelines enclosed for an installation in the Ex zone. Please take extra care during installation/fitting of the sensor to ensure that neither the polymer body/shaft of the sensor nor the thread (Pg 13.5) suffer any injury or damage.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

3 Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro™ 3200/SG/xxx/yyy	Type designation/SG = ancillary electrode/xxx = rod length in mm/yyy = temperature sensor
combination pH	Type of electrode
pH 0...14	pH measurement range
0... 100 °C (140 °C)	Temperature range for operation (sterilization)
Order No. 52 00X XXX	Order number
Ex SEV14 ATEX 0168 X	Identification and number of certificate
Ex IECEx SEV 14.0025X	Identification and number of certificate
Ex CML 22 UKEX 2108X	Identification and number of certificate
See instruction manual	Instruction (Observe the operating instructions)

In addition, each electrode is allocated a serial number on the red part of the plug-in head to enable identification.

4 Installation and commissioning

1. On unpacking, check the sensor for possible mechanical damage of the pH sensitive glass, the diaphragm, the shaft and the VarioPin (VP) connector. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and briefly rinse the sensor with deionized water. Carefully remove the protection droplet of silicone on the diaphragm with the enclosed knife, as shown on the green pictures on the packing.
3. Again rinse the sensor well with deionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH sensitive glass, since this can lead to electrostatic charging and sluggish response times.
4. Carefully check the space behind the pH sensitive glass for the presence of any air bubbles and remove same by gently swinging the sensor in a vertical plane.
5. Please refer to the instruction manual of the relative housing in order to install the sensor correctly.
6. Connect the sensor to the pH transmitter, using a cable with VarioPin (VP) connector. Please refer to the wiring and terminal diagram accompanying the cable.
7. For electrodes with Pt100 temperature sensors, 3-wire technique is recommended in order to obtain best measurement accuracy. Please refer to the cable and transmitter instruction manuals for the corresponding cable wire connections.
8. For InPro 3200 SG electrodes, with an additional Pt ring for solution ground, connect the corresponding blue wire in the VP-cable to your transmitter in accordance with instructions in the connection diagram of VP-cable and/or in the transmitter instruction manual. For ORP measurement (sep. transmitter needed) connect solution ground with ORP input.

5 Operation

5.1 Calibration of the sensor and pH transmitter

For InPro 3200 sensors, a 2-point calibration is recommended. After removing the watering cap, the sensor is dipped in succession into two different buffer solutions with given pH reference values and the transmitter is calibrated to these buffer values.

5.2 Sterilization of sensor

In the case of sterile processes, such as in fermentation applications, the sensors are sterilized prior to the production cycle. Autoclaving of the whole fermenter/reactor unit including pH electrode is possible with this sensor type. Sterilization can also be effected in-situ by means of steam or superheated process medium in contact with the lower sensor parts.

5.3 Temperature sensor

The integrated Pt100/1000 temperature sensor is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation, and **not** for the accurate indication or control of the process temperature.

6 Maintenance

1. Carefully rinse the sensor tip and diaphragm with deionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt™ electrolyte.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with deionized water or alcohol, and afterwards dry carefully.

7 Trouble-shooting

At sluggish response time (e.g. after a long and dry storage period), reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further sensor parts!) into Reconditioning (HF) Solution (article no. 51 319 001) for 5 min. Afterwards, rinse well with deionized water and let sensor tip and diaphragm soak in reference electrolyte for minimum 12 hours.

Rapidly increased resistance of reference electrode and/or drift could be signs of protein precipitation. Use Cleaning Solution Pepsin/HCl (article no. 51 340 068) to remove such precipitation. At repeated problems of this type, consider switching to sensors with Friscolyt™ based gel electrolyte.

8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Guarantee

On manufacturing defects, 12 months after delivery.

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All rights reserved