# METTLER TOLEDO

## DM145-SC

**Silver ring electrode with Sensor Chip**  
**Silberringelektrode mit Sensor Chip**  
**Electrode combinée à anneau d'argent**  
**avec Sensor Chip**  
**Electrodo con anillo de plata combinado**  
**con Sensor Chip**  
**Электрод с серебряным кольцом**  
**с технологией Sensor Chip**  
**具有智能芯片的金属银环电极**

**Plug & Play combined silver ring electrode**  
with movable sleeve junction for argentometric  
titrations in difficult samples.

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## Care instructions

### Attention

Due to transportation and storage the movable sleeve ring of the sleeve junction may stick to the glass and thus stop the flow of electrolyte. This may cause a very instable or drifting signal and may increase the response time.

### Procedure

1. Close the filling aperture with the rubber cap.
2. Warm the movable sleeve ring with warm water and loosen it by hand by turning.
3. Turn the electrode upside down and allow it to cool (loose sleeve junction).
4. Close the movable sleeve by a slight turning motion.
5. Check the fill level of the electrolyte and refill if necessary.
6. The electrode is ready to use.

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
## 1. Start-up maintenance

- 1.1 To equalize pressure remove rubber cap (6) or perforate its nipple (7).
- 1.2 Refill the reference electrolyte METTLER TOLEDO No. 51340047 to approx. 1 cm below the aperture (6).
- 1.3 Place the electrode in deion. water for 15 minutes to dissolve the  $\text{KNO}_3$  crystals in or on the sleeve junction (3).

## 2. Plug & Play performance

The Sensor Chip which stores the sensor parameter is located in the screw cap of the sensor.

The sensor is recognized after connection to the sensor input of the titrator as well as after each restart of the system by reading the sensor parameter on the chip. A popup appears on the titrator terminal when the sensor was successfully recognized. This popup allows you to directly enter the sensor setup. You find further information concerning the functionality of the Plug & Play sensor.

## 3. Testing

The electrode can be tested by measurement of its potential with a 0.1 mol/L  $\text{AgNO}_3$  and its response time or a titration.

- 3.1 The electrode potential for the 0.1 mol/L  $\text{AgNO}_3$  test solution should be within approx. 250 and 350 mV.
- 3.2 Response time of the electrode: After 30 s the indicated electrode potential should not change more than 2 mV within the following 30 s.
- 3.3 A titer determination of 0.1 mol/L  $\text{AgNO}_3$  is best suited for a test titration. The titration has to show a sharp potential jump of at least approx. 130 mV. The equivalence point should be in the range of 50–200 mV.

## 4. Storage and lifetime

- 4.1 Store the electrode in the reference electrolyte – the sleeve junction has to be immersed – and close the aperture (perforated rubber cap is sufficient).
- 4.2 Never let the electrode get dry, otherwise the  $\text{KNO}_3$ -crystals in or on the sleeve junction have to be dissolved.
- 4.3 If testing a used electrode the requirements indicated under 3.3 should hold true.

## 5. Maintenance Instructions

- 5.1 The tolerable temperature range extends from 0 to 60 °C (32 – 140 °F).
- 5.2 The sample should never enter through the sleeve junction into the reference electrode. Therefore the level of the reference electrolyte must always be higher than the sample solution level.
- 5.3 Do not let the cable connector plug corrode.
- 5.4 Ensure that your sample is acidified with diluted nitric or citric acid. Strictly avoid highly concentrated nitric acid, which leads to corrosion of the silver ring!

## 6. Cleaning procedures

- 6.1 Contamination of the electrolyte: drain out the electrolyte solution and refill with a new one METTLER TOLEDO No. 51340047. Place the electrode for several hours in deion. water.
- 6.2 Blocked sleeve junction by protein: place it at least for 1 hour in pepsin solution (5% pepsin in 0.1 mol/L HCl, e.g. Article No. 51340068).
- 6.3 Contamination of the silver surface: depending on the contamination several cleaning procedures can be carried out:
  - a) Place the electrode for a few minutes in an ultrasonic bath with deion. water or ethanol.
  - b) Place the electrode for 60 s in sulfochromic acid or 10% HNO<sub>3</sub> and rinse it thoroughly with deion. water.
  - c) Slightly polish the silver surface with Al<sub>2</sub>O<sub>3</sub> powder.

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