
InPro 6800 Series O₂ Sensors

Instruction manual Bedienungsanleitung Instructions d'utilisation

InPro 6800 Series O₂ Sensors
52 200 953 C

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InPro 6800 Series O₂ Sensors

Instruction manual

InDip, InFit, InPro, ISM, and InTrac are trademarks of the METTLER TOLEDO Group.

Contents

1	Introduction.....	5
2	Important notes	6
2.1	Notes on operating instructions	6
2.2	Intended use	6
2.3	Safety instructions.....	7
2.4	Environmental protection	9
2.5	Examples of some typical applications.....	9
2.6	Use in Ex-zones	10
2.7	Ex-classification ATEX	11
2.7.1	Introduction.....	11
2.7.2	Rated data	11
2.7.3	Special conditions	12
2.8	Ex-classification FM approved.....	13
2.8.1	Introduction.....	13
2.8.2	Rated data	14
2.8.3	Special conditions	14
2.8.4	Applied standards:.....	15
3	Product description.....	16
3.1	General information	16
3.2	Principle.....	16
3.3	Scope of delivery	16
4	Installation.....	19
4.1	Mounting the sensor	19
4.2	Connection	20
4.2.1	Connecting the InPro 6800 to a cable.....	20
4.2.2	Connecting the VP or AK9 cable to the transmitter.....	20
5	Operation.....	21
5.1	Start-up and polarizing.....	21
5.2	Calibration	22
5.2.1	Purpose of calibration	22
5.2.2	What you have to know for calibration	22
5.2.3	Single point calibration.....	23
5.2.4	Dual point calibration.....	24
6	Maintenance	25
6.1	Inspection of the sensor.....	25
6.1.1	Visual inspection	25
6.1.2	Testing the sensor with the METTLER TOLEDO O ₂ Sensor - Master InPro 6800	26
6.1.3	Testing the sensor via a transmitter	27
6.2	Changing the electrolyte, the membrane body or the interior body	28
7	Storage.....	31
8	Product specification	31
8.1	Certificates	31
8.2	Specifications.....	32
9	Ordering information	33
9.1	Sensors.....	33
9.2	Accessories.....	33
9.3	Spare parts	34
9.4	Recommended transmitters.....	34
9.5	Recommended housings	34
10	Theory of the polarographic sensor.....	35
10.1	Introduction.....	35
10.2	Principle of the design of an oxygen electrode.....	36
10.3	Parameters determining current	36
10.4	Polarization voltage	37
10.5	Temperature	37
10.6	Dependence on flow	38
10.7	Oxygen partial pressure – oxygen concentration.....	39

1 Introduction

Thank you for buying the **InPro™ 6800 sensor from METTLER TOLEDO.**

The construction of the InPro 6800 sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, im proper use could lead to hazards for the user or a thirdparty, and/or adverse effects on the plant or other equipment. **Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro 6800.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Important notes

2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the oxygen sensor. The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies additional **information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended use

METTLER TOLEDO InPro 6800 sensors are intended solely for inline measurement of the oxygen partial pressure in liquids and gases, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.

- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with METTLER TOLEDO Process Analytics in the event of any uncertainties.

2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, «Section 9.3»).

- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.
- When using cables for sensor connection not supplied and recommended by the manufacturer the following standards need to be considered:
 - US National Electrical Code® (ANSI/NFPA 70 [NEC®]), where applicable.
 - Canadian Electrical (CE) Code® (CEC Part 1, CAN/CSA-C22.1), where applicable.
- For guidance on US installations, see ANSI/ISA-SP12.06.01, Installation of Intrinsically Safe Systems for Hazardous (Classified) Locations.
- Control room equipment connected to intrinsically safe associated apparatus in the US should not use or generate more than the specified Um of the associated apparatus.
- Care must be taken during installation to avoid impacts or friction that could create an ignition source.
- Tampering and replacement with non-factory components may adversely affect the safe use of the system.
- Insertion or withdrawal of removable electrical connectors or modules is to be accomplished only when the area is known to be free of flammable vapors.
-  **Warning** – Intrinsically safe apparatus can be a source of ignition if internal spacings are shorted or connections opened.
-  **Warning** – Do not open when an explosive atmosphere is present.
-  **Warning** – Substitution of components may impair intrinsic safety.
-  **Warning** – Substitution of components may impair suitability of the equipment.
-  **Warning** – For connection only to non-flammable processes.
-  **Warning** – To maintain the enclosure IP 66 rating, the connector must be fully engaged.

2.4 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

2.5 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

Measurement in liquids:

- Fermentation
- Yeast propagation
- Wort aeration
- Spring water conditioning
- Storage and processing of fruit juices

Measurement in gases:

- Monitoring of the oxygen limit value to protect products from oxidation.
- Monitoring of the oxygen limit value in the exhaust air of fermentation vessels.
- Monitoring of oxygen limit values during inertization processes.
- Measurement of oxygen concentrations in process gas mixtures.

2.6 Use in Ex-zones

Attention!

For an installation in Ex-zones please read the guidelines following hereafter:

Ex-classification ATEX:

⚡ II 1/2G Ex ia IIC T6...T3 Ga/Gb

⚡ II 1/2D Ex ia IIIC T69 °C...T161 °C Da/Db

Number of the test certificate:

SEV 14 ATEX 0169X

IECEX SEV 14.0026X

CML 22 UKEX 2212X

Ex-classification FM approved:

IS / I, II, III / 1/ ABCDEFG / T6 Ta = 60 °C
- 53 800 002; Entity

Project ID: 3021227

FM certificate number: FM16US0034X

FM18CA0021X

2.7 Ex-classification ATEX

2.7.1 Introduction

According to Directive 2014/34/EU (ATEX 114)¹⁾ Appendix I, InPro 6XXX*/*/*/* oxygen sensors are devices group II, category 1/2G and according to RL 99/92/EG (ATEX 137)²⁾ may be used in zones 0/1 or 1/2 or 1 or 2 and gas groups IIA, IIB and IIC that are potentially explosive due to combustible substances in above stated temperature classes.

For use / installation, the requirements of EN 60079-14 must be observed.

According to Directive 2014/34/EU (ATEX 114)¹⁾ Appendix I, InPro6XXX*/*/*/* oxygen sensors are devices group II, category 1/2D and according to RL 99/92/EG (ATEX 137)²⁾ may also be used in zones 20/21 or 21 or 22 that contain combustible dusts.

The digital sensor circuit is part of a common intrinsically safe system and is for operation connected to a separately certified transmitter.

The digital sensor circuit as part of an intrinsically safe system is isolated from the not-intrinsically safe electric circuits up to a maximum rated voltage of 375 V and from grounded parts up to a maximum rated voltage of 30 V.

2.7.2 Rated data

Measuring circuit

Fail-safe ignition protection class Ex ia IIC only for connection to a certified fail-safe circuit.

Analog O₂ oxygen sensor maximum values

$$U_i \leq 16 \text{ V}$$

$$I_i \leq 190 \text{ mA}$$

$$P_i \leq 200 \text{ mW}$$

$$L_i = 0 \text{ (effective internal inductance)}$$

$$C_i = 900 \text{ pF (effective internal capacitance)}$$

Note: The above maximum values are each the total of all individual circuits of the associated intrinsically safe power supply and transmitter.

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR)

2.7.3 Special conditions

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1G applications, zone 0, is shown in the following table:

Temperature class	Max. ambient or media temperature
T 6	68 °C
T 5	80 °C
T 4	108 °C
T 3	160 °C

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1D applications, zone 20, is shown in the following table:

Temperature class	Max. ambient or media temperature
T 69 °C	68 °C
T 81 °C	80 °C
T 109 °C	108 °C
T 161 °C	160 °C


- The capacitance and inductance of the connecting cable has to be considered.
- The oxygen sensor type InPro 6XXX can be used in / with the fittings InFit™ 76*-*** or InTrac™ 7**-***, or in / with other suitable fittings in potentially explosive areas.
- The metal body of the oxygen sensors, or the fittings InFit 76*-*** or InTrac 7**-***, or other appropriate fitting is optionally included in the routine pressure test of the system.
- The independent fitting used for installation of oxygen sensor must be conductively connected <to the equipotential bonding system.

2.8 Ex-classification FM approved

2.8.1 Introduction

The following FM control drawing and the standards listed in section 2.8.4 must be observed, where applicable:

2.8.2 Rated data

For rated data please observe section 2.7.2.

2.8.3 Special conditions

In type of protection intrinsically safe apparatus, the analytical probes (InPro 6000 Series with InFit 76X/Y and InTrac 7XX/YY Series Housings) equipment is designated with the following specific conditions of use.

1. The dissolved oxygen sensor shall be installed in compliance with the enclosure, mounting, spacing and segregation requirements of the ultimate application, including a tool removable cover, and is suitable for use with FM Approved InFit 76X/Y and InTrac 7XX/YY Series housings.
2. Process temperature no greater than +130 °C.
3. Maximum permissible working pressure is 12 bar (174 psig).
4. Potential Electrostatic Charging Hazard – To prevent the risk of electrostatic sparking, the non-metallic surface should only be cleaned with a damp cloth.
5. Enclosures containing titanium constitute a potential risk of ignition by impact or friction. Care must be taken into account during installation and use to prevent impact or friction.

2.8.4 Applied standards:

United States Standards

Title	Number	Issue Date
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements.	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations.	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use.	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements.	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety "i".	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	ANSI/UL 61010-1:2018	2018

Canadian Standards

Title	Number	Issue Date
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0:2019	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety "i".	CAN/CSA-C22.2 No. 60079-11-14:R2018	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1-12:R2017	R2017

3 Product description

3.1 General information

The **sensor InPro 6800** with integrated temperature probe is used for oxygen measurement.

The sensor is **sterilizable** and **autoclavable** and **compatible with CIP** (cleaning in place).

InPro 6800 sensors with ISM™ functionality offer Plug and Measure as well as enhanced diagnostics features. ISM sensors are available with K8S (fully digital) or VP connector.

3.2 Principle

Here is a short summary of the principle of polarographic measurement on which this sensor is based (Clark 1961).

- a) The Clark polarographic sensor basically consist of a working electrode (cathode), a counter/reference electrode (anode), and an oxygen permeable membrane which separates the electrodes from the sample medium.
- b) The transmitter supplies a constant polarization voltage to the cathode, needed to reduce oxygen.
- c) The oxygen molecules which migrate through the permeable membrane are reduced at the cathode. At the same time, oxidation takes place at the anode and oxidized anode metal (silver) is liberated as silver ions into the electrolyte. The electrolyte closes the electric circuit between anode and cathode (ion conductivity).
- d) The current produced by the reactions described above is measured by the transmitter and is proportional to the partial pressure of oxygen (pO₂) in the sample medium.

Please refer to Section „10 Theory of the polarographic sensor“ for further information.


3.3 Scope of delivery


Each sensor is supplied fully assembled and factory-tested for correct function together with:

- a quality control certificate
- inspection certificates 3.1 (complying with EN 10204)

See electrolyte order details at section „9.3 Spare parts“

Equipment features

12 mm Sensor

25 mm Sensor

4 Installation

4.1 Mounting the sensor

Important! Remove the green protective cap before mounting the sensor.

Mounting the sensor in a housing



Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

Mounting the sensor directly on a pipe or a vessel

The sensor can be mounted vertically or in an inclined position. In the case of inclined mounting, the angle of the housing must be equal to or greater than 15° above the horizontal.

Caution! Installation outside the admissible range of mounting positions is not allowed, otherwise correct operation of the sensors / electrodes may be impaired.


The **12 mm sensors** can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.


The **25 mm sensors** can be mounted directly through a standard weld-in socket or the safety weld-in socket from METTLER TOLEDO.

4.2 Connection


4.2.1 Connecting the InPro 6800 to a cable

 Sensors with ISM functionality require the use of a special AK9 cable as well as an ISM compatible O₂ transmitter. The sensor is connected to the transmitter via a AK9 cable. The AK9 cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP68 connector housing guarantees maximum process safety. If you have chosen to use a sensor with digital ISM functionality, the K8S connector will incorporate a preamplifier. It is absolutely necessary to protect this electronic component against any electrical charge

 **Do not touch the sensor at the VP connector plug.** To connect the VP cable to the sensor align the slit of the VP connector with the pin in the plug. Then tightly screw the plug to fasten the two parts.

4.2.2 Connecting the VP or AK9 cable to the transmitter

 **Note:** Cable assignment can be found in the METTLER TOLEDO VP cable instruction manual.

 **Note:** For connecting the cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual.

5 Operation

Attention! Before using the sensors for the first time, the electrolyte should be replaced (see «Chapter 6.2»). Due to possible adverse conditions during transport and storage (e.g. airfreight; pressure and temperature variations), the quality of the electrolyte may become impaired. Poor electrolyte quality can lead to erroneous measurement values.

5.1 Start-up and polarizing

Attention! The protective cap at the tip of the sensor must be removed before putting the sensor into operation.

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O₂ Sensor-Master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O₂ transmitter or to a sensor master. After 6 hours, the sensor is fully polarized and ready for operation.

A shorter polarization period is sufficient if the sensor has been disconnected for only a few minutes. The following table serves to establish the correct polarization time in relation to the depolarization time.

Depolarization time ¹⁾ t_{depol} [Min.]	Minimum required polarization time ²⁾ [Min.]
$t_{\text{depol}} > 30$	360
$30 > t_{\text{depol}} > 15$	$6 * t_{\text{depol}}$
$15 > t_{\text{depol}} > 5$	$4 * t_{\text{depol}}$
$t_{\text{depol}} < 5$	$2 * t_{\text{depol}}$

¹⁾ Depolarization time:

Time span in which the polarization voltage is cut off from the sensor. This is the case during:

- change of electrolyte
- change of membrane body
- the time the cable is disconnected or no transmitter or sensor master is connected to the cable

²⁾ Polarization time:

Time span during which the sensor is under a polarization voltage.

Attention! Setting of the polarization voltage on the transmitter for correct measurements:

- Standard applications (e.g. measurement in biotechnology): **– 675 mV**
- Measurement of permanently low oxygen concentrations < 500 ppb in the presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries): **– 500 mV**

Important! To ensure the supply of the correct polarization voltage the transmitter must be set accordingly.

5.2 Calibration

5.2.1 Purpose of calibration

Each oxygen sensor has its own individual slope and own individual zero point. Both values are subject to change, for example, through electrolyte consumption or after exchange of electrolyte or membrane body. To ensure high measurement accuracy of the sensor, a calibration must be carried out after each change of electrolyte or membrane. Prior to calibration, the sensor has to be polarized for at least 6 hours.


 **A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.**


 **To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100%. If not, then the sensor needs a new calibration.**

5.2.2 What you have to know for calibration

 **Attention!** If dissolved oxygen is being measured at low concentrations (< 500 ppb) in the presence of volatile acidic components (e.g. carbon dioxide (CO₂) measurements in breweries), the calibration procedure described below should be followed in order to achieve the best measurement performance.

No special calibration procedure is required if dissolved oxygen is being measured in standard applications (e.g. measurement in biotechnology).

 **First step:** Adjust the transmitter polarization voltage to – 675 mV. Detailed information is given in the transmitter instruction manual.

 **Second step:** After adjusting the polarization voltage to – 675 mV, equilibration time is required. Wait for about 5 minutes before carrying out step 3.

 **Third step:** Perform calibration according to the transmitter' instruction manual.

 **Fourth step:** Adjust the transmitter polarization voltage back to – 500 mV.

 **Steps 1 to 4** can be accomplished automatically if you use the M 700 transmitter. Thus the routine can be minimized to a few keystrokes. In addition the M 700 transmitter must be equipped with the software function SW-700-011 «high CO₂ compensation». Please ask your local METTLER TOLEDO distributor.

 **General remarks:**

- **For calibration in air, the sensor membrane must be dry**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the **oxygen saturation index** of the calibration medium is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is only very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.
- For correct calibration, a minimum flow rate of the calibration medium is necessary.
- **Calibration in a fermenter should be performed after sterilization** (as sterilization may alter the sensor slope), but prior to inoculation. If it is not possible to perform the calibration after sterilization, the use of an existing membrane body that has been pre-sterilized while mounted in the sensor is recommended. A slope alteration of some per cent can occur with new membrane bodies, particularly after a first sterilization, as the tension of the membrane may be altered by the sterilization process.
- Make sure that all other parameters, such as temperature and pressure, are constant.

For continuous applications, we recommend **periodic recalibration in line with your requirements on accuracy, the type of process in operation and your own experience**. The frequency of the need for recalibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

5.2.3 Single point calibration

By carrying out a single point calibration, the factual slope of the sensor can be established. The calibration medium can be water with known oxygen saturation index (e.g. air-saturated water) or air with known water-vapor saturation (e.g. water-vapor saturated air).

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air,

20.95 % O₂, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Single point calibration should be sufficient for almost all process applications.

5.2.4 Dual point calibration

By carrying out a dual point calibration both slope and zero point of the sensor can be established.

Important! In case of a **dual point calibration, always start by the zero point calibration** before calibrating the slope.

Due to the very low zero current of METTLER TOLEDO sensors, dual point calibration is normally not necessary for standard applications.

As a rule, the zero point should be adjusted to zero manually, or it is automatically performed by the transmitter (see instruction manual for the transmitter).

A zero point calibration is only advisable if very high accuracy is required at low oxygen concentrations.

Attention! Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.995 %.

After the sensor signal has stabilized (after 20...30 minutes), the sensor can be calibrated through the relevant transmitter to the 100 % value of the desired measurable variable, e.g. 0 % air, 0.0 % O₂, or 0.0 ppm (see instruction manual for the transmitter).

6 Maintenance

6.1 Inspection of the sensor

6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

1. The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
2. Check the cable for buckling, brittle areas or ruptures.
3. Before calibration always examine the membrane foil optically for signs of damage. The foil must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.

Note: An undulated membrane has no influence on the sensor performance, assuming the membrane is intact.

4. The membrane body must be replaced if the sensor has too long a response time, the reading is unstable or subject to drift, and if the sensor cannot be calibrated or the membrane shows sign of mechanical damage.
5. Check the cathode area for discoloration, contamination or cracks in the glass. If necessary rinse with demineralized water and clean with a clean soft brush or soft paper tissue.

Attention! Do not use any cleaning agents containing alcohol. This could damage the sensor or lead to fault current.


Attention! The glass body is fragile and sensitive to vibration.

6.1.2 Testing the sensor with the METTLER TOLEDO O₂ Sensor - Master InPro 6800

We recommend to use the METTLER TOLEDO O₂ Sensor-Master InPro 6800 to check the quality of your sensor as follows (not applicable for digital ISM sensors):

- Connect the sensor to the O₂ Sensor-Master.

As soon as the sensor is connected to the O₂ Sensor-Master, the polarization function is automatically activated. Please note: if the sensor was disconnected from the transmitter for longer than 5 minutes, the sensor must be polarized first (polarizing time see «Section 5.1») to get representative test results.

- **Battery Check:**
Push the switch to the left. If the battery is ok and the O₂ Sensor-Master is operational the green LED lights up. Otherwise, please consult the instruction manual of the O₂ Sensor-Master.
- **Sensor Check:**
For this test the O₂ sensor must be **fully polarized** and the **membrane of the sensor must be dry and clean.**

Expose the sensor (connected to O₂ Sensor-Master) to the air. By pushing the switch to the right to the position "2 – Sensor check", the O₂ Sensor-Master checks whether the electrode current for air measurement delivered by the sensor is within the admissible range (40 to 100 nA for InPro 6800). If the green LED lights up the current for air measurement is within the admissible range.

If the LED does not light up, you should check the battery of the O₂ Sensor-Master (see instruction manual «Accessories»). If the battery is working, than there is probably a problem with your sensor. You should change the electrolyte and/or the membrane body of your sensor. If after a membrane change the LED still does not light up, this means that there is maybe something wrong with the interior body of the sensor. You should than change it (see «Section 6.2»).

Important! The Sensor Check function only verifies the correctness of the electrode current for air measurement. In order to be absolutely sure of the functionality of the sensor, the residual signal in an oxygen free medium should also be controlled (see «Section 6.1.3»).

6.1.3 Testing the sensor via a transmitter

A periodic zero current measurement (**no zero point calibration**) is recommended for verification of proper sensor function.

Note: At the time you carry out the zero current measurement, the sensor must be polarized.

Zero current measurement can be done by using zeroing gel (order no. 30 300 435) or nitrogen (N₂) or carbon dioxide (CO₂) calibration gases with a purity of at least 99.995 %, alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 10 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If the measured values are too high, this suggests a depleted electrolyte or a defective membrane. In the first instance replace the electrolyte, and in the second case exchange both the membrane body and the electrolyte accordingly.

If after such procedures the above mentioned values are still not reached, replace the interior body. If this doesn't solve the problem too send the sensor to your local METTLER TOLEDO representative for inspection.

Many sample media contain volatile substances which, even at very low concentrations, have a clearly perceptible smell. Similarly to oxygen, these substances are able to invade the electrolyte through the gas-permeable membrane. Accordingly, they become noticeable when changing the electrolyte. In most cases, such substances have absolutely no influence on the measuring properties of the sensor. Slight discoloration of the electrolyte also has no effect on the measuring properties.

6.2 Changing the electrolyte, the membrane body or the interior body

METTLER TOLEDO DO sensors are supplied with fitted membrane body and have been checked for proper function.

However, if a sensor is to be stored for several months, the electrolyte should be replaced before use.

If the membrane exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body has to be replaced.

Warning! The O₂ electrolyte has an alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.**

If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.

When changing the electrolyte, the membrane body or the interior body, please observe the following instructions (see also the following illustration):

Attention! Make sure that this maintenance step is carried out in **clean place**.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.

Note: steps 4 to 7 may only be carried out when changing the interior body.


4. Unscrew the retainer nut of the interior body with an adjustable wrench or with a 3/8" wrench.
5. Remove the interior body by pulling it out of the sensor shaft. If necessary use a plier.

Warning! Do not twist the interior body.

Otherwise the connection pins can be damaged.

6. Insert the new interior body in the sensor shaft. Turn the interior body in the shaft until the slit of the interior body is aligned with the pin placed in the shaft.
7. Press the body in the shaft and screw the new retainer nut in place.
8. Examine the O-rings visually for mechanical defects, and replace if necessary.
9. Half-fill the new membrane body with O₂ electrolyte and make sure that all bubbles are removed. Air bubbles can be removed by carefully tapping on the membrane body.


10. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue

Attention! No electrolyte, sample media or contamination may be present between the membrane body and the cap sleeve. Please check carefully.

11. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
12. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated.

7 Storage

The sensor can be stored for several months filled with O₂ electrolyte (order no. 30 298 424) as long as the protective cap is in place. In order to avoid the 6-hour polarization time before reuse of the sensor, it can be kept stored connected to the METTLER TOLEDO O₂ Sensor-Master InPro 6800.

Attention! The electrolyte should be replaced before use if the storage period exceeds 3 months.

If the storage period exceeds 6 months, the sensor should be **stored dry**, i.e. without any electrolyte in the membrane body.

Attention! A sensor being stored dry (without electrolyte in the membrane body) **may on no account be connected to the O₂ Sensor-Master InPro 6800**.

8 Product specification

8.1 Certificates

Each sensor is delivered with a set of **3.1 certificates** (complying with EN 10204).

All wetted metal parts (sensor shaft, cap sleeve and membrane body) are identified with a engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and membrane body) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

8.2 Specifications


Attention! Please refer to the supplementary sheet 52 201 146 if the sensor, together with the membrane body T-6800 gas, is been used for measurement of oxygen partial pressure in gases. This supplementary sheet is only included in the packaging of T-6800 gas membrane bodies and kits.

InPro 6800	
Measurement principle	Amperometric/Polarographic
Working conditions	
Pressure resistance measurement	0.2...6 bar absolute [2.9...87.0 psi absolute]
Mechanical pressure resistance	max. 12 bar absolute [174.0 psi absolute]
Measuring temperature range	0...80°C [32...176°F]
Temperature range	- 5...140°C [23...284 °F] (sterilizable/autoclavable)
Construction	
Temperature compensation	Automatic with built-in RTD
Cable connection	VarioPin (IP 68) straight or angled K8S straight (digital ISM sensors)
O-ring material	Silicone FDA and USP VI approved
Membrane material	PTFE/Silicone/PTFE (reinforced with steel mesh)
Wetted metal parts	Stainless steel Special material on request
Surface roughness of wetted metal parts (ISO 1320:1992)	N5 ($R_a < 0.4 \mu\text{m}$ [16 μin])
Quick disconnect interior body	Standard
Cathode	Pt
Anode	Ag
Guard ring	No
Dimensions	
Sensor diameter	12 or 25 mm (0.47" or 0.98")
Immersion length (a) for 12mm sensor	70, 120, 220, 320, 420 mm (2.8, 4.7, 8.66, 12.6, 16.54")
Immersion length (a) for 25mm sensor	80,160, 260, 360 mm (3.15, 6.3, 10.24, 14.17")
Performances	
Detection limit	6ppb
Accuracy	$\leq \pm$ [1% + 6 ppb] of reading in liquids
Response time at 25°C/77°F (air → N ₂)	98% of final value < 90 s
Sensor signal in ambient air (25°C/77°F)	50...110 nA
Residual signal in oxygen-free medium	< 0.1% of the signal in ambient air
Maximum flow error	\leq 5%
Certification	
EHEDG, 3A	Yes
3.1 (EN 10204)	Yes
ATEX, IECEx und UKEX certificate	Yes
FM Approval	Yes
FDA / USP VI	Yes
Quality control	Yes
Compatibility	
with METTLER TOLEDO transmitters	see «Section 9.4»
with METTLER TOLEDO housings	see «Section 9.5»

9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

9.1 Sensors

Sensor without ISM functionality

Immersion length (a)	12 mm diameter		25 mm diameter	
	straight	angled	straight	angled
70 mm (2.8")	52 200 964	52 200 969	—	—
80 mm (3.2")	—	—	52 200 974	52 200 978
120 mm (4.7")	52 200 965	52 200 970	—	—
160 mm (6.3")	—	—	52 200 975	52 200 979
220 mm (8.7")	52 200 966	52 200 971	—	—
260 mm (10.2")	—	—	52 200 976	52 200 980
320 mm (12.6")	52 200 967	52 200 972	—	—
360 mm (14.2")	—	—	52 200 977	52 200 981
420 mm (16.5")	52 200 968	52 200 973	—	—

9.2 Accessories

Accessories	Order no.
O ₂ Sensor-Master InPro 6800	52 200 892
O ₂ Sensor-Simulator	52 200 891
Adapter T82 socket – VP plug	52 200 939
Adapter VP socket – T82 plug	52 200 940
VP cable VP6-ST/ 3 m	52 300 108
VP cable VP6-ST/ 5 m	52 300 109
VP cable VP6-HT/ 3 m	52 300 112
VP cable VP6-HT/ 5 m	52 300 113
VP cable VP8-ST/ 3 m	52 300 354
VP cable VP8-ST/ 5 m	52 300 355
VP cable VP8-HT/ 3 m	52 300 361
VP cable VP8-HT/ 5 m	52 300 362
Oxygen zeroing gel (3 × 25 ml)	30 300 435

For other cable lengths or types, please contact your local METTLER TOLEDO representative.

9.3 Spare parts

Spare part	Order no.
Membrane body, single T-96	52 200 071
Membrane kit T-96 (4 membrane bodies, 1 O-ring set, 25 ml Electrolyte bottle)	52 200 024
O ₂ Membrane bodies T-96 (16 pcs.)	52 206 114
O ₂ Electrolyte pack (3 × 25 ml)	30 298 424
Interior body InPro 6800 (with quick disconnect)	52 200 899
Cap sleeve N-type	52 200 037
Cap sleeve with protective cage P-type	52 200 038

9.4 Recommended transmitters

Note: Ask your local distributor or refer to the technical data sheet.

9.5 Recommended housings

Housings	12 mm Ø	25 mm Ø	
Static housings			Please ask your local METTLER TOLEDO representative.
INGOLD "safety weld-in socket"	–	✓	
InFit 761 CIP	✓	–	
Retractable housings			
InFit 777 e	✓	–	
InFit 797 e	✓	–	
Immersion housing			
InDip™ 550	✓	–	

Note: The housings are available in different versions. Please contact your distributor to get the right ordering information.

10 Theory of the polarographic sensor


10.1 Introduction

Two types of electrodes are employed in analytical work: **potentiometric** and **amperometric** electrodes.

- Potentiometric electrodes develop a voltage generated by the activity of a particular ion. Examples of such electrodes are glass electrodes (like pH electrodes) and most ion-selective electrodes. Their individual potentials cannot be determined. The measurable quantity is the difference of potential between the measuring electrode and an inert reference electrode. The potential of the reference electrode must be constant.
- All potentiometric electrodes are subject to Nernst's law and for this reason electrodes and measuring instruments are in most cases interchangeable. An important requirement of potentiometric measurements is the virtually currentless determination of the electrode voltage. During measurement no chemical reaction occurs and the solution remains in equilibrium.
- In the case of **amperometric electrodes**, such as the **oxygen electrode**, activity measurement is based on a current measurement.

The oxygen electrode consists of a cathode and an anode conductively connected by an electrolyte. A suitable polarization voltage between the anode and the cathode selectively reduces the oxygen at the cathode.

Cathode reaction:

Anode reaction:

These chemical reactions result in an electric current which is proportional to the oxygen partial pressure (PO_2). The oxygen electrode consumes oxygen which is continuously extracted from the solution. The viscosity and flow rate of the solution are therefore important parameters.

The electrode current of an oxygen electrode is determined not only by the oxygen partial pressure but by many other electrode parameters. The electrode currents of different electrode types may differ by several powers of ten. For this reason oxygen electrodes and amplifiers cannot be freely interchanged.


10.2 Principle of the design of an oxygen electrode

There are two main types of oxygen electrodes:

- Electrodes **without** membrane
- Electrodes **with** gas-permeable membrane (Clark Principle)

The membrane electrode according to Clark is today most widely used. As compared to the electrode without membrane it possesses the following advantages:

- Oxygen measurement in gases and solutions
- No mutual contamination of electrode and solution
- No or little dependency on flow

In the case of the Clark electrode, geometrical configuration is very important. In particular, the thickness of the electrolyte film between the cathode and the membrane must be within narrow tolerances so as to ensure good linearity and a low zero current (current in a nitrogen atmosphere). The figure above exemplifies the principle of the design of a Clark-type oxygen electrode.

10.3 Parameters determining current

The quantity of oxygen diffused in and the magnitude of the electrode current are influenced by the following parameters:

- Oxygen partial pressure of the solution
- Membrane material and thickness
- Size of cathode
- Polarization voltage
- Temperature
- Flow conditions in the solution

Fick's law gives the mathematical relationship between these parameters:

$$I = k \times D \times a \times A \times \frac{pO_2}{X}$$

I = electrode current

k = constant

D = coefficient of diffusion of O₂ in the membrane

a = solubility of O₂ in the membrane material


A = cathode surface

pO₂ = oxygen partial pressure in the solution

X = thickness of gas-permeable membrane

10.4 Polarization voltage

The voltage between the anode and the cathode is so selected that the oxygen is fully (> A, see polarogram) reduced while other gases are unaffected (< D). The ideal voltage for the Pt/Ag/AgCl system is between -500 and -750 mV.

The polarization voltage should remain as constant as possible. Besides a constant source of voltage the following prerequisites must be met:

The electrical resistance of electrolyte film must not exceed a particular value so as to avoid a voltage drop. The anode must have a large surface so as to prevent polarization of the anode by the electrode current.

10.5 Temperature

The temperature dependence of the current passing through an oxygen electrode when referred to a constant O₂ partial pressure is determined mainly by the properties of the gas-permeable membrane.

10.6 Dependence on flow

With most oxygen electrodes the electrode current is smaller in stagnant solutions than in agitated ones. The consumption of oxygen by the electrode results in extraction of oxygen from the solution in the close vicinity of the cathode outside the membrane. The oxygen is replaced by diffusion. If the electrode current is strong, the solution cannot fully restore the oxygen by diffusion. This results in an electrode current weaker than would correspond to conditions in the solution. In agitated solutions the oxygen is conveyed to the surface of the membrane not only by diffusion but additionally by the flow (convection). In that case no oxygen impoverishment occurs at the membrane surface.

A high degree of flow dependence occurs mainly with large cathodes, thin and highly permeable membranes, i.e. where electrode currents are large. The problem of flow dependence is often solved by prescribing a minimum flow rate.

In METTLER TOLEDO InPro 6800 electrodes, the thin PTFE membrane determining the electrode current (i.e. the actual measuring signal) is separated from the sample solution by a relatively thick silicone membrane. This latter is highly permeable to oxygen molecules and thus acts as an oxygen reservoir. The diffusion of oxygen out of the sample solution into the silicone membrane is spread over a wide area. Since this results in less oxygen being extracted from the sample solution per unit area, the double PTFE/Silicone membrane forms an effective buffer against disturbances due to hydrodynamic flow.

10.7 Oxygen partial pressure – oxygen concentration

The electrode current depends on the oxygen partial pressure and the O₂ permeability of the membrane – but not on the O₂ solubility in the solutions. The oxygen concentration in mg O₂/l (CL) cannot therefore be determined directly with an electrode.

According to Henry's law the oxygen concentration is proportional to its partial pressure (PO₂).

$$CL = pO_2 \times a$$

a = solubility factor

If " a " is constant, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

Medium, sat. with air	Solubility at 20 °C (68 °F) and 760 mm Hg
Water	9.2 mg O ₂ /l
4 mol/l KCl	2 mg O ₂ /l
50 % Methanol-water	21.9 mg O ₂ /l

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions.

Thus, determination of the oxygen concentration is only possible with constant and known solubility factors " a ".

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

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US Normen

Titel	Nummer	Ausgabe-Datum
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements.	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations.	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use.	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements.	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	ANSI/UL 61010-1:2018	2018

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Titel	Nummer	Ausgabe-Datum
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0:2019	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	CAN/CSA-C22.2 No. 60079-11- 14:R2018	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1- 12:R2017	R2017

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Titre	Numéro	Dates de publication
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements.	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations.	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use.	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements.	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	ANSI/UL 61010-1:2018	2018

Normes canadiennes

Titre	Numéro	Dates de publication
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0:2019	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	CAN/CSA-C22.2 No. 60079-11- 14:R2018	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1- 12:R2017	R2017

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All rights reserved