METTLER TOLEDO

InPro 8000 Series

Turbidity sensors
Trübungssensoren
Sondes de turbidité

Instruction manual
Bedienungsanleitung
Instructions d'utilisation

English page 3

German page 23

French page 43

InPro 8000 Series

Turbidity sensors

Instruction manual

Subject to technical changes without prior notice.

© It is strictly forbidden to reprint this instruction manual or any parts thereof without the written permission of Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, 8902 Urdorf, Switzerland. No section or excerpt whatsoever may be reproduced or with the assistance of electronic systems be edited, duplicated or distributed, in particular in the form of photocopies, photographs, magnetic media or other recording methods. All rights reserved, especially the right of duplication and translation as well as in regard to patent and registration rights.

InFit, InTrac, InPro, InDip and CaliCap are trademarks of the METTLER TOLEDO Group.

Viton and Kalrez are registered trademarks of DuPont Performance Elastomers LLC.

All other trademarks are the property of their respective holders.

Contents

1. Product description	5
1.1 Introduction	5
1.2 Equipment and scope of delivery	6
1.3 Technical data	6
2. Safety	9
2.1 Application compatibility	9
2.2 Proper utilization	9
2.3 Safety measures	9
2.4 Ex-Declaration	10
2.4.1 Use in Ex-zones	10
2.4.2 Ex-Classification	11
3. Installation and start up	12
3.1 Handling	12
3.2 Installation in vessels/reactors	12
3.3 Installation in pipes	12
3.4 Fiber optic cable extensions	14
3.5 Calibration in open laboratory vessels	14
4. Maintenance and troubleshooting	15
4.1 Conditions of warranty	15
4.2 Maintenance	15
5. Removal from operation, storage, disposal	16
5.1 Removal from operation	16
5.2 Storage	16
5.3 Disposal	16
6. Dimensional drawings	17
7. Ordering information	18
7.1 Sensors	18
7.2 Accessories	19
8. Certificates	21
8.1 Material certificate	21

1. Product description

1.1 Introduction

We thank you for your decision to purchase a METTLER TOLEDO InPro™ 8000 Series fiber optic measuring sensor with operation based on the principle of backward scattered light.

Please keep this instruction manual in a safe place where it is readily to hand for users of this sensor. If you have any questions regarding function of the sensor and ancillary equipment, or if you consider that the descriptions/instructions do not tally or are lacking, please contact your nearest METTLER TOLEDO Authorized Distributor, who will be pleased to provide assistance.

Markings

The markings labeled on each sensor state the following information:

METTLER TOLEDO

Serial no. xxxxxx Order no. 52 xxx xxx

The code following the model designation of the InPro 8050, InPro 8100 and InPro 8200 Sensors as printed on the package label provides information on the shaft material, sealing material for the optical window (if present), and the length of the shaft in mm.

e.g. InPro 8200/S/Kalrez®-FDA/205

with S = stainless steel shaft DIN 1.4435 or 1.4404 (316L)
Kalrez-FDA 6230 window seal
shaft length = 205 mm (8.07").

In addition:

H = Hastelloy C-22 shaft

Epoxy = Epoxy bonded window

See also chapter 7.1: "Ordering information sensors".

1.2 Equipment and scope of delivery

Scope of delivery

The sensors of the InPro 8000 Series are delivered ready for use including fiber optic cable. Each sensor is accompanied by this instruction manual.

A METTLER TOLEDO Process System for the measurement of turbidity or of suspended particles is completed by a METTLER TOLEDO Turbidity Transmitter Type M800 1-channel (order no. 30 026 633). Detailed information for calibration of the measuring system can be found in the instruction manual of the M800 transmitter (order no. 52 121 825).

Packaging

The packaging consists of cardboard and foamed plastic. Safekeep the packaging for later use during storage or transportation of the sensor. Should you wish to dispose of the packaging material, please observe your local regulations as well as those data and instructions given in **chapter 5.3** of this manual.

Unpacking and inspection


Please check the sensors immediately during unpacking in order to determine possible damage or missing items. Any irregularities should immediately be notified to your carrier and to your supplier.

1.3 Technical data

Function

Sensors of the InPro 8000 Series are optical sensors which measure the concentration of suspended particles in a liquid medium on the basis of the principle of backscattered light. The Transmitter M800 1-channel is equipped with a LED which beams a light – almost in the infra-red range (wavelength 880 nm) – via a fiber optic cable into the liquid medium. If the light hits on particles, then it is scattered in all directions. The light scattered at an angle of 180° (backscattered light) is captured and led back via a fiber optic cable to a photodiode in the transmitter, processed as photocurrent and the signal transformed into a value for the transmitter display.

The sensors are divided into two types, single and dual optical fiber sensors. The difference relates to the number of optical fibers present at the end (nose) of the sensor shaft.

At low particle concentrations, the measuring range is expressed in FTU (Formazin Turbidity Units). At high concentrations, the measurements are given in grams of suspended solids per liter of sample volume, using diatomaceous earth as reference substance.

Specifications InPro 8000 series turbidity sensor

	InPro 8050	InPro 8100	InPro 8200
Measuring principle	backscatter, 1-fiber	backscatter, 1-fiber	backscatter, 2-fiber
Sensor material	PSU (Polysulfone)	DIN 1.4435 (316L)	DIN 1.4435 (316L) or HA-C22
Sensor diameter	12 mm	12 mm	12 mm
Insertion length	120 mm	120, 205, 297 or 407 mm	120, 205, 297 or 407 mm
Process connection	Pg 13.5	Pg 13.5	Pg 13.5
Measuring range:			
Formazin Turbidity Units	10...4000 FTU	10...4000 FTU	5...4000 FTU
Suspended Solids (Diatomaceous earth as reference)	0...250 g/L	0...250 g/L	0...30 g/L
Design:			
Cable connectors	SMA	SMA	SMA
Cable lengths	6 m (20 ft) optical	3 m (10 ft) optical	3 m (10 ft) optical
Cable type	duplex HCS fiber, fixed	duplex HCS fiber, fixed	duplex HCS fiber, fixed
Wetted parts:			
Metals	none	DIN 1.4435 (316L)	DIN 1.4435 (316L) or Hastelloy® (HA-C22)
Plastic	PSU	none	none
O-ring	Viton®-FDA	Viton-FDA	Viton-FDA
Sapphire spigot for fiber protection	yes	none	none
Sapphire window for fiber protection	none	none	Kalrez®-type: Kalrez-FDA 6230 sealed Epoxy-type: Epoxy bonded
Working conditions:			
Pressure range	0...2 bar (0...29 psi)	0...6 bar (0...87 psi)	Epoxy-type: 0...16 bar (0...232 psi) With Swagelok® adapter: 60 bar (0...870 psi) Kalrez-type: 0...16 bar (0...232 psi)
Measuring temperature range	0...60 °C (32...140 °F)	-30...130 °C (-22...266 °F)	Epoxy-type: -30...130 °C (-22...266 °F) Kalrez-type: -10...130 °C (14...266 °F)
Temperature range (sterilization)	n/a	0...130 °C (32...266 °F)	Kalrez-type: 0...130 °C (32...266 °F) Epoxy type: n/a
Certificate and Approval:	n/a	☞ Ex II 1/2G Ex op is/ op pr IIA 95 °C Ga/Gb ☞ Ex II 1/2G Ex op is/ op pr IIB 38 °C Ga/Gb	☞ Ex II 1/2G Ex op is/ op pr IIA 95 °C Ga/Gb ☞ Ex II 1/2G Ex op is/ op pr IIB 38 °C Ga/Gb

2. Safety

2.1 Application compatibility

The wetted material parts of the sensor (shaft, window, window sealing material, O-ring) can under circumstance be non-compatible with the particular composition of the process medium and/or of the operating conditions. Responsibility to verify application compatibility lies wholly with the user.

The compatibility of different types of material are outlined on

<http://www.coleparmer.com/techinfo/chemcomp.asp>.

Mettler-Toledo GmbH, Process Analytics accepts no responsibility whatsoever for the correctness or accuracy of such details.

2.2 Proper utilization

METTLER TOLEDO InPro 8000 Series sensors are intended solely for the measurement of suspended particles or turbidity in aqueous solutions in industrial applications.

Any other use, or any operation over and above that intended by the manufacturer, are deemed to be non-permissible and incorrect, and can lead to harm or injury to material/equipment and persons. This is also relevant for applications which do not comply with the technical data of the sensor. For any damage possibly arising from such misuse, the user assumes full and sole responsibility.

2.3 Safety measures

The sensors InPro 8000 Series have been manufactured in line with state-of-the-art technology and in accordance with accepted technical safety regulations. Nevertheless, the sensors can still represent a source of risk and danger:

- if the sensors are operated by insufficiently trained personnel,
- if the sensors are employed incorrectly or not as intended by the manufacturer
- if the sensors are not regularly maintained or serviced.

Local legislation and regulations must be observed at all times. Such stipulations do not form an integral part of this instruction manual.

It is on principle necessary for persons handling or using the sensors to wear personal safety outfit such as protective goggles and protective clothing.

The user is responsible for the instruction and training of personnel. In this respect, additional copies of the instruction manual can be ordered from your supplier. This instruction manual is an essential element of the sensor equipment and must at all times be readily at hand for operators directly at the location of employment of the sensors.

Before the sensor is removed from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

No modification whatsoever may be carried out on the sensors. Any unauthorized modification or manipulation of the sensors results in immediate expiry of the full scope of warranty granted by the manufacturer.

2.4 Ex-Declaration

2.4.1 Use in Ex-zones

Note: Before installation in Ex-zones, please read and observe the following guidelines:

Ex-Classification:

⊕_{Ex} II 1/2G Ex op is/ op pr IIA 95 °C Ga/Gb

⊕_{Ex} II 1/2G Ex op is/ op pr IIB 38 °C Ga/Gb

Designation and number of the test certificate:

SEV 12 ATEX 0169X

CML 22 UKEX 9210X

2.4.2 Ex-Classification

 II 1/2G Ex op is/ op pr IIA 95 °C Ga/Gb and
 II 1/2G Ex op is/ op pr IIB 38 °C Ga/Gb according to
**EU type Examination certificate SEV 12 ATEX 0169X and
 CML 22 UKEX 9210X**

Introduction

The turbidity sensors METTLER TOLEDO InPro 81XX/82XX */*/*/* may be used in accordance with EN IEC 60079-0:2018 and EN60079-28:2015 in the categories 1 to 3 (Zones 0 to 2) and gas groups IIA and IIB only according to above (see Ex-classification) stated max. temperature specifications.

Special conditions X

1. The turbidity probes InPro 8100/** and InPro 8200/** may only be used with the transmitter Trb 8300 (phased out) or M800 1-channel. These transmitters, not situated in the explosive atmosphere but containing the light source (LED) are associated apparatuses which have a specified ambient temperature range from -10 °C to $+50\text{ °C}$.
2. Outside the tank/vessel in which the probe is located, the ambient temperature range is: -20 °C ... $+60\text{ °C}$.
3. The maximum permissible temperature (a.k.a. limiting temperature) of the process medium in which the probe is located (inside the tank/vessel) determines the probe maximum surface temperature. These temperatures should be considered as follows:

Gas group	Medium limiting temperature → probe maximum surface temperature
IIA	95 °C
IIB	38 °C

4. The metal body of the turbidity probes (resp. safety shrink-wrapped nozzle, fittings InFit™ 76*_* or InTrac™ 7**_*) must be conductively connected to the equipotential bonding system of the installation.
5. The metal body of the turbidity probes (resp. safety shrink-wrapped nozzle, fittings InFit 76*_* or InTrac 7**_*) must be included in the periodic pressure test of the installation.

6. The conduit around the optical fibre provides an additional “op pr” EPL Gb type of protection of the optical fibre on a length of 1.0 m for probe InPro 8100/***(plastic conduit) and on a length of 52 cm InPro 8200/***(metallic armouring).

3. Installation and start up

3.1 Handling

The fiber optic cables have a core of Hard Clad Silica (HCS) fibers. When laying fiber optic cables, the bend radius must always exceed at least 10 cm (4"). Fiberoptic cables may never be shortened independently. The transparent plastic cap on the sensor is to be removed before use. Any formation of deposits of solid particles on the surface is to be avoided. Gas bubbles in the process medium are also able to produce a signal and therefore must be avoided.

3.2 Installation in vessels/reactors

The turbidity sensors are to be installed at a distance of at least 10 cm (4") from the walls of the vessel, agitator and other internals, in order to prevent disturbances through reflections. When employing special METTLER TOLEDO retractable housings InTrac with modified immersion tubes, the installation guideline valid for each individual housing is to be observed. This step likewise, is to avoid measurement disturbances arising from irrelevant light reflections (**see chapter 7.2: “Accessories/Housings”**).

3.3 Installation in pipes

The turbidity sensors are to be installed such that the distance from the pipe walls is at least 10 cm (4"). If this is not feasible in pipes of small diameter, then the sensor can be installed sloped, at an angle of 15...45° or, alternatively, in a knee bend. The pipe must always be completely full of liquid.

Connecting InPro 8000 Series Sensor to M800 Process 1-Channel Transmitter


Fig. 1: Special slotted gaskets for cable gland fitting with fiber optic.

Fig. 2: Gaskets mounted on fiber optic cable.

Two special slotted gaskets for cable glands (M25 x 1.5) are included in M800 Process 1-Channel transmitter delivery (Fig 1). Correct installation of these two gaskets is important to ensure proper transmitter ingress protection and to ensure the fiber optic cables are fixed properly to the transmitter.

The gaskets are pre-installed to the fiber optic cables, as Fig 2.

Note: Any vibration/movement of the fiber optic cables will affect the stability of the measurement.

Fig. 3: Gaskets are fit through M25 x 1.5 cable glands.

Fig. 4: SMA-connection to M800 transmitter "Receiver/Emitter".

The gaskets are fit through M25 x 1.5 cable glands. (Fig.3). Fiber optic cables are connected to "Receiver" and "Emitter" on M800 transmitters (Fig 4). Tighten the cable glands' screws by hand.

Note: Blue and orange fiber optic cables are identical. Which means commissioning engineer can connects the two fiber optic cables without considering the fiber optic cable color.

Note: Do not use any tools for tightening SMA to transmitter. Using tools could cause irreversible damage to the threads. Hand tightening is sufficient.

Note: Fiber optic cables are fragile. Do not bend the fiber optic cables.

Refer to M800 transmitter instruction manual for transmitter wall mounted/pipe mounted installation.

3.4 Fiber optic cable extensions

Fiber optic cables are available according to the list in **chapter 7.2 "Accessories"**. The duplex cable extensions are supplied complete with four fitted SMA plugs and two couplings. The two SMA plugs of the sensor fixed cable are to be screwed hand-tight onto the couplings. Two SMA plugs of the extension cable can then also be screwed by hand onto the couplings in any order. To protect the cable connection from dust and splashwater, it is recommended to use the METTLER TOLEDO Connection Box IP 65 (Nema 4X) (order no. 52 800 241).

If extension cables are installed subsequently at a later date, the measuring system has to be newly calibrated (see instruction manual for the METTLER TOLEDO Transmitters M800 1-channel, **chapter: "Calibration"**).

3.5 Calibration in open laboratory vessels

The sensor must be calibrated during first start-up. Multi-point/Process/in-Situ calibration are available on M800 1-channel transmitter. Recalibration is required when installing a new sensor/extension cables are installed. Refer to M800 transmitter instruction manual for sensor calibration procedures.

If the turbidity measuring system is calibrated off-line, it is recommended to use the METTLER TOLEDO calibration tool, CaliCap™ (order no. 52 800 210), for fitting on the sensor shaft. If CaliCap is not used, a minimum distance of 10 cm (4") to the vessel walls has to be observed as well. During calibration, the sample should be stirred in order to prevent the sedimentation of solid particles. Avoid direct sun light/exposing the sensor end (nose) to direct light.

4. Maintenance and troubleshooting

4.1 Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for the period of one year from date of delivery ex works. If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, carriage paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiry of the period of warranty, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

4.2 Maintenance

Depending on plant Standard Operating Procedure (SOP), it is recommended to verify the measurement system with a CaliCap (dry check) before each new production batch. When deviation of more than $\pm 10\%$ is observed, it is recommended to recalibrate the sensor. For continuous production, it is recommended to verify the measurement system with a CaliCap every three months.

Dirty or contaminated sensors can deliver incorrect measurement values. If fouling is presumed, the sensor is to be removed and cleaned. Suitable cleaning solutions are mild detergents or strongly diluted acids ($< 0.5\%$ by wt.) such as hydrochloric acid.

When handling acids, precautionary measures are to be taken at all times.

5. Removal from operation, storage, disposal

5.1 Removal from operation

The sensor is only conditionally subject to aging. When employed correctly as intended and appropriately maintained and serviced, the lifetime of the sensor can extend to several years.

Before removing the sensor from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

After demounting, the sensor should first be flushed with clear water.

If the sensor incurs a defect, it can as a rule not be repaired and must be disposed under observance of prevailing regulations.

5.2 Storage

If the sensor is not needed for use, it can be stored dry.


5.3 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.


6. Dimensional drawings


InPro 8200/S(H)

InPro 8100/S

InPro 8050

All dimensions in millimeters

7. Ordering information

7.1 Sensors

Specifications duplex fiber optic cable

Cladding material	HCS (Hard Clad Silica)
Core diameter	600 ±10 µm
Attenuation @ 850 nm	≤8 dB/km
Long-term bend radius	≥94 mm
Outer jacket	PVC

Ordering information

Designation code	InPro 8100/S/xxx	S = shaft material DIN 1.4435 (316L)
	InPro 8200/S(H)/...../xxx	H = Hastelloy C22 = sealing for sapphire window xxx = insertion length in mm

Designation	Insertion lengths	Body material	Order no.
InPro 8050	120 mm (4.7")	PSU	52 800 209
InPro 8100/S/120	120 mm (4.7")	DIN 1.4435 (316L)	52 800 205
InPro 8100/S/205	205 mm (8.1")	DIN 1.4435 (316L)	52 800 206
InPro 8100/S/297	297 mm (11.7")	DIN 1.4435 (316L)	52 800 207
InPro 8100/S/407	407 mm (16.0")	DIN 1.4435 (316L)	52 800 208
InPro 8200/S/Epoxy/120	120 mm (4.7")	DIN 1.4435 (316L)	52 800 216
InPro 8200/S/Epoxy/205	205 mm (8.1")	DIN 1.4435 (316L)	52 800 217
InPro 8200/S/Epoxy/297	297 mm (11.7")	DIN 1.4435 (316L)	52 800 218
InPro 8200/S/Epoxy/407	407 mm (16.0")	DIN 1.4435 (316L)	52 800 219
InPro 8200/H/Epoxy/120	120 mm (4.7")	HA-C22 (Hastelloy)	52 800 220
InPro 8200/H/Epoxy/205	205 mm (8.1")	HA-C22 (Hastelloy)	52 800 221
InPro 8200/H/Epoxy/297	297 mm (11.7")	HA-C22 (Hastelloy)	52 800 222
InPro 8200/H/Epoxy/407	407 mm (16.0")	HA-C22 (Hastelloy)	52 800 223
InPro 8200/S/Kalrez-FDA/120	120 mm (4.7")	DIN 1.4435 (316L)	52 800 224
InPro 8200/S/Kalrez-FDA/205	205 mm (8.1")	DIN 1.4435 (316L)	52 800 225
InPro 8200/S/Kalrez-FDA/297	297 mm (11.7")	DIN 1.4435 (316L)	52 800 226
InPro 8200/S/Kalrez-FDA/407	407 mm (16.0")	DIN 1.4435 (316L)	52 800 227

7.2 Accessories

O-rings

Designation	Size	Order no.
Viton-FDA	10.77 x 2.62 mm	59 905 680
Silicon-FDA	10.77 x 2.62 mm	59 905 678
Kalrez-FDA	10.78 x 2.62 mm	59 905 789
EPDM-FDA	10.77 x 2.62 mm	59 905 768

Fiber optic cable extensions including two couplings, for all sensors:

Designation	Lengths	Order no.
Fiber optic extension 3 m	3 m (10 ft)	52 800 228
Fiber optic extension 5 m	5 m (15 ft)	52 800 229
Fiber optic extension 6 m	6 m (18 ft)	52 800 230
Fiber optic extension 10 m	10 m (30 ft)	52 800 231
Fiber optic extension 15 m	15 m (45 ft)	52 800 232
Fiber optic extension 20 m	20 m (60 ft)	52 800 233
Fiber optic extension 25 m	25 m (75 ft)	52 800 234
Fiber optic extension 30 m	30 m (100 ft)	52 800 235
Fiber optic extension 75 m	75 m (230 ft)	52 800 177
Fiber optic extension 100 m	100 m (330 ft)	52 800 154

Additional accessories, for all sensors:

Designation	Order no.
CaliCap (calibration tool)	52 800 210
Couplings (2 included in each fiber kit)	52 800 240
Coupling box IP 65/NEMA 4x	52 800 241
Swagelok adapter NPT 1/2" ¹⁾	52 800 242

¹⁾ Only to be used with InPro 8200/S/Epoxy Sensors for pressure range 0...60 bar (0...870 psi).

Transmitter

Designation	Order no.
M800, 1-channel, Process	30 026 633

Housings for InPro 8050

Type of housing installation	Designation	Wetted part material	Insertion length (mm)	Sensor length (mm)	Process connection	O-ring	Order no.
Direct	InFit 761/NS/0070/	PVDF	70	120	DN25	Viton-FDA	52 400 311
Inline	PVDF/D00/Vi9-						
Insertion	InFit 761/NS/0070/	PP	70	120	DN25	Viton-FDA	52 400 316
	PP_/D00/Vi9-						
Open	InDip®550/1000 PVC	PVC	1000	120	N/A	Viton-FDA	52 400 320
Basin	Pg13.5						
	InDip550/1500 PVC	PVC	1500	120	N/A	Viton-FDA	52 400 582
	Pg13.5						
	InDip550/2000 PVC	PVC	2000	120	N/A	Viton-FDA	52 400 588
	Pg13.5						
	InDip550/2500 PVC	PVC	2500	120	N/A	Viton-FDA	52 400 594
	Pg13.5						
	InDip550/3000 PVC	PVC	3000	120	N/A	Viton-FDA	52 400 600
	Pg13.5						

Housings for InPro 8100 and InPro 8200

Type of housing installation	Designation	Wetted part material	Insertion length (mm)	Sensor length (mm)	Process connection	O-ring	Order no.
Direct	InFit761/WS/0070/	DIN 1.4435	70	120	DN25	Viton	59 900 753
Inline	4435/D10/Vi2-						
Insertion	InFit761/NS/0070/	DIN 1.4435	70	120	DN25	Silicon-FDA	59 900 796
	4435/D00/Si9-						
	InFit761/NC/0070/	DIN 1.4435	70	120	DN25	Silicon-FDA	52 400 491
	4435/D00/Si9-						
	InFit761/NC/0033/	DIN 1.4435	33	120	Varivent®	Silicon-FDA	52 400 502
	4435/V02/Si-						
	InFit761/NC/0070/	DIN 1.4435	70	120	TC 1.5"	Silicon-FDA	52 400 494
	4435/T01/Si-						
	InFit761/NC/0070/	DIN 1.4435	70	120	TC 2"	Silicon-FDA	52 400 495
	4435/T02/Si-						
Retractable from the process	InTrac779M/070/	DIN 1.4404	70	205	DN25	Viton-FDA	52 403 236
	4404/D00/Vi/A00						
	InTrac779M/100/	DIN 1.4404	100	297	DN25	Viton-FDA	52 403 237
	4404/D00/Vi/A00						

Note: Please consult your local METTLER TOLEDO representative for more option.

8. Certificates

8.1 Material certificate

Each sensor with a stainless steel shaft made of DIN 1.4435 or 1.4404 (316L) is delivered with a set of 3.1B certificates (complying with EN10204.3/1.B). All wetted metal parts are identified with an engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part is polished in order to get a surface roughness lower than $0.4 \mu\text{m}$ ($16 \mu\text{in}$). This represents a roughness grade number of N5 (according to ISO1320:1992)

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