Puncture-Type Combination pH Electrodes LoT406-M6-DXK-S7/25

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fi	Käyttöohje	22
fr	Instructions d'utilisation	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

LoT 406 Sensor Series
59 906 021 B

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
Operating Instruction for puncture-type combination pH electrodes

LoT406-M6-DXK-S7/25

Contents

1	Introduction	15
2	Important notes.....	15
3	Connecting the electrode	15
4	Operation	15
4.1.	Getting ready	15
4.2.	Calibrating the electrode.....	15
4.3.	Measuring.....	16
5	Maintenance.....	16
5.1.	Cleaning	16
5.2.	Storing the electrode	16
5.3.	Environmental protection.....	16
6	Trouble-shooting	17
7	Electrode specification	17
8	Factors that shorten electrode life	17
9	Literature.....	17

Description of the electrode

Argenthal, Xerolyt and Friscoylt are trademarks of the METTLER TOLEDO Group.

1 Introduction

METTLER TOLEDO puncture-type pH the electrodes with Xerolyt polymer electrolyte are used for direct measurement at pH in foods, for example fresh or cured meat, cheese and fish. Before using the electrode, read instructions carefully.

2 Important notes

On receipt, check that the pH electrode is undamaged and in good condition. For further information, ordering and service (advice), please refer to your supplier.

Safety


- The electrodes tip is made of glass. If the electrode breaks, there is a risk of injury and of glass splinters in the measured material.
- Wear gloves and goggles when handling solutions, e.g. when attaching the watering cap filled with electrolyte, and when calibrating with pH buffer solutions.

3 Connecting the electrode


Use an AS7 cable to connect the puncture pH electrode directly to a pH meter or pH transmitter (Fig. 4).

4 Operation

4.1. Getting ready


4.2. Calibrating the electrode

Please read the instructions for the pH meter or transmitter. We recommend calibrating at the temperature of the measured material. A pH/temperature table is printed on the buffer label. Automatic temperature compensation is possible by using a separate temperature sensor. If highly accurate measurement is required, insert the electrode the material 3–4 minutes before starting to take readings.

4.3. Measuring


Follow the instructions for the pH/mV meter or transmitter.

5 Maintenance


5.1. Cleaning

We recommend regular cleaning

Slight contamination

Fat residues

5.2. Storing the electrode

Friscolyt™ or KCl 3 mol/l

Important:

- Do not keep pH electrodes in demineralized or distilled water, as this will shorten their life.
- Do not store pH electrodes dry. If the electrode has been stored dry, regenerate it overnight in the watering cap with KCl 3 mol/l or Friscolyt.

5.3. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

6 Troubleshooting

Problem	Cause	Remedy
Air bubbles inside electrode tip	Bubbles formed while stored horizontally	<ul style="list-style-type: none"> – With tip pointing down, shake the electrode with short flicking movements (like a clinical thermometer) – Store electrode in upright position
Air bubbles on aperture diaphragm;	Temperature fluctuations	Place electrode in tap water at 40 °C until bubbles disappear
Slope not adjustable at pH 4	Temperature setting on pH meter	Make sure pH meter is set to the buffer temperature
	Buffer solution is dirty or old	Fresh buffer solution
	Electrode oily or dirty	If heavily fouled, clean electrode with ethanol/ether 1:1 or hand-hot (40 °C) soap solution
The pH meter shows unstable readings	Electrode oily or dirty	If heavily fouled, clean electrode with ethanol/ether 1:1 or hand-hot (40 °C) soap solution
	Problems with cable	Check that cable is properly connected, or change the cable
Unusual measured values	pH meter, electrode, cable	<ul style="list-style-type: none"> – Check contacts at electrode/cable and cable/pH meter – Check cable or use a different one – Change electrode and recalibrate your pH meter

7 Electrode specification

Designation	LoT406-M6-DXK-S7/25
Order No.	10 406 3123
pH measurement range	2 ... 11
Temperature range	0 ... 80 °C
Reference system	Argenthal
Reference electrolyte	Xerolyt
Diaphragm type	Aperture
Plug-in connector type	S7
Penetration depth	25 mm
Diameter	6 mm
Shaft	PBT (FDA approved)

PBT = Polybutyleneterephthalate

8 Factors that shorten electrode life

- Dry storage
- pH < 2
- Prolonged immersion in organic solvents
- Storage in distilled or demineralized water
- Temperature fluctuations

9 Literature

- Practice and theory of the pH measurement
20 606 1531-E
- Application pH Measurement in meat
52 000 455
- Application pH Measurement in cheese making
52 000 459
- Application pH Measurement in raw meat processing
52 000 467

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All rights reserved