

Reference Manual

METTLER TOLEDO

METTLER TOLEDO


Standard Interface Command Set

MT-SICS 0 version 2.3x

MT-SICS 1 version 2.2x

MT-SICS 2 for Excellence Balances

MT-SICS 3 for Excellence Balances



Contents

1	Introduction	4
2	Overview of all commands	6
3	Data interchange with the balance	9
3.1	Command formats	9
3.2	Response formats	9
3.2.1	Format of the response with weight value	10
3.2.2	Format of the response without weight value	11
3.2.3	Error messages	12
3.2.4	Tips for the programmer	12
4	Commands and responses	13
4.1	Commands and responses MT-SICS level 0	13
4.2	Commands and responses MT-SICS level 1	23
4.3	Commands and responses MT-SICS level 2 for Excellence balances	30
4.4	Commands and responses MT-SICS level 3 for Excellence balances	84
4.5	Additional commands – FastHost for Excellence balances	94
5	Special features	101
6	An example	102
7	What if...?	104

1 Introduction


In weight measurements the demands on the readability and maximum capacity of balances and scales range from less than one microgram up to several hundred tonnes. To meet these and other requirements, METTLER TOLEDO offers an extensive range of balances and scales. Many of the balances and scales used have to be capable of integration in a complex computer or data acquisition system. To enable you to integrate balances in your system in a simple manner and utilize their capabilities to the full, most balance functions are also available as appropriate commands via the data interface.

Standardization of the commands

All new METTLER TOLEDO balances launched on the market support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS), which is divided into 4 levels, depending on the functionality of the balance:

- MT-SICS level 0 Command set for the simplest balance, e.g. weighing cell.
- MT-SICS level 1 Extension of the command set for standard balances, i.e. balances without integrated applications.
- MT-SICS level 2 Extension of the command set by the commands specific for a balance family, e.g. MT-SICS level 2 for the Excellence balance line.
- MT-SICS level 3 Application-specific commands as an extension of the command set, e.g. MT-SICS level 3 for piece counting or FreeWeigh application with Remote -x software.

A particular distinguishing feature of this concept is that the commands combined in MT-SICS level 0 and 1 are identical for all balances. Both the simplest weighing balance and a fully expanded weighing work station recognize the commands of MT-SICS level 0 and 1. Investigations of various applications have shown that the vast majority of all system solutions can be handled with the commands of MT-SICS level 0 and 1. This means for you: if you restrict yourself to the commands of MT-SICS level 0 and 1, you can expand your system with additional balances from METTLER TOLEDO without having to change your application programs.



The commands written in this booklet contain MT-SICS commands for Excellence balances. As the balances can differ based on model and software version, not all the MT-SICS level 2 and 3 commands are usable on every model. We therefore recommend to use the "iO" command to get an overview of all commands that are supported by a particular balance. Commands that are listed with the "iO" command, but not described in this booklet, are not available for the user.

What do the commands of MT-SICS level 0 and 1 offer?

You can use the commands of MT-SICS level 0 and 1 to perform the following operations via the interface:

- Request weighing results,
- Tare the balance and preset the tare weight,
- Zero the balance,
- Identify MT-SICS implementation,
- Identify the balance,
- Reset the balance,
- Control the display,
- Control the keys for operation of the balance.

The commands of MT-SICS level 2 and 3

You can naturally use the data interface to exploit all functions available with your current balance or application. These additional functions are collected in the commands of MT-SICS level 2 and 3.

When creating your software application, please note that whereas the commands of MT-SICS level 2 have been specially tailored to your balance family, the commands of MT-SICS level 3 apply to your weighing application and can not be supported by every METTLER TOLEDO balance. In the enclosure with this manual, you can see what commands are supported by your balance and application.

Additional documentation on data interface

Settings of the interface such as baud rate, number of data bits, parity, handshake protocols and connector pin assignment are described in the operating instructions of the peripheral instrument or cable in question.

Version number of the MT-SICS

Each level of the MT-SICS has its own version number which can be requested with the command I1 from level 0.

This manual describes

- MT-SICS level 0, version 2.3x
- MT-SICS level 1, version 2.2x
- MT-SICS level 2 for Excellence balances
- MT-SICS level 3 for Excellence balances

You can use the command **I1** via the interface to request the MT-SICS level and MT-SICS versions implemented on your balance.

Please make sure that the versions implemented on your balance agree with those listed above.

2 Overview of all commands

Commands MT-SICS level 0		Page
I0	Inquiry of all implemented MT-SICS commands	13
I1	Inquiry of MT-SICS level and MT-SICS versions	14
I2	Inquiry of balance data	15
I3	Inquiry of balance SW version and type definition number	15
I4	Inquiry of serial number	16
I5	SW-Identification number	16
S	Send stable weight value	17
SI	Send weight value immediately	18
SIR	Send weight value immediately and repeat	19
Z	Zero	20
ZI	Zero immediately	21
@	Reset	22

Commands MT-SICS level 1		Page
D	Balance display	23
DW	Weight display (Display show Weight)	23
K	Key control	24
SR	Send weight value on weight change (Send and Repeat)	26
T	Tare	27
TA	Inquiry/setting of tare weight value	28
TAC	Clear tare value	29
TI	Tare Immediately	29

Commands MT-SICS level 2 for Excellence balances		Page
C0	Inquiry/setting of calibration setting	30
C1	Initiate calibration according to current setting	32
C2	Initiate calibration with external weight	33
C3	Initiate calibration with internal weight	34
COM	Inquiry/Setting the communication parameters of the serial interface	35
COPT	Command to configure interface options	37
DAT	Date	41
I10	Balance ID – Inquiry of balance identification	42
I11	Balance type	42
I14	Inquiry of balance info	43
MO1	Inquiry/setting of weighing mode	44
MO2	Inquiry/setting of environment	44
MO3	Inquiry/setting of AutoZero	45
MO4	Inquiry/setting of SmartSens functions	45

M05	Inquiry of user list/method	46
M06	Inquiry/setting of current user/method number	46
M08	Inquiry/setting of display brightness	47
M09	Inquiry/setting of display contrast	47
M11	Inquiry/setting of beeper volume	48
M12	Creating beeper tone	48
M13	Inquiry/setting of Touch function	49
M14	Inquiry of available languages	49
M15	Inquiry/setting of language	50
M16	Inquiry/setting of standby mode	51
M17	Inquiry/setting of ProFACT time criteria	52
M18	Inquiry/setting of ProFACT/FACT temperature criterion (D temp.)	53
M19	Inquiry/setting of adjustment weight	54
M20	Inquiry/setting of test weight	55
M21	Inquiry/setting of unit	56
M22	Inquiry/setting of custom unit definitions	57
M23	Inquiry/setting of readability, 1d/xd	58
M24	Inquiry/setting of print key function	58
M25	Inquiry/setting of application selection	59
M26	Inquiry/setting of current application	59
M27	Inquiry of adjustment history	60
M28	Inquiry of temperature probe	61
M29	Inquiry/setting of value release	61
M36	Inquiry/setting of LevelControl	62
M37	Inquiry/setting of draft shield door opening	63
M38	Execute reset	64
P100	Print out text on the printer	65
P101	Send stable weight value to printer channel	65
P102	Send weight value to printer channel immediately	66
P120	Reset SmartTrac according to application	66
P121	Set SmartTrac as +/- display	67
P122	Activate individual pointers of SmartTrac	68
P123	Activate SmartTrac by segments	69
P124	Switch off SmartTrac	69
PWR	Power on/off	69
SIS	Inquiry of current net weight values	70
SIRU	Send weight value with currently displayed unit immediately and repeat	72
SIU	Send weight value with currently displayed unit immediately	72
SNR	Send stable weight value and repeat after each deflection	73
SNRU	Send stable weight value with currently displayed unit and repeat after each deflection	74

MT-SICS Basic Information

SRU	Send stable weight value with currently displayed unit after deflection	75
ST	Send stable weight value after pressing F (transfer) key	76
SU	Send stable weight value with currently displayed unit	77
TIM	Time	78
TST0	Inquiry/setting of the test function	79
TST1	Initiate test function in the current setting	80
TST2	Initiate test function with external weight	81
TST3	Initiate test function with internal weight	82
UPD	Inquiry/setting of the update rate of the host interface	82
WS	Inquiry/setting of draft shield	83

Commands MT-SICS level 3 for Excellence balances

		Page
PW	Piece Counting: Inquiry/setting of the piece weight	84
A01	Percent Weighing: Inquiry/setting of reference in %	85
A06	Dynamic Weighing: Inquiry/setting of dynamic weighing filter	86
A07	Dynamic Weighing: Inquiry/setting of dynamic weighing AutoStart	87
A08	Dynamic Weighing: Inquiry/setting of dynamic weighing AutoTara	88
A09	Dynamic Weighing: Inquiry/setting of data acquisition	89
SM0	Dynamic Weighing: Reset of all SMx commands	89
SM1	Dynamic Weighing: Start immediately and transfer the result	90
SM2	Dynamic Weighing: Start after a minimum load is exceeded and transfer the result once	91
SM3	Dynamic Weighing: Start after a minimum load is exceeded, transfer the result and repeat	92
SM4	Dynamic Weighing: Inquiry/setting of time interval	93

Additional commands – FastHost for Excellence balances

		Page
B00	FastHost-list of commands	94
B01	Inquiry of FastHost individual value	94
B02	FastHost start/stop continuous mode	95
B03	FastHost Inquiry of time interval of value counter	95
B04	FastHost Inquiry and setting of format specification	96
B05	FastHost Inquiry and setting of stability criteria	98
B06	FastHost Inquiry and setting sampling reduction	99
B07	FastHost Inquiry and setting switch-on mode	99
B08	Set to zero with FastHost stability criteria	100

3 Data interchange with the balance

Each command received by the balance via the data interface is acknowledged by a response of the balance to the transmitter.

Commands and responses are data strings with a fixed format, and will be described in detail in chapter 3.

3.1 Command formats

Commands sent to the balance comprise one or more characters of the ASCII character set. Here, the following must be noted:

- Enter commands only in uppercase.
- The possible parameters of the command must be separated from one another and from the command name by a space (ASCII 32 dec., in this description represented as \square).
- The possible input for "text" is a sequence of characters of the 8-bit ASCII character set from 32 dec to 255 dec.
- Each command must be closed by $C_{R}L_{F}$ (ASCII 13 dec., 10 dec.).

The characters $C_{R}L_{F}$, which can be inputted using the Enter or Return key of most entry keypads, are not listed in this description, but it is essential they be included for communication with the balance.

Example

Command to balance which writes Hallo into the balance display:

$D\square$ "Hallo" The command terminator $C_{R}L_{F}$ is not shown.

Comment

The quotation marks " " must be inserted in the entry.


3.2 Response formats

All responses sent by the balance to the transmitter to acknowledge the received command have one of the following formats:

- Response with weight value
- Response without weight value
- Error message

3.2.1 Format of the response with weight value

A general description of the response with weight value is the following.



ID	Response identification.
▯	Space (ASCII 32 dec.).
Status	Status of the balance, see description of the commands and responses.
WeightValue	Weighing result; shown as number with 10 digits, incl. decimal point and sign – directly in front of the first digit if value negative. The weight value appears right-aligned. Preceding zeros are not shown with the exception of the zero to the left of the decimal point. With METTLER TOLEDO DeltaRange balances, outside the fine range the last decimal place is shown as a space.
Unit	Weight unit actually set under unit 1.
C_R	Carriage Return (ASCII 13 dec.).
L_F	Line Feed (ASCII 10 dec.).

Comment

C_RL_F will not be shown in this description.

Examples

Response with stable weight value of 0.256 g:

S▯S▯S▯S▯S▯S▯S▯0.256▯g

Response with stable weight value outside the fine range:

S▯S▯S▯S▯S▯S▯4875.2▯g

3.2.2 Format of the response without weight value

A general description of the response without weight value is the following.



1 - 4 characters 1 character

ID	Response identification.
┘	Space (ASCII, 32 dec.).
Status	Status of the balance, see description of the commands and responses.
Parameters	Command-dependent response code.
C_R	Carriage Return (ASCII 13 dec.).
L_F	Line Feed (ASCII 10 dec.).

Comment

C_RL_F will not be shown in this description.

Example

Response to **D┘"Halla"** when Hallo appears unbridged in the display: **D┘A**.

3.2.3 Error messages

ID	C _R	L _F
----	----------------	----------------

There are three different error messages. The identification always comprises two characters.

ID	Error identification Possible error messages are
ES	Syntax error The balance has not recognized the received command.
ET	Transmission error The balance has received a "faulty" command, e.g. owing to a parity error or interface break.
EL	Logical error The balance can not execute the received command.

C_R Carriage Return (ASCII 13 dec.).

L_F Line Feed (ASCII 10 dec.).

Comment

C_RL_F will not be shown in this description.

3.2.4 Tips for the programmer

Command and response

You can improve the dependability of your application software by having your program evaluate the response of the balance to a command. The response is the acknowledgement that the balance has received the command.

Reset

To be able to start from a determined state, when establishing the communication between balance and system, you should send a reset command to the balance. When the balance or system is switched on or off, faulty characters can be received or sent.

Quotation marks " "

Quotation marks included in the command must always be entered.
Set/get configuration for options.

4 Commands and responses

4.1 Commands and responses MT-SICS level 0

IO Inquiry of all implemented MT-SICS commands

Command	IO	Send list of all implemented MT-SICS commands.
Response	IO␣B␣x1␣"1.Command"	x1 = Number of the MT-SICS level where the 1. Command belongs to.
	IO␣B␣x1␣"2.Command"	: 2nd (next) command implemented.
	:	:
	IO␣A␣x1␣"last Command"	Last command implemented.
	IO␣I	The list cannot be sent at present as another operation is taking place.

Example

Command	IO	Send list of commands.
Response	IO␣B␣0␣"IO"	Level 0 command "IO" implemented.
	IO␣B␣0␣"I1"	Level 0 command "I1" implemented.
	:	:
	:	:
	:	:
	IO␣B␣0␣"S"	Level 0 command "S" implemented.
	:	:
	:	:
	IO␣B␣0␣"Z"	Level 0 command "Z" implemented.
	IO␣B␣0␣"@"	Level 0 command "@" (reset) implemented.
	IO␣B␣1␣"D"	Level 1 command "D" implemented.
	IO␣B␣1␣"DW"	Level 1 command "DW" implemented.
	:	:
	IO␣A␣3␣UPD	(last command).

Comments

- The IO command lists all commands implemented in the present software. Thus, there is no need of the supplement sheet delivered with the previous versions of this manual.
- All level 0 commands are listed in alphabetical order before all commands of level 1 etc. This order corresponds to the order how the commands are described in this manual.

I1 Inquiry of MT-SICS level and MT-SICS versions

Command	I1	Inquiry of MT-SICS level and MT-SICS versions.
Response	I1└┬┘└┬┘"x1"└┬┘"x2"└┬┘"x3"└┬┘"x4"└┬┘"x5"	
	x1 = 0	Balance with MT-SICS level 0 (simplest balance).
	x1 = 01	Balance with MT-SICS level 0 and 1 (standard balance).
	x1 = 012	Balance with MT-SICS level 0, 1 and 2 (standard balance with extensions).
	x1 = 03	Balance with MT-SICS level 0 and 3 (simplest balance with a special application).
	x1 = 013	Balance with MT-SICS level 0, 1 and 3 (standard balance with a special application).
	x1 = 0123	Balance with MT-SICS level 0, 1, 2, and 3 (standard balance with extensions and a special application).
	x1 = 3	Application device with MT-SICS level 3 (not necessarily a balance).
	X2	Version of the implemented MT-SICS 0 commands.
	X3	Version of the implemented MT-SICS 1 commands.
	X4	Version of the implemented MT-SICS 2 commands.
	X5	Version of the implemented MT-SICS 3 commands.
	I1└┬┘	Command understood, not executable at present.

Example

Command	I1	Inquiry of MT-SICS level and versions.
Response	I1└┬┘└┬┘"01"└┬┘"2.30"└┬┘"2.22"└┬┘"└┬┘"	
		01 Level 0/1 implemented.
		2.30 Level 0, version V2.30
		2.22 Level 1, version V2.22

Comments

- In the case of the MT-SICS level, only fully implemented levels are listed. In other words, if it is not possible to implement all commands from a certain level, the level is not specified.
- In the case of the MT-SICS version, all levels are specified even those only partially implemented.

I2 Inquiry of balance data

Command	I2	Inquiry of balance data.
Response	I2LA "text"	Balance data as "text".
	I2/I	Command understood, not executable at present.

Example

Command	I2	Inquiry of balance type.
Possible responses	I2LA "PR5002DRUR-StandardL5100.90Lg"	
	I2LA "X205BDU ExcellenceL220.0090Lg"	

Comments

- With DeltaRange balances, the last decimal place is available only in the fine range.
- The number of characters of "text" depends on the balance type.

I3 Inquiry of balance SW version and type definition number

Command	I3	Inquiry of balance SW version and type definition number.
Responses	3LA "text"	Balance SW version and type definition number as "text".
	I3/I	Command understood, not executable at present.

Examples

Command	I3	Inquiry of SW version number(s) and type definition number.
Response	I3LA "1.05L1.1.1.17.7"	1.05 Software version number 1.1.1.17.7 Type definition number

Comments

The first number (digits prior to the first space in the text string) is the SW version number. The second SW version number is optional, and depends on the balance type. The last number (following the last space) is the type definition number for service purposes.

I4 Inquiry of serial number

Command **I4** Inquiry of serial number.
 Responses **I4LAL"text"** Serial number as "text".
I4LI Command understood, not executable at present.

Example

Command **I4** Inquiry of serial number.
 Response **I4LAL"0123456789"**

Comments

- The serial number agrees with that on the model plate and is different for every balance.
- The serial number can be used, for example, as a device address in a network solution.
- The response to **I4** appears unsolicited after switching on and after the reset command (@).

I5 SW-Identification number

Command **I5** Inquiry of SW-Identification number.
 Responses **I5LAL"x"** SW-Identification number as Text.
 x: SW-Identification number.
I5LI Command understood, not executable at present.

Example

Command **I5** Inquiry of SW-Identification number.
 Response **I5LAL"12345678A"**
 SW-Identification number with index.

Comment

The SW-Identification number is unique for every Software.

S	Send stable weight value
----------	---------------------------------

Command	s	Send the current stable net weight value.
Response	SLSLWeightValueLUnit	Current stable weight value in unit actually set under unit 1.
	SLI	Command not executable (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).
	SL+	Balance in overload range.
	SL-	Balance in underload range.

Example

Command	s	Send a stable weight value.
Response	SLSLSLSLSLS100.00Lg	The current, stable weight value is 100.00 g.

Comments

- The duration of the timeout depends on the balance type.
- To send the stable weight value in actually displayed unit, see "sT" command in level 2
- The draft shield closes with this command, when the "Door function" is set on "Automatic". It opens after sending a stable weight.

SI	Send weight value immediately
-----------	--------------------------------------

Command	SI	Send the current net weight value, irrespective of balance stability.
Response	SUSUWeightValueUUnit	Stable weight value in unit actually set under unit 1.
	SUDUWeightValueUUnit	Nonstable (dynamic) weight value in unit actually set under unit 1.
	SI	Command not executable (balance is currently executing another command, e.g. taring).
	SU+	Balance in overload range.
	SU-	Balance in underload range.

Example

Command	SI	Send current weight value.
Response	SUDUUUUUU129.07Ug	The current weight value is unstable (dynamic) and is 129.07 g.

Comments

- The response to the command **SI** is the last internal weight value (stable or dynamic) before receipt of the command **SI**.
- To send weight value immediately in actually displayed unit, see "**SIU**" command in level 2.

SIR	Send weight value immediately and repeat
------------	---

Command	SIR	Send the net weight values repeatedly, irrespective of balance stability.
Response	SUSUWeightValueUnit	Stable weight value in unit actually set under unit 1.
	SUDUWeightValueUnit	Nonstable (dynamic) weight value in unit actually set under unit 1.
	SI	Command not executable (balance is currently executing another command, e.g. taring).
	SU+	Balance in overload range.
	SU-	Balance in underload range.

Example

Command	SIR	Send current weight values at intervals.
Response	SUDUUUUUU129.07ug	
	SUDUUUUUU129.08ug	
	SUSUUUUUU129.09ug	
	SUSUUUUUU129.09ug	
	SUDUUUUUU114.87ug	
	...	The balance sends stable or nonstable weight values at intervals.

Comments

- **SIR** is overwritten by the commands **S**, **SI**, **SR**, **@** and hardware break and hence cancelled.
- The number of weight values per second depends on the balance type.
- To send weight value in actually displayed unit, see "**SIRU**" command in level 2

Z	Zero
----------	-------------

Command	z	Zero the balance.
Response	zL A	The following then holds: gross = net + tare = 0. Zero setting performed, i.e. stability criterion and zero setting range complied with.
	zL I	Zero setting not performed (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).
	zL +	Upper limit of zero setting range exceeded.
	zL -	Lower limit of zero setting range exceeded.

Example

Command	z	Zero.
Response	zL A	Zero setting performed.

Comments

- The tare memory is cleared during zero setting.
- The zero point determined during switching on is not influenced by this command, i.e. the measurement ranges remain unchanged.
- The duration of the timeout depends on the balance type.
- The draft shield closes with this command, when the "Door function" is set on "Automatic". It opens after sending a stable weight.

ZI	Zero immediately
-----------	-------------------------

Command	ZI	Zero the balance immediately regardless the stability of balance.
Response	ZI\perpD	Re-zero performed under non-stable (dynamic) conditions.
	ZI\perpS	Re-zero performed under stable conditions.
	ZI\perpI	Zero setting not performed (balance is currently executing another command, e.g. taring).
	ZI\perp+	Upper limit of zero setting range exceeded.
	ZI\perp-	Lower limit of zero setting range exceeded.

Example 1

Command	ZI	Zero immediately.
Response	ZI\perpS	Zero setting performed, weight value was stable.

Example 2

Command	ZI	Zero immediately.
Response	ZI\perpD	Zero setting performed, weight value was dynamic (nonstable).

Comments

- The tare memory is cleared after zero setting.
- The zero point determined during switching on is not influenced by this command, i.e. the measurement ranges remain unchanged.

@	Reset
----------	--------------

Command @ Resets the balance to the condition found after switching on, but without a zero setting being performed.

Response I4LAL"text" Serial number of the balance, the balance is ready for operation.

Example

Command @

Response I4LAL"1114350697"

Balance is reset, its serial number is 1114350697.

Comments

- All commands awaiting responses are cancelled.
- Key control is set to the default setting KL1.
- The tare memory is reset to zero.
- The "reset" command is always executed.
- If the balance is on standby, it is switched on.

4.2 Commands and responses MT-SICS level 1

D Balance display

Write into balance display

Command	D "text"	Write text into balance display.
Response	D A	Text appears unabridged left-aligned in the balance display marked by a symbol, e.g. *.
	D I	Command not executable.
	D L	Command understood, parameter wrong or balance with no display.

Example

Command	D "HALLO"	Write "HALLO" into the balance display.
Response	D A	The full text "HALLO" appears in the balance display.

Clear balance display

Command	D " "	Clear balance display.
Response	D A	Balance display cleared, marked by a symbol, e. g. *.
	D I	Clear balance display.

Comments

- A symbol in the display, e.g. * indicates that the balance is displaying an invalid weight value.
- The maximum number of characters of "text" visible in the display depends on the balance type.

DW Weight display (Display show Weight)

Command	DW	Switch main display to weight mode.
Response	DW A	Main display shows the current weight value.
	DW I	The command has been understood, but is not executable.

K	Key control
----------	--------------------

Commands	K_L1	When a key is pressed, execute the corresponding function, but do not send.
	K_L2	When a key is pressed, do not execute the corresponding function and send nothing.
	K_L3	When a key is pressed, do not execute the key function, but send the corresponding key code.
	K_L4	When a key is pressed, execute the corresponding function and send its function code. If the corresponding function can not be executed immediately, the function code K_LB_Uy for the start of the function and K_LA_Uy or K_LI_Uy for the end of the function are sent. This behavior applies to taring, zeroing, calibrating, testing, printing, etc. If a function may not be executed, the function code K_LI_Uy is sent.
Response	K_LA	Key control command understood and successfully executed.
	K_LI	Key control command understood but not executable at present, e.g. balance actually in menu or input mode.
	K_LL	Key control command understood, but command parameter wrong.

The key commands of the activate **K_L3** command are defined as follows:

Key commands of the Excellence XP balances

1		9
2		8
3	4	7
5	10	5

Response when pressed long.

When a code with a long press is sent, new key commands will not be accepted.

Example with an activated K_{L3} command:

$K_{LR}5$	Key 5 was pressed and held around 2 seconds.
$K_{LC}5$	Key 5 was released.

Response when K_{L4} is active

$K_{LA}y$	Function y was released by pressing the correspondent key and successfully executed.
$K_{LI}y$	Function y was released by pressing the correspondent key, but it could not be successfully executed, e.g. calibration was aborted by user.
$K_{LB}y$	Function y was released and started, the execution needs time to complete. These functions are marked with an asterix (*). After this response, either $K_{LA}y$ or $K_{LI}y$ follows. The balance functions are coded as follows: Calibration* $y = 0$ Tare* $y = 1$ Re-zero* $y = 2$ Data transfer to printing device* $y = 3$ Test* $y = 7$

Command	K_{L4}	When a key is pressed, execute the corresponding function and send the function code as an acknowledgement.
Responses	K_{LA}	Each time a key is pressed, immediate acknowledgement with the corresponding function code will be sent.
	$K_{LB}1$	The taring function has been started -> taring active.
	$K_{LA}1$	Taring completed successfully.
	$K_{LB}1$	The taring function has been started -> taring active.
	$K_{LI}1$	Taring not completed successfully, taring aborted.

Comments

- K_{L1} is the factory setting (default value).
- K_{L1} active after balance switched on and after the reset command
- Only one K command is active at any one time.
- A distinction must be made between key code K_{L3} and function code K_{L4} . The key code is specific to the balance type, the function code corresponds to the above table.

SR Send weight value on weight change (Send and Repeat)

Command	SR \downarrow PresetValue \downarrow Unit	Send the current stable weight value and then continuously after every weight change greater or equal to the preset value a nonstable (dynamic) value followed by the next stable value, range = 1d to max. load.
	SR	If no preset value is entered, the weight change must be at least 12.5 % of the last stable weight value, minimum = 30d.
Response	S \downarrow S \downarrow WeightValue \downarrow Unit	Current, stable weight value in unit actually set under until 1. Weight change.
	S \downarrow D \downarrow WeightValue \downarrow Unit	Dynamic weight value in unit actually set under until 1.
	S \downarrow S \downarrow WeightValue \downarrow Unit	Next stable weight value in unit actually set under until 1.
	S \downarrow I	Command not executable (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).
	S \downarrow L	Command understood, parameter wrong.
	S \downarrow +	Balance in overload range.
	S \downarrow -	Balance in underload range.

Example

Command	SR \downarrow 10.00 \downarrow g	Send the current stable weight value followed by every load change • 10 g.
Response	S \downarrow S \downarrow 100.00 \downarrow g	Balance stable.
	S \downarrow D \downarrow 115.23 \downarrow g	100.00 g loaded.
	S \downarrow S \downarrow 200.00 \downarrow g	Balance again stable.

Comments

- **SR** is overwritten by the commands **S**, **SI**, **SIR**, **@** and hardware break and hence cancelled.
- If, following a nonstable (dynamic) weight value, stability has not been reached within the timeout interval, the response "**S** \downarrow **I**" is sent and then a nonstable weight value. Timeout then starts again from the beginning.
- The preset value must be entered in unit actually set under until 1.

T	Tare
----------	-------------

Command	T	Tare, i.e. store the next stable weight value as a new tare weight value.
Response	TLSLWeightValueLUnit	Taring performed, i.e. stability criterion and taring range complied with. The tare weight value returned corresponds to the weight change on the balance in the unit actually set under unit 1 since the last zero setting.
	TLI	Taring not performed (balance is currently executing another command, e.g. zero setting, or timeout as stability was not reached).
	TL+	Upper limit of taring range exceeded.
	TL-	Lower limit of taring range exceeded.

Example

Command	T	The balance is tared and has a value of 100.00 g in the tare memory.
Response	TLSLLLLLLL100.00Lg	

Comments

- The tare memory is overwritten by the new tare weight value.
- The duration of the timeout depends on the balance type.
- The function of the combined tare and zero setting key corresponds to the zero setting (Z) command of the interface.
- Clearing tare memory: see command **TAC**.
- Unit 1 is the weight unit displayed after the balance has been switched on.
- The draft shield closes with this command, when the "Door function" is set on "Automatic". It opens after sending a stable weight.

TA	Inquiry/presetting of tare weight value
-----------	--

Inquiry of tare weight value

Command	TA	Inquiry of the tare weight value.
Response	TA┐┐┐TareWeightValue┐Unit	Current tare weight value in unit actually set under until 1.
	TA┐I	Current tare weight value can not be transferred at present as another operation is taking place.

Setting of tare preset value

Command	TA┐TarePresetValue┐Unit	Entry of a tare preset value in unit actually set under unit 1.
Response	TA┐┐┐WeightValue┐Unit	Entry accepted, returned value rounded to actual readability. The balance display shows the net value referred to the inputted tare value.
	TA┐I	Taring not performed (balance is currently executing another command, e.g. zero setting, or timeout as stability was not reached).
	TA┐L	Command understood, parameter wrong.

Example

Command	TA┐100.00┐g	Tare.
Response	TA┐┐┐┐┐┐┐┐┐100.00┐g	The balance has 100.00 g in the tare memory.

Comments

- The tare memory will be overwritten by the preset tare weight value.
- The inputted tare value will be automatically rounded by the balance to the current readability.
- The preset value must be entered in the unit actually set under unit 1.
- The taring range is specified to the balance type.

TAC Clear tare value

Command	TAC	Clear tare value.
Response	TACLA	Tare value cleared, 0 is in the tare memory.
	TACLI	Command not executable (balance is currently executing another command, e.g. zero setting, or timeout as stability was not reached).

TI Tare Immediately

Command	TI	Tare immediately, i.e. store the current weight value, which can be stable or non stable (dynamic), as tare weight value.
Response	TIUSWeightValueUnit	Taring performed, stable tare value. The new tare value corresponds to the weight change on the balance since the last zero setting.
	TIUDWeightValueUnit	Taring performed, non-stable (dynamic) tare value.
	TIUI	Taring not performed (balance is currently executing another command, e.g. zero setting).
	TIUL	The command is not executable, e.g. certified version of balance.
	TIU+	Upper limit of taring range exceeded.
	TIU-	Lower limit of taring range exceeded.

Example

Command	TI	Tare immediately.
Response	TIUDUUUUUU117.57ug	The tare memory holds a non-stable (dynamic) weight value.

Comments

- The tare memory will be overwritten by the new tare weight value.
- After a non-stable (dynamic) stored tare weight value, a stable weight value can be determined. However, the absolute value of the stable weight value determined in this manner is not accurate.
- The stored tare weight value is sent in the unit actually set under unit 1.
- The taring range is specified to the balance type.balances

4.3 Commands and responses MT-SICS level 2 for Excellence balances

C0 Inquiry/setting of calibration setting

Inquiry of calibration setting

Command **C0** Inquiry of the calibration setting.

Response **C0Lx1Lx2L"WeightValueLUnit"**
 Weight value and unit specify the value of the weight for an external calibration requested from the user via the display (see command **C2**). The unit corresponds to the factory setting of unit 1, e.g. gram (g) with standard balances or carat (ct) with carat balances respectively. With internal calibration, neither weight value nor unit appears.

x1 Calibration mode

x1=0 Mode = Manual
 The calibration can only be triggered manually. A change in the ambient conditions has no influence on the initiation of the calibration procedure.

x1=1 Mode = Auto, status display AutoCal or Cal not activated.
 The sensors built into the balance monitor the ambient conditions; however, the change is so small that a calibration is not necessary.

x1=2 Mode = Auto, status display "AutoCal" or "Cal" flashes.
 The sensors built into the balance have determined a considerable change in the ambient conditions. The balance requests a calibration or at least a test (see **TST** command).

x2 Calibration weight

x2=0 Internal weight (factory setting)
x2=1 External weight
 The current value of the external weight can be seen in the menu of the balance under "Calibration" (see Operating instructions).

C0LI The calibration status and the current setting of the calibration can not be transferred at present as another operation is taking place.

Example

Command **C0** Inquiry of status and setting of the calibration.

Response **C0LAL2L1L"LLLL100.000Lg"**
 Current setting of mode is "Auto". The ambient conditions of the balance have changed so much that the balance requests a calibration (x1=2) with the external weight (x2=1). For a calibration initiated with the command **C2**, a weight of 100.000 g is needed.

Setting the calibration setting

Command **C0Lx1Lx2** Set calibration setting.

x1 Calibration mode

x1=0 Mode = Manual
 A change in the ambient conditions has no influence on the initiation of the calibration procedure.

x1=1 Mode = Auto, the sensors built into the balance monitor the ambient conditions. When a considerable change in the ambient conditions is determined, the status display AutoCal or Cal will be activated; this means the balance will ask for calibration.

x2 Calibration weight

x2=0 Use internal weight (factory setting)

x2=1 Use external weight
 The current value of the external weight can be seen in the menu of the balance under "Calibration" (see Operating instructions).

Response **C0LA** Calibration setting set.

C0LL Calibration setting can not be set, e.g. parameter wrong or certified version of the balance.

C0LI Command not executable as the balance is, e.g. being tared.

Example

Command **C0L0L1** Set calibration setting to manual and external.

Response **C0LA** Calibration setting set.

Comment

Setting x1=1 and x2=0 corresponds to the menu setting "FACT" under "Calibration".

C1	Initiate calibration according to current setting
-----------	--

Command	C1	Start calibration in the current setting.
First response	C1LB	The calibration procedure has been started. Wait for second response (see Comment).
	C1LI	A calibration can not be performed at present as another operation is taking place. No second response follows.
	C1LL	Calibration operation not possible, e.g. with certified balance. No second response follows.
Further responses	C1L"text"	Weight request with external calibration.
	C1LA	Calibration has been completed successfully.
	C1LI	The calibration procedure was aborted as, e.g. stability not attained or wrong weights loaded.

Example

Command	C1	Start calibration.
Response	C1LB	Calibration operation started.
	C1L"LLLLLL0.00Lg"	Prompt to unload the balance.
	C1L"LLLL2000.00Lg"	Prompt to load calibration weight 2000.00 g.
	C1L"LLLLLL0.00Lg"	Prompt to unload the balance.
	C1LA	Calibration completed successfully.

Comment

Commands sent to the balance during the calibration operation are not processed and responded to in the appropriate manner until the calibration is at an end.

C2	Initiate calibration with external weight
-----------	--

Command	C2	Initiate external calibration. Inquiry of the weight used by means of the C0 command.
First response	C2LB	The calibration procedure has been started.
	C2LI	A calibration can not be performed at present as another operation is taking place. No second response follows.
	C2LJ	Calibration operation not possible, e.g. as a calibration with an external weight is not admissible (certified balance). No second response follows.
Further responses	C2L"text"	Prompt to unload or load the balance.
	C2LA	Calibration has been completed successfully.
	C2LI	The calibration procedure was aborted as, e.g. stability not attained or wrong weight loaded.

Example

Command	C2	Start calibration.
Response	C2LB	Calibration operation started.
	C2L"LLLLLLLL0.00Lg"	Prompt to unload the balance.
	C2L"LLLL2000.00Lg"	Prompt to load calibration weight 2000.00 g.
	C2L"LLLLLLLL0.00Lg"	Prompt to unload the balance.
	C2LA	Calibration completed successfully.

Comment

Commands sent to the balance during the calibration operation are not processed and responded to in the appropriate manner until the calibration is at an end.

C3	Initiate calibration with internal weight
-----------	--

Command	C3	Initiate internal calibration.
First response	C3LB	The calibration procedure has been started. Wait for second response.
	C3LI	A calibration can not be performed at present as another operation is taking place. No second response follows.
	C3LL	Calibration operation not possible, e.g. as internal weight missing. No second response follows.
Further responses	C3LA	Calibration has been completed successfully.
	C3LI	The calibration was aborted as, e.g. stability not attained or the procedure was aborted with the C key.

Example

Command	C3	Initiate internal calibration.
Response	C3LB	Calibration operation started.
	C3LA	Calibration completed successfully.

Comment

Commands sent to the balance during the calibration operation are not processed and responded to in the appropriate manner until the calibration is at an end.

COM Inquiry/setting the communication parameters of the serial interface

We recommend to use these parameters only for weighing platforms (without terminal).

Inquiry of the parameters of the serial interface

Command **COM** Inquiry of the current settings for all currently installed interfaces.

Response **COMLAL0LBaudLBitLHS**

Baud Baud rate (transmission speed)

Baud = 4 2400 baud

Baud = 5 4800 baud

Baud = 6 9600 baud (factory setting)

Baud = 7 19200 baud

Baud = 8 38400 baud

Bit Number of bits per character, parity, and number of stop bits

Bit = 0 7 bits / even parity / 1 stop bit

Bit = 1 7 bits / odd parity / 1 stop bit

Bit = 2 7 bits / no parity / 1 stop bit

Bit = 3 8 bits / no parity / 1 stop bit (factory setting)

Bit = 4 7 bits / even parity / 2 stop bits

Bit = 5 7 bits / odd parity / 2 stop bits

Bit = 6 7 bits / no parity / 2 stop bits

Bit = 7 8 bits / no parity / 2 stop bits

HS handshake (data-flow control)

HS = 0 No handshake (factory setting)

HS = 1 Software handshake (Xoff – Xon protocol)

HS = 2 Hardware handshake (CTS – RTS protocol)

Command **COM** Inquiry of the current settings for serial interface

Response **COMLAL0LBaudLBitLHS**

For an explanation of baud, bit, and HS see the above explanation of the COM command.

COMLIL Command understood, parameter incorrect, e.g. specified COM port does not exist.

Example 1 (weighing platform with fix interface)

Command **COM** Transmit the current setting of the interface parameters for the installed interface.

Response **COMLAL0L7L3L2**

Fix interface (RS-232) is set to 19200 baud, 8

Setting of the interface parameter

Command `COM_LA_0_Baud_LBit_LHS`

Set the parameters of the specified interface to the desired values.

Baud Baud rate (transmission speed)

Baud = 4 2400 baud
Baud = 5 4800 baud
Baud = 6 9600 baud (factory setting)
Baud = 7 19200 baud
Baud = 8 38400 baud

Bit Number of bits per character, parity, and number of stop bits

Bit = 0 7 bits / even parity / 1 stop bit
Bit = 1 7 bits / odd parity / 1 stop bit
Bit = 2 7 bits / no parity / 1 stop bit
Bit = 3 8 bits / no parity / 1 stop bit (factory setting)
Bit = 4 7 bits / even parity / 2 stop bits
Bit = 5 7 bits / odd parity / 2 stop bits
Bit = 6 7 bits / no parity / 2 stop bits
Bit = 7 8 bits / no parity / 2 stop bits

HS handshake (data-flow control)

HS = 0 No handshake (factory setting)
HS = 1 Software handshake (Xoff – Xon protocol)
HS = 2 Hardware handshake (CTS – RTS protocol)

Response `COM_LA`

Command successfully executed.

`COM_LL`

Command understood but cannot be executed (e.g. incorrect parameter or specified interface (COM port) not installed).

Example (weighing module with two interfaces)

Command `COM_0_8_3_0`

Setting the parameters for the serial interface to 38400 baud, 8 data bits, no parity, 1 stop bit, no handshake.

Response `COM_LA`

The parameters have been set to the desired values.

Comments

- The response takes place with the current settings, the settings are changed after the response.

COPT Command to configure interface options

Only for weighing platforms (without terminal) with interface option.

1 COPT command in general

Command	<code>COPT x1 x2 x3 x4</code>	Set/get configuration for options. x1 = Function: "start", "get", "set", "end" x2 = Identifier of the datum (get/set) x3 = Type of the datum (for "set" only) x4 = Value of the datum (for "set" only)
Start	<code>COPT start</code>	Starts configuration. Host connection is suspended.
End	<code>COPT end</code>	Ends the configuration. The host connection is resumed.
Read	<code>COPT get id</code>	Read a datum. id = Identifier of the datum
Set	<code>COPT set id typ val</code>	Sets configurations for options. id = Identifier of the datum type = Type of the datum value = Value of the datum
Responses	<code>COPT a1</code>	Command is successful. a1 = for "get" only
	<code>COPT L</code>	Illegal parameter.
	<code>COPT I</code>	Command not executable: – No option inserted. – COPT command already active. – "start" command not executed for "set" commands

Example

We change the baud rate of an RS232 option:

Command `COPT␣"start"`

Response `COPT␣A␣" "`

Command `COPT␣"get"␣"System.Infos.DeviceName"`

Response `COPT␣A␣"Value='RS232␣Option' "`

Command `COPT␣"set"␣"Connection.ConnectionList.1.Baudrate"␣"i4"␣"5"`

Response `COPT␣A␣" "`

Command `COPT␣"end"`

Response `COPT␣A␣" "`

Comments

Identifier: These are not case sensitive. For details, see option specific section.

Data type: For details, see option specific sections.

2 Specific interface option information

2.1 For all options

Following read only parameters ("get") apply to all options:

Identifiers:	Type:	Length:
" System.Infos.DeviceName "	"string"	MinLength: 0, MaxLength: 20
" System.Infos.SWNumber "	"string"	MinLength: 0, MaxLength: 20
" System.Infos.SWVersion "	"string"	MinLength: 0, MaxLength: 20

Example

Command `COPT␣"get"␣"System.Infos.DeviceName"`

Response `COPT␣A␣"Value='Ethernet␣Option' "`

2.2 RS232 option

Additional read only parameters:

Identifier:	Type:	Possible values:
" System.Infos.Status "	"i4"	"-4" ("busy") "-3" ("wrongly/not configured") "-2" ("starting up") "-1" ("out-of-order") "0" ("ready")

Additional read/write parameters:

Identifiers:	Type:	Possible values:
"Connection.ConnectionList.1.Baudrate"	"i4"	"2" ("600") "3" ("1200") "4" ("2400") "5" ("4800") "6" ("9600") "7" ("19200") "8" ("38400")
"Connection.ConnectionList.1.Parity"	"i4"	"0" ("7Bit/None") "1" ("7Bit/Even") "2" ("7Bit/Odd") "3" ("8Bit/None")
"Connection.ConnectionList.1.Handshake"	"i4"	"0" ("None") "1" ("XON/XOFF") "2" ("RTS/CTS")
"Connection.ConnectionList.1.EndOfLine"	"i4"	"1" ("<CR><LF>") "2" ("<CR>") "3" ("<LF>")
"Connection.ConnectionList.1.CharSet"	"i4"	"0" ("ANSI/WIN") "1" ("IBM/DOS")

Example

Commands:

```
COPT␣"set"␣"Connection.ConnectionList.1.Baudrate"␣"i4"␣"5"
```

```
COPT␣"set"␣"Connection.ConnectionList.1.Parity"␣"i4"␣"1"
```

Set the baud rate to 4800 Bd and the format to 7 Bit odd parity.

2.3 LC options, PS/2 option, BTS option

Additional read only parameters (the devices have no further parameters):

Identifier:	Type:	Possible values:
"System.Infos.Status"	"i4"	"-4" ("busy") "-3" ("wrongly/not configured") "-2" ("starting up") "-1" ("out-of-order") "0" ("ready")

Example

Command: COPT␣"get"␣"System.Infos.Status"

Response: COPT␣A␣"Value='0'"

2.4 Ethernet option

Additional read only parameters:

Identifier:	Type:	Possible values / Length:
"System.Infos.Status"	"i4"	"-5" ("lost DHCP lease") "-4" ("busy") "-3" ("wrongly/not configured") "-2" ("starting up") "-1" ("out-of-order") "0" ("ready")
"System.Infos.SerialNumber"	"string"	MinLength: 0, MaxLength: 10

Additional read/write interface option parameters:

Identifiers:	Type:	Possible values / Length:
"Interface.DHCP"	"i4"	"0" ("Off") "1" ("On")
"Interface.IPAddress"	"string"	MinLength: 0, MaxLength: 15
"Interface.SubnetMask"	"string"	MinLength: 0, MaxLength: 15
"Interface.StandardGateway"	"string"	MinLength: 0, MaxLength: 15
"Interface.DomainNameServer"	"string"	MinLength: 0, MaxLength: 15
"Interface.Hostname"	"string"	MinLength: 0, MaxLength: 41
"Interface.Apply"		n.a. there is no parameter.

This function has to be called for changes in parameters with an identifier starting with "Interface" to take effect. This command be called to make the previous interface parameters permanent.

Type: "" two quotes with no space between, i.e. none, its a function

Values: "" two quotes with no space between, i.e. none, its a function

Call: COPT"set" "Interface.Apply" "" ""

"Connection.ConnectionList.1.Client"	"i4"	"0" ("Off") "1" ("On")
"Connection.ConnectionList.1.ClientAddress"	"string"	MinLength: 0, MaxLength: 40
"Connection.ConnectionList.1.ClientPort"	"i14"	Min: 1, Max: 65535, Increment: 1
"Connection.ConnectionList.1.Server"	"i4"	"0" ("Off") "1" ("On")
"Connection.ConnectionList.1.ServerPort"	"i14"	Min: 1, Max: 65535, Increment: 1
"Connection.ConnectionList.1.EndOfLine"	"i4"	"0" ("<CR>") "1" ("<LF>") "2" ("<CR><LF>")
"Devices.Host.MuxID"	"string"	MinLength: 0, MaxLength: 16

Example

```
COPT␣"set"␣"Interface.DHCP"␣"i4"␣"0"
COPT␣"set"␣"Interface.IPAddress"␣"string"␣"172.24.113.7"
COPT␣"set"␣"Interface.SubnetMask"␣"string"␣"255.255.248.0"
COPT␣"set"␣"Interface.StandardGateway"␣"string"␣"172.24.112.1"
COPT␣"set"␣"Interface.Apply"␣"␣"␣"
COPT␣"set"␣"Connection.ConnectionList.1.ServerPort"␣"i4"␣"8001"
COPT␣"set"␣"Connection.ConnectionList.1.Server"␣"i4"␣"1"
```

These commands switch DHCP off and configure an IP address permanently (Apply). Then we make a server at port 8001 visible to clients.

DAT	Date
------------	-------------

Inquiry of date

Command	DAT	Inquiry of current date of the balance.
Response	DAT␣A␣dd␣mm␣yyyy	"dd␣mm␣yyyy" represents the date in the format day␣month␣year.
	DAT␣I	Inquiry of the date not possible at present as another operation is taking place.

Set date

Command	DAT␣dd␣mm␣yyyy	Set date in the format "dd␣mm␣yyyy".
Response	DAT␣A	Date has been set.
	DAT␣I	Date can not be set at present as another operation is taking place.
	DAT␣L	Command not executed as the date format was not correct.

Example

Command	DAT	Inquiry of date of the balance.
Response	DAT␣A␣01␣10␣2003	Current date of the balance is 1 October 2003.

Comment

The set date is retained even after the reset command "@".

I10 Balance ID – Inquiry of balance identification

Inquiry of balance identification

Command I10 Inquiry of balance identification.

Response I10 \square A \square "text"
 "text" represents the current balance identification (max. 20 alphanumeric characters).

Set balance identification

Command I10 \square "text" Set balance identification text accordingly.

Response I10 \square A Balance identification has been set.

I10 \square I The balance identification can not be set at present as another operation is taking place.

I10 \square L Command not executed as the name is too long (max. 20 characters).

Example

Command I10 Inquiry of balance identification.

Response I10 \square A \square "My Balance"
 Current balance identification is My Balance.

Comments

- A sequence of maximum 20 characters is possible as "text".
- The set balance identification is retained even after the reset command "@".

I11 Balance type

Command I11 Inquiry of model designation of the balance.

Response I11 \square A \square "text" "text" represents the model designation.
 I11 \square I The model designation can not be transferred at present as another operation is taking place.

Example

Command I11 Inquiry of model designation of the balance.

Response I11 \square A \square "XS204"
 The balance is a XS204.

Comment

A sequence of maximum 20 characters is possible as "text".

114 Inquiry of balance info

Command	I14LNo	Inquiry of balance info.
		No: 0 = Balance configuration 1 = Balance description 2 = SW-Identification number 3 = SW version 4 = Serial number 5 = TDNR number
Response	I14LABNoLIndexL "Info"	Index: Number of module Info: Required information
	I14LI	Parameter wrong.
	I14LI	Command not executable.

Example

Command	I14L2	Inquiry of SW-Identification number.
Response	I14LB2L1L "11670123 "	SW identification number of "bridge".
	I14LB2L2L "11670456 "	SW identification number of "terminal".
	I14LAL2L3L "11670789 "	SW identification number of "option".

M01 Inquiry/setting of weighing mode

Command	M01	Inquiry of weighing mode.
Response	M01┐A┐x	x: Weighing 0 = Normal weighing / Universal 1 = Dosing 2 = Sensor mode 3 = Check weighing
	M01┐L	Parameters are missing, the command can thus not be executed.
	M01┐I	Command not executable.
Command	M01┐x	Setting of weighing mode. See Inquiry.
Response	M01┐A	Command executed.
	M01┐L	Parameters wrong (value range, ...).
	M01┐I	Command not executable.

Example

M01 0 → **M01 A** Setting of weighing mode to normal.

M02 Inquiry/setting of environment

Command	M02	Inquiry of environment.
Response	M02┐A┐x	x: Environment 1 = Stable 2 = Standard 3 = Unstable
	M02┐L	Parameters are missing, the command can thus not be executed.
	M02┐I	Command not executable.
Command	M02┐x	Setting of environment. See Inquiry.
Response	M02┐A	Command executed.
	M02┐L	Parameters wrong (value range, ...).
	M02┐I	Command not executable.

Example

M02 3 → **M02 A** Setting of environment to "unstable".

M03 Inquiry/setting of AutoZero

Command	M03	Inquiry of AutoZero.
Response	M03 $\underline{\text{L}}$ A $\underline{\text{L}}$ x	x: Weighing 0 = AutoZero is switched off 1 = AutoZero is activated
	M03 $\underline{\text{L}}$ L	Parameters are missing, the command can thus not be executed.
	M03 $\underline{\text{L}}$ I	Command not executable.
Command	M03 $\underline{\text{L}}$ x	Setting of AutoZero. See Inquiry.
Response	M03 $\underline{\text{L}}$ A	Command executed.
	M03 $\underline{\text{L}}$ L	Parameters wrong (value range, ...).
	M03 $\underline{\text{L}}$ I	Command not executable.

Example

M03 1 → **M03 A** Setting on AutoZero function.

M04 Inquiry/setting of SmartSens functions

Command	M04	Inquiry of SmartSens functions.
Response	M04 $\underline{\text{L}}$ B $\underline{\text{L}}$ x1 $\underline{\text{L}}$ x2	x1: SmartSens number 0 = left SmartSens 1 = right SmartSens
	...	
	M04 $\underline{\text{L}}$ A $\underline{\text{L}}$ x1 $\underline{\text{L}}$ x2	x2: Function of SmartSens 0 = no function 1 = activate 1st function 2 = activate 2nd function .. = other
	M04 $\underline{\text{L}}$ L	Parameters are missing, the command can thus not be executed.
	M04 $\underline{\text{L}}$ I	Command not executable.
Command	M04 $\underline{\text{L}}$ x1 $\underline{\text{L}}$ x2	Setting of SmartSens function. See Inquiry.
Response	M04	Command executed.
	M04 $\underline{\text{L}}$ L	Parameters wrong (value range, ...).
	M04 $\underline{\text{L}}$ I	Command not executable.

Examples

M04 $\underline{\text{L}}$ **0 2** → **M04 A** Setting of first (left) SmartSens to 2nd function (e.g. door).

M04 → **M04 B 0 0** 1. SmartSens = no function.
M04 A 1 2 2. SmartSens = function 2 (e.g. zeroing).

M05 Inquiry of user list/method

Command	M05	Inquiry of user list method.
Response	M05 ␣ B ␣ No ␣ "Name"	No: Number of user method
	M05 ␣ B ...	
	M05 ␣ A ␣ No ␣ "Name"	Name: Name of user method
	M05 ␣ L	Parameters are missing, the command can thus not be executed.
	M05 ␣ I	Command not executable.

Examples

M05 → **M04 B 1 "Meier"**
M05 A 2 "Method EX-2"

M06 Inquiry/setting of current user/method number

Command	M06	Inquiry of activ user method number.
Response	M06 ␣ A ␣ x	x: Number of user method.
	M06 ␣ L	Parameters are missing, the command can thus not be executed.
	M06 ␣ I	Command not executable.
Command	M06 ␣ x	Setting of user method. See Inquiry.
		User number: 0..max. number of users.
Response	M06 ␣ A	Command executed.
	M06 ␣ L	Parameters wrong (value range, ...).
	M06 ␣ I	Command not executable.

Example

M06 3 → **M06 A** User method 3 is activated.

M08 Inquiry/setting of display brightness

Command	M08	Inquiry of display brightness.
Response	M08┐A┐x	x: 0..100, display brightness in %
	M08┐L	Parameters are missing, the command can thus not be executed.
	M08┐I	Command not executable.
Command	M08┐x	Setting of display brightness. See Inquiry.
Response	M08┐A	Command executed.
	M08┐L	Parameters wrong (value range, ...).
	M08┐I	Command not executable.

Examples

M08 55	→ M08 A	Setting of display brightness.
M08	→ M08 A 60	Inquiry of display brightness.

M09 Inquiry/setting of display contrast

Command	M09	Inquiry of display contrast.
Response	M09┐A┐x	x: 0..100, display contrast in %
	M09┐L	Parameters are missing, the command can thus not be executed.
	M09┐I	Command not executable.
Command	M09┐x	Setting of display contrast. See Inquiry.
Response	M09┐A	Command executed.
	M09┐L	Parameters wrong (value range, ...).
	M09┐I	Command not executable.

Example

M09	→ M09 A 60	Inquiry of display contrast.
------------	-------------------	------------------------------

M11 Inquiry/setting of beeper volume

Command	M11	Inquiry of beeper volume.
Response	M11L Ax	x: Beeper volume 0...100 %
	M11LL	Parameters are missing, the command can thus not be executed.
	M11LI	Command not executable.
Command	M11Lx	Setting of beeper volume. See Inquiry.
Response	M11LA	Command executed.
	M11LL	Parameters wrong (value range, ...).
	M11LI	Command not executable.

Example

M11 → **M11 A 50** Inquiry of beeper volume.

M12 Creating beeper tone

Command	M12Lx	Creating beeper tone. x: 0 = Variant 1 (e.g. 1x beep) 1 = Variant 2 2 = Variant 3
Response	M12LA	Command executed.
	M12LL	Parameters wrong (value range, ...).
	M12LI	Command not executable.

Example

M12 1 → **M12 A** Creating beeper tone.

Comment

This command creates an acoustic signal independent of the beeper volume setting (**M11** command).

M13 Inquiry/setting of Touch function

Command	M13	Inquiry of the Touch function.
Response	M13┐A┐x	x: 0 = Touch is switched off, inactivated 1 = Touch function is switched on
	M13┐L	Parameters are missing, the command can thus not be executed.
	M13┐I	Command not executable.
Command	M13┐x	Setting of Touch function. See Inquiry.
Response	M13┐A	Command executed.
	M13┐L	Parameters wrong (value range, ...).
	M13┐I	Command not executable.

Example

M13 1 → **M13 A** Switching on Touch function.

Comment

The functions of all other keys (hot keys, etc.) are unaffected.

M14 Inquiry of available languages

Command	M14	Inquiry of available languages.
Response	M14┐B┐No┐"language"	No: Number of language
	M14┐B┐...	
	M25┐A┐No┐"language"	Language: Name of language
	M14┐L	Parameters are missing, the command can thus not be executed.
	M14┐I	Command not executable.

Examples

M14 → **M14 B 0** **"English"**
M14 B 1 **"Deutsch"**
M14 A 2 **"Français"**

M15 Inquiry/setting of language

Command	M15	Inquiry of actual language.
Response	M15┐A┐x	x: Number of language
	M15┐L	Parameters are missing, the command can thus not be executed.
	M15┐I	Command not executable.
Command	M15┐x	Setting of language.
		x: Number according to available languages (command M14).
Response	M15┐A	Command executed.
	M15┐L	Parameters wrong (value range, ...).
	M15┐I	Command not executable.

Example

M15 3 → **M15 A** Language 3 is activated.

Comment

Language number: Number of the language according to the available languages (command **M14**).

M16 Inquiry/setting of standby mode
--

Command	M16	Inquiry of standby mode.
Response	M16┐A┐x	x: 0 = Mode switched off, inactivated 3 = 30 min. 4 = 60 min. 5 = 120 min. 6 = 240 min.
	M16┐L	Parameters are missing, the command can thus not be executed.
	M16┐I	Command not executable.
Command	M16┐x	Setting of standby mode. See Inquiry.
Response	M16┐A	Command executed.
	M16┐L	Parameters wrong (value range, ...).
	M16┐I	Command not executable.

Example

M16 1 → **M16 A** Setting of standby mode (5 min.).

Comments

- The balance switches to standby mode if it is not operated within 5 min.
- Operating the balance includes pressing a key, significant changes in weight, item counter > 0 or interface commands.

M17 Inquiry/setting of ProFACT time criteria

Command **M17** Inquiry of ProFACT time criteria.

Response **M17** **A** **hh** **mm** **ss** **x**

- hh: hours (00..23h)
- mm: minutes (00..59)
- ss: seconds (00..59)
- x: initiating time of ProFACT in hours, minutes and days
 - 0 = time criterion is inactivated
 - 1 = Monday
 - 2 = Tuesday
 - 4 = Wednesday
 - 8 = Thursday
 - 16 = Friday
 - 32 = Saturday
 - 64 = Sunday

M17 **L** Parameters are missing, the command can thus not be executed.

M17 **I** Command not executable.

Command **M17** **U** **hh** **mm** **ss** **x**

Setting of ProFACT time criterion.
See Inquiry.

Response **M17** **A** Command executed.

M17 **L** Parameters wrong (value range, ...).

M17 **I** Command not executable.

Example

M17 12 00 00 5 → M17 A

As 5 = 4 + 1, adjustment takes place on Mondays and Wednesday at 12:00.

M18 Inquiry/setting of ProFACT/FACT temperature criterion (Δ temp.)
--

Command	M18	ProFACT/FACT temperature criterion.
Response	M18LALx	<p>x:</p> <p>0 = Temp. criterion is switched off</p> <p>1 = 1st setting is activated</p> <p>2 = 2nd setting is activated</p> <p>3 = 3rd setting is activated</p> <p>4 = 4th setting is activated</p>
	M18LL	Parameters are missing, the command can thus not be executed.
	M18LI	Command not executable.
Command	M18Lx	Setting of ProFACT/FACT temperature criterion. See Inquiry.
Response	M18LA	Command executed.
	M18LL	Parameters wrong (value range, ...).
	M18LI	Command not executable.

Example

M18 1 → **M18 A** 1. Setting is activated.

Comments

Temperature difference is defined as the criterion. The balance is automatically adjusted if the temperature of the balance increases by the defined temperature difference.

M19 Inquiry/setting of adjustment weight

Command	M19	Inquiry of adjustment weight.
Response	M19 <code>└A┐Value└Unit</code>	Value: Value of adjustment weight Unit: Weight unit of adjustment weight = Defined unit of balance
	M19 <code>└L</code>	Parameters are missing, the command can thus not be executed.
	M19 <code>└I</code>	Command not executable.
Command	M19 <code>└Value└Unit</code>	Setting of external adjustment weight. This must be entered in the defined unit of the balance. See Inquiry.
Response	M19 <code>└A</code>	Command executed.
	M19 <code>└L</code>	Parameters wrong (value range, ...).
	M19 <code>└I</code>	Command not executable.

Example

M19 → **M19 A 100.123 g** Initiate adjustment weight.

Comments

The adjustment weight must be entered in the defined unit of the balance. This unit can be found by entering an inquiry command.

M20 Inquiry/setting of test weight

Command	M20	Inquiry of external test weight.
Response	M20┐A┐Value┐Unit	
		Value: Value of test weight
		Unit: Weight unit of test weight
	M20┐L	Parameters are missing, the command can thus not be executed.
	M20┐I	Command not executable.
Command	M20┐Value┐Unit	Setting of external test weight.
		See Inquiry.
Response	M20┐A	Command executed.
	M20┐L	Parameters wrong (value range, ...).
	M20┐I	Command not executable.

Example

M20 → **M20 A 100.123 g** Inquiry of test weight.

M21 Inquiry/setting of unit

Command	M21	Inquiry of unit.
Response	M21┐B┐Des┐Unit	
	M21┐B┐...	Des: Designation of unit
	M21┐A┐Des┐Unit	0 = Unit 1, to MT-SICS
		1 = Display unit
		2 = Info unit
		Unit: 0 = g
		1 = kg
		2 = t
		3 = mg
		4 = Microgram
		5 = Carat
		8 = Ounces
		9 = Troy ounces
		10 = Grain
		11 = Penny weight
		12 = Momme
		14 = Tael Hong Kong
		15 = Tael Singapore
		16 = Tael Taiwan
		17 = Tical
		18 = Tola
		19 = Baht
		20 ... 24 = Reserved
		25 = no unit
		26 ... Custom unit1
	M21┐L	Parameters are missing, the command can thus not be executed.
	M21┐I	Command not executable.
Command	M21┐Des┐Unit	Setting of unit(s).
		See Inquiry.
Response	M21┐A	Command executed.
	M21┐L	Parameters wrong (value range, ...).
	M21┐I	Command not executable.

Examples

M21	0 1	→	M21 A	Setting of unit 1 to "kg".
M21		→	M21 B 0 0	Inquiry of unit, unit 1 = "g".
		→	M21 B 1 3	Display unit = "mg".
		→	M21 B 2 5	Info unit = "carat".

Comments

All S commands are given in Unit 1 according to the definition of the MT-SICS. Only weight units are accepted as Unit 1.

M22 Inquiry/setting of custom unit definitions

Command **M22** Inquiry of custom unit definitions.

Response **M22┐B┐x1┐x2┐x3┐"x4"┐x5**

...

M22┐A┐x1┐x2┐x3┐"x4"┐x5

- x1: Number of custom units
1 = custom unit 1
2 = custom unit 2
x = other
- x2: Formula
0 = (weight offset) x factor
1 = factor / (weight offset)
- x3: Factor
- x4: Name of unit
- x5: Rounding step

M22┐L Parameters are missing, the command can thus not be executed.

M22┐I Command not executable.

Command **M22┐x1┐x2┐x3┐"x4"┐x5**

Setting of custom units. See Inquiry.

Response **M22┐A** Command executed.

M22┐L Parameters wrong (value range, ...).

M22┐I Command not executable.

Examples

M22	→	M22┐B┐1┐0┐15.5┐"sfr"┐0.05
		The custom unit is set ...
		M22┐A┐2┐1┐25.4┐"h1"┐0.1

M23 Inquiry/setting of readability, 1d/xd

Command	M23	Setting of readability.
Response	M23┐A┐x	x: Readability 0 = 1 d 1 = 10 d
	M23┐L	Parameters are missing, the command can thus not be executed.
	M23┐I	Command not executable.
Command	M23┐x	Setting of readability. See Inquiry.
Response	M23┐A	Command executed.
	M23┐L	Parameters wrong (value range, ...).
	M23┐I	Command not executable.

Example

M23 1 → **M23┐A** Setting readability = 10 d.

M24 Inquiry/setting of print key function

Command	M24	Inquiry of print key function.
Response	M24┐A┐x	x: Function of print key 0 = Print stable weight 1 = Print weight immediately 2 = No print function
	M24┐L	Parameters are missing, the command can thus not be executed.
	M24┐I	Command not executable.
Command	M24┐x	Setting of function of print key. See Inquiry.
Response	M24┐A	Command executed.
	M24┐L	Parameters wrong (value range, ...).
	M24┐I	Command not executable.

Example

M24 1 → **M24┐A** Setting of print key function.

M25 Inquiry of application selection

Command	M25	Inquiry of application selection.
Response	M25 ␣ B ␣ No ␣ "Name"	No: Number of application
	M25 ␣ B ␣...	
	M25 ␣ A ␣ No ␣ "Name"	Name: Name of application
	M25 ␣ L	Parameters are missing, the command can thus not be executed.
	M25 ␣ I	Command not executable.

Examples

```

M25 → M25 B 0 "Weighing"
      M25 B 2 "...
      M25 A x "Density"

```

M26 Inquiry/setting of current application

Command	M26	Inquiry of actual current application.
Response	M26 ␣ A ␣ x	x: Number of application
	M26 ␣ L	Parameters are missing, the command can thus not be executed.
	M26 ␣ I	Command not executable.
Command	M26 ␣ x	Setting application number.
Response	M26 ␣ A	x: Number according to appl. list (command M25). Command executed.
	M26 ␣ L	Parameters wrong (value range, ...).
	M26 ␣ I	Command not executable.

Example

```

M26 3 → M26 A Application 3 is activated.

```

Comment

Application number: Number of the application according to the application list (command **M25**).

M27 Inquiry of adjustment history

Command	M27	Inquiry of adjustment history.
Response	M27 <code>LB</code> <code>Nr</code> <code>tt</code> <code>mm</code> <code>jjjj</code> <code>hh</code> <code>mm</code> <code>Mode</code> <code>Wgt</code> "	
	M27 <code>LB</code> <code>...</code>	No: Number of adjustment entry
	M27 <code>LA</code> <code>Nr</code> <code>tt</code> <code>mm</code> <code>jjjj</code> <code>hh</code> <code>mm</code> <code>Mode</code> <code>Wgt</code> "	
		tt: Date, day of adjustment
		mm: Date, month
		jjjj: Date, year
		hh: Time, hour
		mm: Time, minute
		Mode: Type of adjustment
		0 = internal adjustment
		1 = external adjustment
		Wgt: Adjustment weight used "100.234 g"
	M27 <code>LI</code>	Parameters are missing, the command can thus not be executed.
	M27 <code>LI</code>	Command not executable.

Examples

M27 → **M27 A 3 1 1 2000 8 26 0 ""**
 1rd adjustment, internal

M27 B 2 14 12 1999 14 30 1 "200.1234 g"
 2nd adjustment, external

M27 B 1 14 12 1999 8 26 1 "200.1234 g"
 3st adjustment, external

M28 Inquiry of temperature probe

Command	M28	Inquiry of measured values of temperature probe.
Response	M28 ▯ B ▯ No ▯ Temp	
		No: Number of the temperature probe.
	M28 ▯ B ▯...	Temp: Temperature of the probe in °C.
	M28 ▯ A ▯ No ▯ Temp	
	M28 ▯ L	Parameters are missing, the command can thus not be executed.
	M28 ▯ I	Command not executable.

Examples

M28 → **M28 B 1 21.5** Temperature of 1st sensor in °C.
M28 A 2 23.2 Temperature of 2nd sensor in °C.

Comment

No information is available on the location, accuracy or correct temperature offset of the probe.

M29 Inquiry/setting of value release

Command	M29	Inquiry of value release.
Response	M29 ▯ A ▯ x	x: Value release
		0 = Very fast
		1 = Fast
		2 = Reliable and fast
		3 = Reliable
		4 = Very reliable
	M29 ▯ L	Parameters are missing, the command can thus not be executed.
	M29 ▯ I	Command not executable.
Command	M29 ▯ x	Setting of value release. See Inquiry.
Response	M29 ▯ A	Command executed.
	M29 ▯ L	Parameters wrong (value range, ...).
	M29 ▯ I	Command not executable.

Example

M29 3 → **M29 A** Setting of value release to "reliable".

M36 Inquiry/setting of LevelControl

Command	M36	Inquiry of LevelControl.
Response	M36 LA sts U text U ton	
		sts: Status 0 = LevelControl switched off 1 = LevelControl switched on
		text: Warning text 0 = Off 1 = Once 2 = Repeat
		sound: Warning Beep 0 = Off 1 = Once 2 = Repeat
	M36 LI	Parameters are missing, the command can thus not be executed.
	M36 LI	Command not executable.
Command	M36 U sts U text U ton	Setting of LevelControl.
		sts: Status 0 = LevelControl switched off 1 = LevelControl switched on
		text: Warning text 0 = Off 1 = Once 2 = Repeat
		sound: Warning Beep 0 = Off 1 = Once 2 = Repeat
Response	M36 LA	Command executed.
	M36 LI	Parameters are missing, the command can thus not be executed.
	M36 LI	Command not executable.

Example

M36 → **M36 A 111** LevelControl switched on. Text once and Beep once.

Example

This command is not supported by weighing platforms without terminal

M37 Inquiry/setting of draft shield door opening

Command	M37	Inquiry of door opening.
Response	M37┐A┐int	int: 0...100, door opening in % (100 = door completely open)
	M37┐L	Parameters are missing, the command can thus not be executed.
	M37┐I	Command not executable.
Command	M37┐int	Setting of door opening. int: 0..100, door opening in % (100 = door completely open)
Response	M37┐A	Command executed.
	M37┐L	Parameters are missing, the command can thus not be executed.
	M37┐I	Command not executable.

Example

M37 → **M37 A 75** Door opening is set on 75%.

Comment

Only the values available in the Set-Up of the balance are excepted.

M38 Execute reset

Command	M38 <u>int</u>	Execute Reset. int: execute Reset 0 = Reset, Clear Window 1 = Application Reset 2 = User Reset 3 = Master Reset
Response	M38 <u>A</u>	Command executed.
	M38 <u>L</u>	Parameters are missing, the command can thus not be executed.
	M38 <u>I</u>	Command not executable.

Example

M38 0 → **M38 A** Execute Reset.

Comments

- Reset: Activated process, such as Zero-ing, Taring, Calibrations, all open entry of result windows will be cancelled or closed.
- Application Reset: All applications will be reset to default.
- User Reset: All user settings will be reset to default.
- Master Reset: All balance parameters are reset to default

P100 Print out text on the printer

Command	P100L"text"	Print out "text" on the printer (24 characters/line).
Response	P100LA	Command executed.
Command	P100LI	Command can not be executed at present as another operation is taking place, e.g. taring.
Response	P100LL	Text could not be printed as, e.g. printer not switched on or connected.

Example

Command	P100L"Hallo"	Print out Hallo on the printer.
Response	P100LA	Printout has been started.

Comments

- A sequence of maximum 24 characters is admissible as text.
- The command **P100L"text"** has no influence on the function of the printer as a record printer.

P101 Send stable weight value to printer channel

Command	P101	Send stable weight value to the printer, corresponds to command S .
Response	P101LA	Command executed.
	P101LI	Command not executable as another operation is taking place.
	P101LL	Command not executable as no printer connected.

Example

Command	P101	A stable weight value appears on the printout, e.g. 200.01 g.
Response	P101LA	Command executed.

Comments

- The command **P101** has no influence on the function of the printer as a record printer.
- The draft shield closes with this command, when the "Door function" is set on "Automatic". It opens after sending a stable weight.

P102 Send weight value to printer channel immediately

Command	P102	Send weight value to the printer immediately, corresponds to command SI .
Response	P102LA	Command executed.
	P102LI	Command not executable as another operation is taking place.
	P102LL	Command not executable no printer connected.

Example

Command	P102	A stable or dynamic weight value appears on the printout, e.g. "D 200.01 g".
Response	P102LA	Command executed.

Comment

The command **P102** has no influence on the function of the printer as a record printer.

P120 Reset SmartTrac according to application

Command	P120	Reset SmartTrac according to application.
Response	P120LA	SmartTrac controlled according to last command.
	P120LI	SmartTrac can not be written to as another operation is taking place.
	P120LL	SmartTrac can not be activated as it is not present.

P121 Set SmartTrac as +/- display

Command	P121 \downarrow x1 \downarrow x2 \downarrow x3	<p>Set SmartTrac as +/- display with tolerance markers.</p> <p>x1: Value/Unit for the target value (pointer in 12 o'clock position).</p> <p>x2: Value/Unit for the +tolerance value (pointer to + marker).</p> <p>x3: Value/Unit for the -tolerance value (pointer to - marker).</p> <p>The specified weight value with unit must lie in the weighing range of the balance.</p> <p>Any weighing unit within Unit2 is admissible as unit (see operating instructions of the balance).</p>
Response	P121 \downarrow A P121 \downarrow I P121 \downarrow L	<p>SmartTrac controlled according to last command.</p> <p>SmartTrac can not be written to as another operation is taking place.</p> <p>The command could not be executed as the value range of one of the parameters has been exceeded.</p>

Example

Command	P121 \downarrow 200.00 \downarrow g \downarrow 15.00 \downarrow g \downarrow 20.00 \downarrow g	<p>Set SmartTrac with target value 200.00 g, +tolerance 15.0 g and -tolerance 20.0 g.</p>
Response	P121 \downarrow A	<p>Setting performed.</p>

P122 Activate individual pointers of the SmartTrac

Command	P122 \square x1 \square x2 \square x3	<p>Set up to 2 individual pointers of the SmartTrac.</p> <p>x1 = 0 SmartTrac without tolerance markers.</p> <p>x1 = 1 SmartTrac with tolerance markers.</p> <p>x2 ... x3 Integers from 0 to 59.</p> <p>x2 Enter position of 1st pointer.</p> <p>x3 Enter position of 2nd pointer.</p>
Response	P122 \square A	SmartTrac controlled according to last command.
	P122 \square I	SmartTrac can not be written to as another operation is taking place.
	P122 \square L	The command could not be executed as the value range of one of the parameters has been exceeded or the number of parameters is wrong.

Example

Command	P122 \square 0 \square 15 \square 30	Set pointers of the SmartTrac, similar to cross hairs.
Response	P122 \square A	Pointer set.

Comment

At least one pointer parameter must be preset.

P123 Activate SmartTrac by segments

Command	P123 \square x 1 \square x 2 \square x 3	Switch on up to four pointer segments of the SmartTrac. x 1 = 0 SmartTrac without tolerance markers. x 1 = 1 SmartTrac with tolerance markers. x 2, x 3 Integers from 0 to 59. x 2 Start of 1st segment. x 3 End of 1st segment.
Response	P123 \square A	SmartTrac controlled according to last command.
	P123 \square I	SmartTrac can not be written to as another operation is taking place.
	P123 \square L	The command could not be executed as the value range of one of the parameters has been exceeded or no SmartTrac.

Example

Command	P123 \square 0 \square 15 \square 30	Switch on the second quarter of the SmartTrac.
Response	P123 \square A	Command executed.

Comment

At least one segment parameter must be preset.

P124 Switch off SmartTrac display

Command	P124	Switch off SmartTrac display.
Response	P124 \square A	Command executed.
	P124 \square I	SmartTrac display can not be switched off as another operation is taking place.
	P124 \square L	The command can not be executed as there is no SmartTrac.

PWR Power on/off

Command	PWR \square x	Switch balance on or off. x = 0 Set balance to standby mode. x = 1 Switch balance on.
Response	PWR \square A	Balance has been switched off successfully.
	PWR \square A	Balance with the serial number according to text has been
	I4 \square A \square "text"	switched on successfully (see also I4 command).
	PWR \square I	Command not executable as the balance is, e.g. being

SIS Inquiry of current NetWeight values

Command **SIS** Inquiry of NetWeight with actual unit and weighing status.

Response **SIS** Command can't be executed now.

Response **SIS** **Net weight** **unit**
readability **step** **approval** **info**
 → at status 0 to 3

Response **SIS** **Error**
 → at status 4 to 6

- status:
- 0 = stable weight value
 - 1 = dynamic weight value
 - 2 = stab. inaccurate weight (MinWeight)
 - 3 = dyn. inaccurate weight (MinWeight)
 - 4 = overload
 - 5 = underload
 - 6 = error, not valid

Net weight: Net weight value (string)

- Unit: Actual Unit
- | | |
|---------|-----------|
| 0 = g | 10 = GN |
| 1 = kg | 11 = dwt |
| 2 = t | 12 = mo |
| 3 = mg | 13 =msg |
| 4 = ug | 14 = tl H |
| 5 = ct | 15 = tl S |
| 6 = N | 16 = tl T |
| 7 = lb | 17 = tcl |
| 8 = oz | 18 = tola |
| 9 = ozt | 20 = bath |

Readability: Amount of decimal places
 Format = 0...x (integer)

Step	Display step (integer)
	1 = "1" step 10 = "10" step
	2 = "2" step 20 = "20" step
	5 = "5" step 50 = "50" step
	100 = "100" step
Approv	Approval status (integer)
	0 not approved
	1 e = d
	10 e = 10d
	100 e = 100d
	-1 not approved with * in display
Info	Weight info (integer)
	0 without tare
	1 net with weighted tare
	2 net with tare store

Example

Command	SIS	Inquiry of weight value with actual unit and weighing status.
Response	SISLAL0L"100.00"LOL2L1L10L0	100.0(0) g
Response	SISLAL1L"10.0"LOL5L2L50L0L2	10.0 ct, carat value, with step 50, in coarse range, with tare store and unstable
Response	SISLAL6L"Error7" or SISLAL6L"	Error, not valid
Response	SISLAL4L"	Overload

Comment

Can't be used with custom unit, piece counting (PCS) or percent weighing (%).

SIRU Send weight value with currently displayed unit immediately and repeat

Command	SIRU	As the "SIR" command, but with currently displayed unit.
Response	SUSUWeightValueUnit	Command executed.
	SUDUWeightValueUnit	Command executed.
	SU+	Balance in overload range.
	SU-	Balance in underload range.
	SUI	Command not executable as balance is, e.g. being tared.

Example

Command	SIU
Response	SUDUUUUUUU12.34Ulb

SIU Send weight value with currently displayed unit immediately

Command	SIU	As the "SI" command, but with currently displayed unit.
Response	SUSUWeightValueUnit	Command executed, stable.
	SUDUWeightValueUnit	Command executed, dynamic.
	SU+	Balance in overload range.
	SU-	Balance in underload range.
	SUI	Command not executable as balance is, e.g. being tared.

Example

Command	SIU
Response	SUDUUUUUUU12.34Ulb

SNR Send stable weight value and repeat after each deflection

Command	SNR \square PresetValue \square Unit	Send current stable weight value and repeat after each deflection greater or equal to the preset value (see Comment).
Response	S \square S \square WeightValue \square Unit	Current stable weight value (1. value)
	S \square S \square WeightValue \square Unit	
	: :	Next stable weight value after preset deflection (2 value) etc.
	S \square I	Command not executable (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).
	S \square L	Command understood, parameter wrong.
	S \square +	Balance in overload range.
	S \square -	Balance in underload range.

Example

Command	SNR \square 50 \square g
Response	S \square S \square \square \square \square \square \square \square \square \square 12.34 \square g
	S \square S \square \square \square \square \square \square \square \square \square 67.89 \square g

Comment

The preset value is optional. If no value is defined, the deflection limit depends on balance readability as follows:

readability	min. deflection
0.01 mg	0.01 g
0.1 mg	0.1 g
0.001 g	1 g
0.01 g	1 g
0.1 g	1 g
1 g	5 g

SNRU Send stable weight value and repeat after each deflection

Command	SNRU \square PresetValue \square Unit	As the SNR command, but with currently displayed unit (see Comment).
Response	S \square S \square WeightValue \square Unit	Current stable weight value (1. value)
	S \square S \square WeightValue \square Unit	
	: :	Next stable weight value after preset deflection (2 value) etc.
	S \square I	Command not executable (balance is currently executing another command, e.g. taring, or timeout as stability was not reached).
	S \square L	Command understood, parameter wrong.
	S \square +	Balance in overload range.
	S \square -	Balance in underload range.

Example

Command	SNRU \square 50 \square g
Response	S \square S \square 12.34 \square g
	S \square S \square 67.89 \square g

Comment

The preset value is optional. If no value is defined, the deflection limit depends on balance readability as follows:

readability	min. deflection
0.01 mg	0.01 g
0.1 mg	0.1 g
0.001 g	1 g
0.01 g	1 g
0.1 g	1 g
1 g	5 g

SRU	Send stable weight value with currently displayed unit after deflection
------------	--



Command	SRU	As the "SR" command, but with currently displayed unit.
	SRU WeightValueUnit	
Response	S SWeightValueUnit	Command executed.
	S DWeightValueUnit	Deflection.
	S +	Balance in overload range.
	S -	Balance in underload range.
	S I	Command not executable as balance is, e.g. being tared.

Example


Command	SRU
Response	S SLLLLLLLL12.34Llb
	S DLLLLLLLL13.88Llb
	S LLLLLLLL15.01Llb

ST Send stable weight value after pressing  (transfer) key


Inquiry of the status

Command	ST	Inquiry of actual status of the ST function.
Responses	STLLO	Function inactive, no weight value is sent when  (transfer key) is pressed.
	STLAL1	Function active, weight value is sent each time when  (transfer key) is pressed.
	STLI	The current status can not be transferred at present as another operation is taking place.

Set ST function

Command	STL1	Send the current stable net weight value each time when  (transfer key) is pressed (see "S" command with MT-SICS level 0).
Responses	STL0	Stop sending weight value when transfer key is pressed.
	STLA	Command understood and successfully executed.
	STLI	Command understood, but not executable at present, e.g. balance is currently executing another function.
	STLL	Command understood, parameter wrong.

Example

Command	STL1	Activate ST function.
Responses	STLA	Command executed. When  (transfer key) is pressed:
	SLSLLLLL123.456Lg	Current net weight is 123.456g.

Comments

- **STL0** is the factory setting (default value).
- **ST** function is not active after switching on and after the reset command.

SU	Send stable weight value with currently displayed unit
-----------	---

Command	SU	As the "S" command, but with currently displayed unit.
Response	SUSUWeightValueUUnit	Command executed.
	SU+	Balance in overload range.
	SU-	Balance in underload range.
	SUI	Command not executable as balance is, e.g. being tared.

Example

Command	SU
Response	SUSUUUUUUUU12.34U1b

Comment

The draft shield closes with this command, when the "Door function" is set on "Automatic". It opens after sending a stable weight.

TIM	Time
------------	-------------

Inquiry of time

Command **TIM** Send current time of the balance.

Response **TIM└A└hh└mm└ss**
 "hh└mm└ss" represents the time in the 24-hour format (hours└minutes└seconds).

TIM└I No inquiry of the time is possible at present as another operation is taking place.

Set time

Command **TIM└hh└mm└ss** Set time in 24-hour format (hours└minutes└seconds).

Response **TIM└A** Time has been set, clock running.

TIM└I The time can not be set at present as another operation is taking place.

TIM└L Command not executed as the time format is not correct (e.g. 22└67└25) or no timer is built in.

Example

Command **TIM** Inquiry of time.

Response **TIM└A└22└56└11**
 The current time of the balance is 22 hours, 56 minutes and 11 seconds.

Comment

The time setting is retained even after the reset command "@".

TSTO Inquiry/setting of the test function

Inquiry of the test function setting

Command	TSTO	Inquiry of the setting for the test function.
Responses	TSTO┐A┐x┐"WeightValue┐Unit"	<p>x = 0 The internal weight is used for the test.</p> <p>x = 1 The external weight is used for the test.</p> <p>Weight value┐Unit Value of the external weight currently set that is requested in the test from the balance user via the display.</p>
	TSTO┐I	The current setting of the test function can not be transferred at present as another operation is taking place.

Set test configuration

Command	TSTO┐x	<p>Set test configuration of the balance.</p> <p>x = 0 Test with internal weight.</p> <p>x = 1 Test with external weight.</p>
Responses	TSTO┐A	Test configuration set.
	TSTO┐L	Wrong parameter.
	TSTO┐I	Command not executable as the balance is, e.g. being tared.

Example

Command	TSTO	Inquiry of current setting for the test and the value of the external test weight.
Response	TSTO┐A┐1┐┐┐┐┐┐┐┐┐┐2000.00┐g	<p>The current setting corresponds to the test with an external weight.</p> <p>For a test initiated with the TST2 command (see below), an external weight of 2000.00 g is needed.</p>

Comments

- The current value of the external weight can be seen in the menu under "Test" (see Operating instructions).
- With an internal test, no weight value appears.

TST1 Initiate test function in the current setting

Command	TST1	Start test function in the current setting.
First response	TST1LB	The test procedure has been started. Wait for next response (see Comment).
	TST1LI	The test function can not be executed at present as another operation is taking place. No second response follows.
Further responses	TST1LL	Test not possible. No second response follows.
	TST1L "text"	Prompt to unload and load the balance.
	TST1LAL "WeightValueLUnit"	Test procedure completed successfully. Value with unit corresponds to the deviation from the specified value displayed after the test (difference "D"). No unit is specified if the test has been performed with the internal weight.
	TST1LI	The test procedure has been aborted as, e.g. stability was not attained or wrong weights were loaded.

Comment

Commands sent to the balance during the test procedure are not processed and responded to in the appropriate manner until the test procedure is at an end.

TST2 Initiate test function with external weight

Command	TST2	Start test function with external weight. Inquiry of the weight used by means of the TST command (see above).
First response	TST2LB	The test procedure has been started. Wait for next response (see Comment).
	TST2LI	The test function can not be executed at present as another operation is taking place. No second response follows.
	TST2LL	Test not possible. No second response follows.
Further responses	TST2L"text"	Prompt to unload and load the balance.
	TST2LAL"WeightValueLUnit"	Test procedure completed successfully. Weight value with unit corresponds to the deviation from the specified value displayed in the top line after the test. (difference "D", see operating instructions of balance). The test procedure has been aborted as, e.g. stability was not attained or wrong weights were loaded.

Example

Command	TST2LI	Initiate test with external weight.
Response	TST2LB	The test procedure could be started.
	TST2L"LLLLLLLL0.00Lg"	Prompt to unload the balance.
	TST2L"LLLLLLLL200.0Lg"	Prompt to load the test weight.
	TST2L"LLLL0.00Lg"	Prompt to unload the balance.
	TST2LAL"LLLL0.01Lg"	External test completed successfully.

Comment

Commands sent to the balance during the test procedure are not processed and responded to in the appropriate manner until the test procedure is at an end.

TST3 Initiate test function with internal weight

Command	TST3	Start test function with built-in weight.
First response	TST3LB	The test procedure has been started. Wait for next response (see Comment).
	TST3LI	The test function can not be executed at present as another operation is taking place. No second response follows.
	TST3LL	Test not possible. No second response follows.
Further responses	TST3LAL	"WeightValue"
	TST3LI	Test procedure completed successfully. Value corresponds to the deviation from the specified value displayed after the test (difference "D").
		The test procedure has been aborted as, e.g. stability was not attained or wrong weights were loaded.

Example

Command	TST3	Initiate test with internal weight.
Response	TST3LB	The test procedure could be started.
	TST3LAL	"LLLLLLLLL0.01" The difference to the specified value is 0.01.

Comment

The commands received immediately after the first response are not processed and responded to in the appropriate manner until after the second response.

UPD Inquiry/setting of the update rate of the host interface

Command	UPD	Inquiry of the update rate of the host interface.
Response	UPDLALint	"int" represents the update rate per second.
	UPDLL	No parameters available.
	UPDLI	Command not executable.
Command	UPDLint	Set update rate per second of balance.
Response	UPDLA	Command executed.
	UPDLI	Command not executable.

Example

Command	UPD	Inquiry of update rate of the interface.
Response	UPDLAL20.2	The update rate per second of the interface is 20.2.

WS – Inquiry/setting of position of draft shield doors

Command	WS	Inquiry of position of draftshield doors.
Response	WSLALx	x: Position 0 = door(s) closed 1 = door open (right) 2 = door open (left) 8 = error 9 = intermediate
	WSLI	Parameters are missing, the command can thus not be executed.
	WSLI	Command not executable.
Command	WSLx	Setting of position of draftshield doors. int: See Inquiry.
Response	WSLA	Command executed.
	WSLI	Parameters wrong (value range, ...).
	WSLI	Command not executable.

Example

WS 1 → **WS A**
Doors opening or door opening to the right.

Comments

The user must ensure that the doors are in the correct position. If the doors are blocked when closing, then they return to their original position. The position can be monitored by an inquiry command.

4.4 Commands and responses MT-SICS level 3 for Excellence balances

PW	Piece Counting: Inquiry/setting of the piece weight
-----------	--

Command	PW	Inquiry of the piece weight for the piece counting application.
Response	PW┐A┐PieceWeightValue┐Unit	Current piece weight value in unit actually set under unit 1.
	PW┐I	Command not executable as the piece counting application is not active or a parameter is incorrect.
	PW┐I	Command not executable as another operation is taking place.
Command	PW┐A┐PieceWeightValue┐Unit	Setting of the piece weight value to the according preset value. The unit should correspond to the unit actually set under unit 1.
Response	PW┐A	Command executed.
	PW┐I	Command not executable as the piece counting application is not active or a parameter is incorrect.
	PW┐I	Command not executable as another operation is taking place.

Comments

- This command can only be used with the application "piece counting".
- The range of the piece weight value is specified to the balance type.

A01	Percent Weighing: Inquiry/setting of reference in %
------------	--

Command	A01	Inquiry of reference for percent weighing.
Response	A01LAX	x: Reference (in %)
	A01LL	Command not executable as the percent weighing application is not active.
	A01LI	Command not executable.
Command	A01LX	Setting of reference (in %).
Response	A01LB	Start to set reference (waiting for stable weight).
	A01LA	Command executed.
	A01LL	Incorrect parameter.
	A01LE	Setting reference aborted (not stable, over- or underload, abortkey,...).
	A01LI	Command not executable.

Example

Command	A01L100.00	Set the reference for percent weighing to 100.00 %
Response	A01LB	Reference is set, waiting for stable weight
	A01LA	Command executed.

Comment

- This command can only be used with the application "percent weighing".

A06 Dynamic Weighing: Inquiry/setting of dynamic weighing filter

Command	A06	Inquiry of dynamic weighing filter.
Response	A06LALx	x: 1 = Stable 2 = Standard 3 = Unstable
	A06LL	Command not executable as the dynamic weighing application is not active
	A06LI	Command not executable.
Command	A06LX	Setting of dynamic weighing filter.
		x: 1 = Stable 2 = Standard 3 = Unstable
Response	A06LA	Command executed.
	A06LL	Incorrect parameter.
	A06LI	Command not executable.

Example

Command	A06L2	Set dynamic weighing filter on standard
Response	A06LA	Dynamic weighing filter setting has been executed.

Comment

- This command can only be used with the application “dynamic weighing”.

A07	Dynamic Weighing: Inquiry/setting of dynamic weighing AutoStart
------------	--

Command	A07LX	Inquiry of the dynamic weighing AutoStart.
Response	A07LAx	x: 0 = off 1 = on
	A07LL	Command not executable as the dynamic weighing application is not active
	A07LI	Command not executable.
Command	A07LX	Setting of the dynamic weighing AutoStart.
		x: 0 = off 1 = on
Response	A07LA	Command executed.
	A07LL	Incorrect parameter.
	A07LI	Command not executable.

Example

Command	A07L0	Set dynamic weighing without AutoStart.
Response	A07LA	AutoStart setting has been executed.

Comment

- This command can only be used with the application "dynamic weighing".

A08 Dynamic Weighing: Inquiry/setting of dynamic weighing AutoTare

Command	A08	Inquiry of the dynamic weighing AutoTare.
Response	A08LALx	x: 0 = off 1 = on
	A08LL	Command not executable as the dynamic weighing application is not active
	A08LI	Command not executable.
Command	A08LX	Setting of dynamic weighing AutoTare.
		x: 0 = off 1 = on
Response	A08LA	Command executed.
	A08LL	Incorrect parameter.
	A08LI	Command not executable.

Example

Command	A08L0	Set dynamic weighing without AutoTare.
Response	A08LA	AutoTare setting has been executed.

Comment

- This command can only be used with the application "dynamic weighing".

A09 Dynamic Weighing: Inquiry/setting of Data Acquisition

Command	A09	Inquiry of the data acquisition for dynamic weighing.
Response	A09LALx	x: 0 = Dynamic Behaviour 1 = Time Interval
	A09LL	Command not executable as the dynamic weighing application is not active.
	A09LI	Command not executable.
Command	A09LX	Setting of the data acquisition for dynamic weighing.
		x: 0 = Dynamic Behaviour 1 = Time Interval
Response	A09LA	Command executed.
	A09LL	Incorrect parameter.
	A09LI	Command not executable.

Example

Command	A09L1	Set the data acquisition for dynamic weighing = time interval.
Response	A09LA	The data acquisition for dynamic weighing has been set to time interval.

Comment

- This command can only be used with the application "dynamic weighing".

SM0 Dynamic Weighing: Reset of all SMx commands

Command	SM0	Reset of all SMx commands.
Response	SM0LA	Command executed.
	SM0LL	Command not executable as the dynamic weighing application is not active or a parameter is incorrect.
	SM0LI	Command not executable as another operation is taking place.

Comment

- This command can only be used with the application "dynamic weighing".

Comments

- This command can only be used with the application „dynamic weighing“.
- Set AutoStart on (see command A07)
- The minimum load can be defined in the dynamic weighing application setup in the menu AutoStart.
- The recurring establishment of the start standby is cancelled by the SMO, SM1, SM2 and @ commands.

SM4 Dynamic Weighing: Inquiry/setting of Time Interval

Command	SM4	Inquiry of the time interval for dynamic weighing.
Response	SM4┐A┐Time Interval	Time interval in seconds for the data acquisition of the dynamic weighing.
	SM4┐L	Command not executable as the dynamic weighing application is not active.
	SM4┐I	Command not executable at present as another operation is taking place.
Command	SM4┐Time Interval	Setting of the time interval (1 ... 99 seconds) for dynamic weighing.
	SM4┐A	Command executed.
	SM4┐L	Incorrect parameter.
	SM4┐I	Command not executable at present as another operation is taking place.

Example

Command	SM4┐10	Set the time interval to 10 seconds.
Response	SM4┐A	The actual time interval has been set to 10 seconds.

Comments

- This command can only be used with the application "dynamic weighing".
- The time interval preset by the SM4 command is automatically set to the actual value indicated in the balance's display when the time interval parameter is selected manually.

4.5 Additional commands – FastHost for Excellence balances

B00 FastHost lists of commands

Command	B00	Inquiry of list of FastHost commands.
Response	B00┐B┐"text" ... B00┐A┐"text" B00┐L	Text: Command String. Inadmissible parameter.

Example

Command	B00	Inquiry of list of FastHost commands.
Response	B00┐B┐"B00" B00┐B┐"B01" B00┐B┐"B02" B00┐B┐"B03" B00┐B┐"B04" B00┐B┐"B05" B00┐B┐"B06" B00┐B┐"B07" B00┐A┐"B08"	Commands B00 to B08 are available.

B01 Inquiry of FastHost individual value

Command	B01	Inquiry of FastHost value.
Response	<output> B01┐L	<output>: as per format specification (cf B04). Inadmissible parameter.

Example

Command	B01	Inquiry of FastHost value.
Response	B┐S┐┐┐┐┐┐┐┐123.45┐unit	As per current format specification.

B02 FastHost start/stop continuous mode

Command	B02Lstart	Start/stop continuous mode. Start: 0 = Stop (Default) 1 = Start
Response	<output>	<output> : continuous as per sampling reduction (cf B06) and format specification (cf B04).
	B02LL	Inadmissible parameter.

Example

Command	B02L1	Inquiry of Start continuous mode.
Response	BUDL123.45Lg BLSL123.54Lg ...	Continuous and as per current format specification.

B03 FastHost inquiry of time interval of value counter

Command	B03	Inquiry of time interval of value counter.
Response	B03LALIntval	Intval: time interval (in seconds).
	B03LL	Inadmissible parameter.

Example

Command	B03	Inquiry of time interval of value counter.
Response	B03LAL0.010923	Time interval is 10.923 ms.

B04 FastHost inquiry and setting of format specification

Command	B04	Inquiry of FastHost format specification.
Response	B04␣A␣"Form"	Form: format string.

Command	B0␣A␣"Form"	Setting of FastHost format specification.
Response	B04␣A	FastHost format specification is set.
	B04␣L	Fault in format specification.
	B04␣I	Parameter not permanently saved.

Example 1

Command	B04␣"S␣%S%␣W:11:3%␣U%␣C:100%"	Setting FastHost format specification.
Response	B04␣A	FastHost format specification is set.

Example 2

Command	B04␣"Weight␣=␣%W: :2%␣U%"	Setting FastHost format specification.
Response	B04␣A	FastHost format specification is set.

Example 3

Command	B04␣"%W%␣mal␣%A%␣directly␣%WA%"	Setting FastHost format specification.
Response	B04␣A	FastHost format specification is set.

Example 4

Command	B04␣"\ "%W%\ "%␣is␣directly␣%W%"	Setting FastHost format specification.
Response	B04␣A	FastHost format specification is set.

Comments

- The format specification is permanently stored under user data.

Format table Information	Format specification	Options
Weight value (net value) with selectable representation (number of places, number of decimal places). The unit of the weight value is always in grams.	%W%	<p>%W : 10 : 5 % Total number of places and number of decimal places. Formatting is right-justified as long as the defined number of places for the output display is sufficient.</p> <p>Notes: The maximum possible number of decimal places is limited to the resolution of the balance (Default). Output is always as full range with rounding (Deltarange is not supported).</p>
Unit (as string)	%U%	<p>None Unit is fixed to gram</p>
<Absolute weight value> starting from fabrication zero point	%A%	<p>%A : 12 : 3 % Total number of places and number of decimal places. Formatting is right-justified as long as the defined number of places for the output display is sufficient.</p>
Stability information(stable / unstable)	%S%	<p>%S : D : S % Specification of the Stability value: 1st indication for instable (default D), 2nd indication for stable (default S). Only the indications A-Z, A-z, 0-9 and <space> are available</p>
Signal width as measure for the measured value stability. Note: The signal width can be used for a user specific detection stable/unstable.	%Q%	<p>%Q : 5 : 3 % Total number of places and number of decimal places. Default: 10:0. Formatting is right-justified as long as the defined number of places for the output display is sufficient.</p>
Time information Note: It concerns measured value counter with fixer counting rate and not around one time.	%C%	<p>%C : 256 % Supplies the measured value counter modulo with indicated value. Note: By default the counter runs from 0..4294967295 (4 Byte unsigned long)</p>
Temperature level of the measuring cell. Note: The temperature value in the cell isn't identical to the balancing ambient temperature! This can't be derived from the temperature of the measuring cell.	%T%	<p>%T : 5 : 2 % Total number of places and number of decimal places. Formatting is right-justified as long as the defined number of places for the output display is sufficient.</p> <p>Note: Maximum number of decimal places is limited to 3 (default)</p>
Constant Text	Insert directly	<p>% % Percent sign as constant text at the start of format or after a <Blank> \%Inverted commas are a text component</p>
Information delimitation	<Blank>	

Remarks on the format data

With overload the results %W% and %A% the value 99999999 are given.

With underload the results %W% and %A% the value -99999999 are given.

%C% concerns an internal measured value counter of signal processing. This is initialised after each PowerOn of the weighbridge with 0. With the count 4294967295 (4Byte unsigned long) an overflow occurs and the counter starts again with 0. With Lotus/Magellan balances with a time interval of measured value counter of 10.923 ms this overflow takes place on the 543rd day after PowerOn. With a „modulo-1000' formatting, the overflow counter runs only to 999 and then returns to 0 and starts again.

B05 FastHost inquiry and setting of stability criteria

Command	B05	Inquiry of FastHost stability criteria.
Response	B05┐A┐Width	Width: Signal width as stable criteria in grams.

Command	B05┐Width	Setting of of FastHost stability criteria.
Response	B05┐A	FastHost stability criteria is set.
	B05┐L	Inadmissible parameter.
	B05┐I	Parameter not permanently saved.

Example

Command	B05	Inquiry of FastHost stability criteria.
Response	B05┐A┐Width	Width of signal for stability is 0.0006 g

Comments

- The condition of stability is permanently stored.
- Only the signal width can be adjusted to determine the Fasthost stability by means of B05. The observation time for the determination of the signal width is given by the balance system. This can be different, depending on the adjustment of the balance parameters and on the type.
- If the value of the Fasthost signal width (B05) is more than 1000 times smaller than the smallest resolution of the balance, then the internal stability information of the balancing system will be sent to the stability marking, i.e. same behaviour concerning stability information as with the S-commands.

B06 FastHost inquiry and setting data sampling reduction

Command	B06	Inquiry of FastHost data sampling reduction.
Response	B06└┬┬Rate	Rate: Sampling reduction 1 = give all values (max. rate). 2 = give every second value 3 = give every third value ...
Command	B06└Rate	Setting of of FastHost data sampling reduction.
Response	B06└A	FastHost data sampling reduction is set.
	B06└L	Inadmissible parameter.
	B06└I	Parameter not permanently saved.

Example

Command	B06	Inquiry of FastHost data sampling reduction.
Response	B06└┬┬4	Width continuous sending every fourth value will be given.

Comment

Settings are permanently stored.

B07 FastHost inquiry and setting switch-on mode

Command	B07	Inquiry of FastHost switch-on mode.
Response	B07└┬┬Mode	Mode: switch-on mode. 0 = normal start 1 = automatic send after switch-on is activated.
Command	B07└Mode	Setting of FastHost switch-on mode.
Response	B07└A	FastHost switch-on mode is set.
	B07└L	Inadmissible parameter.
	B07└I	Parameter not permanently saved.

Example

Command	B07	Inquiry of FastHost switch-on mode.
Response	B07└┬┬1	Automatic send after switch-on is activated.

Comment

The switch-on mode is permanently stored.

B08	Set to zero with FastHost stability criteria
------------	---

Command	B08	Set to zero with FastHost stability criteria.
Response	B08┐A	Set to zero FastHost successfully completed.
Response	B08┐+	Command stopped due to overload.
	B08┐-	Command stopped due to underload.
	B08┐I	Command not available, stopped or timed-out.

Example

Command	B08	Set to zero with FastHost stability criteria.
Response	B08┐A	Set to zero FastHost successfully completed.

5 Special features

Parameter values after switching balance On/Off

The commands of the standard command are saved on the permanent memory of the balance. This means that all values changed via the interface are saved when the balance is switched off.

Several commands in succession

If several commands are sent in succession without waiting for the corresponding responses, it is possible that the balance confuses the sequence of command processing or ignores entire commands.

Weight unit of weight value

In response strings with a weight value, unit always signifies the unit actually set under unit 1 in the menu of the balance (exceptions see **SU**, **SIU**, **SIRU** and **SRU** commands (MT-SICS level 2)).

METTLER TOLEDO DeltaRange balances

If the fine range of DeltaRange balances has been exceeded at the time of transmission, the balance sends a weight value as response in which the tenth character is a space.

Repeat rate and timeout


The repeat rate with repeat commands and the duration of the timeout (time-limit function) depend on the balance type, see technical data of the balance in question.

6 An example

The following simple formula weighing application shows the data interchange between the computer with the formula weighing program and the balance.

A substance ($S = 55 \text{ g}$) comprising components $K1 = 100 \text{ g}$ and $K2 = 21 \text{ g}$ needs to be weighed into a beaker.

If too much or too little of the first component is weighed in, the target weight of the second component should be adjusted so that the ratio of the two components remains the same. The user is guided by the balance display and acknowledges his actions with the tare key.



- @ Reset balance.

← I2LAL"1114350697"
- KL3 Disable key function and report each keystroke.

← KLA

→ DL"BEAKER"


 Prompt "(load) BEAKER" appears in the display.
- ← DLA Response to the prompted text.

← KLC10 Acknowledges pressing of the tare key.
- T Tare weight on the balance.

← TUSUUUUU70.0000g

 Beaker weighs 70.0000 g.
- DL"C1 100g"

 Prompt, add component 1 = 100 g.



- ← D L A Response to the prompted text.
- ← K L C L 1 0 Acknowledges pressing of the tare key.

- s Send target weight of component 1.
- ← S L S L L L L L L L 1 0 5 . 0 0 0 0 L g Target weight of component 1 missed by 5 g.
- T Tare weight on the balance.
- ← T L S L L L L L L L 5 5 . 0 0 0 0 L g Contents of the tare memory, now corresponds to gross weight.

- D L " C 2 L 2 1 L g " Prompt, add component 2 = 21 g.
- ← D L A Response to the prompted text.
- ← K L C L 1 0 Acknowledges pressing of the tare key.

- D L " S u b L 7 6 L g " Display "76 g substance weighed in".
- ← D L A Response to the "Display" command.

7 What if...?

Tips from actual practice when the communication between the system (computer) and the balance does not function.

Establishing the communication

Switch the balance off with the corresponding "Off" key and then on again with the "On" key. The balance must now send identification string **I4**, e.g. **I4LAL"0123456789"**. If this is not the case, check the following points.

Connection

For bidirectional communication, at least three connecting lines are needed:

- Data line from the balance (TxD signal with RS232 interface).
- Data line to the balance (RxD signal with RS232 interface).
- Signal ground line (SG with RS232 interface).

Make sure that all these connections are in order. Check the connector pin assignment of the connection cables.

Interface parameters

For the transmission to function properly, the settings of the following parameters must match at both the computer and the balance:

- Baud rate (send/receive rate)
- Number of data bits
- Parity bit

Check the settings at both devices.

Handshake

For control of the transmission, in part separate connection lines are used (CTS/DTR). If these lines are missing or wrongly connected, the computer or balance can not send or receive data.

Check whether the balance is prevented from transmitting by handshake lines (CTS or DTR). Set the parameter "protocol" for the balance and the peripheral device to "No Handshake" or "none". The handshake lines now have no influence on the communication.

**To protect your METTLER TOLEDO product's future:
METTLER TOLEDO service assures the quality, measuring
accuracy and preservation of value of all METTLER TOLEDO
products for years to come.
Please send for full details about our attractive terms of service.
Thank you.**



* P 1 1 7 8 0 7 1 1 *

Subject to technical changes and to the availability
of the accessories supplied with the instruments.

© Mettler-Toledo GmbH 2005 11780711B Printed in Switzerland 0512/2.12

Mettler-Toledo GmbH, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland
Phone+41-44-9442211, Fax+41-44-9443060, Internet:<http://www.mt.com>