**Operating Instructions
Bedienungsanleitung
Mode d'emploi
Manual de Instrucciones
Istruzioni d'uso**

**METTLER TOLEDO
Printer P25, P26, P28
Drucker P25, P26, P28
Imprimante P25, P26, P28
Impresora P25, P26, P28
Stampante P25, P26, P28**

Page 3

English

Seite 21

Deutsch

Page 39

Français

Página 57

Español

Pagina 75

Italiano

METTLER TOLEDO

Contents

Overview	4	6 Maintenance	13
1 Introduction	5	6.1 Changing the Paper Roll	14
1.1 Capabilities of the P25, P26, P28 Printer	5	6.2 Replacing the Ribbon.....	14
1.2 Overview of Functions.....	5	6.3 Cleaning	14
2 Startup Procedure	5	6.4 Disposal	14
2.1 Safety Measures.....	5	7 Troubleshooting	15
2.2 Power Supply	6	8 Recommended Settings and Printer Models	16
2.3 Installing the Printer	6	9 Consumables and spare parts	18
3 Keystroke Functions	7	10 RS232C Interface	18
4 Menu	10	11 Technical data	19
5 Printer Test	13	12 Declaration of Conformity	20

Overview

Controls:



RS-P26

Connections:

- 1 On/Off Switch
- 2 Pilot lamp – status indicator
- 3 Control panel
- 4 AC power cord connection
- 5 Interface connector

Functions:

Key	P25	P26	P28
	<ul style="list-style-type: none"> Paper feed Menu option (+ On Switch) 	<ul style="list-style-type: none"> Paper feed Scroll up menu option 	<ul style="list-style-type: none"> Paper feed Scroll up menu option
	—	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit 	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit
	—	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option 	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option
	—	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings 	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings
	—	—	<ul style="list-style-type: none"> Totaling application
	—	—	<ul style="list-style-type: none"> Statistics application

1 Introduction

1.1 Capabilities of the P25, P26, P28 Printer

This Printer is a simple to operate dot matrix printer for METTLER TOLEDO instruments. The printout on normal paper meets the requirements of modern quality assurance systems (GLP, GMP, ISO 9001, etc.). These norms define the principles and user requirements related to measuring and test instruments. The Printers comply with the requirements of measurement traceability by precisely documenting the origin of data, as well as calibration procedures. It meets the requirement of archiving measurement results by providing a thermally stable and light resistant printout.

1.2 Overview of Functions

	RS-P25	USB-P25	RS-P26	RS-P28
Auto baudrate detection ¹⁾	✓	*	✓	✓
Balance feature "on/off" ²⁾	✓	*	✓	✓
USB Interface	–	✓	–	–
RS232 Interface	✓	–	✓	✓
Zero / Tare button	–	–	✓	✓
Date / Time	–	–	✓	✓
Statistics / Totaling	–	–	–	✓

* USB has standard settings

¹⁾ Auto baudrate detection: The Printer detects the settings of the connected device every time it is switches on/off by the key. Switch this feature off if the connected device has this feature not implemented (see section 8).

²⁾ Balance feature "on/off": This feature disables the keys on the Printer (see section 8).

English

2 Startup Procedure

2.1 Safety Measures

The Printer has been tested for the connections and intended purposes documented in these Operating Instructions. However, it is your responsibility to perform the necessary tests related to the methods and purposes of its intended use. Therefore, you should observe the following safety measures.

- The Printer must be used indoors only, and never in hazardous areas.
- The electronics of the Printer are protected against the ingress of splashing water and dust. However, since the printer housing is not watertight, the printer should not be used in the presence of liquids.
- Comply with cleaning instructions (see section 6.3), protect power plug against wet conditions.
- Never open the printer housing. If you have any problems with your printer, please contact your responsible METTLER TOLEDO service.

2.2 Power Supply

The power supply of the printer is integrated in the housing (100-240 V, 50-60Hz). This allows you to connect the Printer to your local line voltage.

Attention

- Connect the Printer only to AC Power receptacles with a connection to ground!
- Use an extension cable with an equipment grounding conductor!

2.3 Installing the Printer

The Printer detects the settings of the Instrument. Make sure that you proceed as follow:

① Connecting the printer to the instrument


- Attach the interconnect cable to the printer interface of the instrument and to the printer.
- Switch the instrument on

② Switching the printer on


- Connect the Printer to the mains power.
- Use the main switch </I/O> to power on.
- RS232 interface only: the pilot lamp will flash until an instrument is connected. If the lamp is still flashing check chapter 7.

③ Inserting paper

- Remove paper cover (pull on the back) (1)
- Insert the paper rotary axis through the core of the paper roll (2).
- Feed paper through the slot in the printing unit in the direction of the arrow (3).
- Press and hold until enough paper has been fed (4).
- Lead paper through the paper cover slot
- Replace the paper cover (1).

Notes:

- If, after the first printouts the connection settings (Baudrate, Bit, Handshake) of the balance change, the printer must be restarted (switch off / on) to start the auto-baudrate detection again.
- Autobaude rate can be switched off.



3 Keystroke Functions

Printer P25




	 Press briefly	 Press and hold down
Operation Level	<div style="display: flex; align-items: center; margin-bottom: 10px;">  <div style="margin-left: 10px;">one line feed</div> </div> <div style="display: flex; align-items: center; margin-bottom: 10px;">  <div style="margin-left: 10px;">Printer test with setting printout</div> </div> <div style="display: flex; align-items: center;">  <div style="margin-left: 10px;">On</div> </div>	<div style="display: flex; align-items: center; margin-bottom: 10px;">  <div style="margin-left: 10px;">continuous feed (until release)</div> </div> <div style="display: flex; align-items: center; margin-bottom: 10px;">  <div style="margin-left: 10px;">open menu</div> </div> <div style="display: flex; align-items: center;">  <div style="margin-left: 10px;">On</div> </div>
Menu Level	<div style="display: flex; align-items: center;">  <div style="margin-left: 10px;">next menu item</div> </div>	<div style="display: flex; align-items: center;">  <div style="margin-left: 10px;">store</div> </div>

English

Printer P26

	 Press briefly	 Press and hold down
Operation Level	 zero / tare the balance  print Date and Time  one line feed  print	 continuous feed (until release)  open menu
Menu Level	 <p>exit down up next</p>	 store

Printer P28

	 Press briefly	 Press and hold down
Operation Level	<p>  zero / tare the balance  print Date and Time  one line feed  print </p>	<p>  continuous feed (until release)  open menu </p>
Menu Level	<p>  exit  down  up  next </p>	<p>  store </p>
Function Level	<p>  Start totaling application  Start statistics application </p>	<p>  Conclude totaling application  Conclude statistics application </p>

English


4 Menu

Printer P25

 Press briefly	 Press and hold down	Printout
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; text-align: center; padding: 2px;">Settings</div> <div style="text-align: center; padding: 2px;">Reset</div> <div style="text-align: center;">↓</div> <div style="border: 1px dashed gray; padding: 5px; margin: 5px 0;"> <div style="text-align: center; padding: 2px;">Auto Bauderate</div> <div style="text-align: center; padding: 2px;">On</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">Auto Bauderate</div> <div style="text-align: center; padding: 2px;">Off</div> </div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; text-align: center; padding: 2px;">Char Set</div> <div style="text-align: center; padding: 2px;">IBM / DOS</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">IBM / DOS</div> <div style="text-align: center; padding: 2px;">Cyrillic</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">IBM / DOS</div> <div style="text-align: center; padding: 2px;">Japanese</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">ANSI / WIN</div> <div style="text-align: center; padding: 2px;">Latin 1</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">ANSI / WIN</div> <div style="text-align: center; padding: 2px;">Cyrillic</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">ANSI / WIN</div> <div style="text-align: center; padding: 2px;">Japanese</div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; text-align: center; padding: 2px;">Balance Feature</div> <div style="text-align: center; padding: 2px;">Enable</div> <div style="text-align: center;">↓</div> <div style="text-align: center; padding: 2px;">Balance Feature</div> <div style="text-align: center; padding: 2px;">Disable</div> </div>	<p style="text-align: center;">Reset the printer</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p>	<p>Settings Reset</p> <p style="text-align: center;">---- Reset done ----</p> <p>Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS</p> <p>Auto Bauderate: On -----STORED----- ... "current settings"</p> <p>Auto Bauderate: Off -----STORED----- ... "current settings"</p> <p>IBM/DOS -----STORED----- ... "current settings"</p> <p>IBM/DOS Cyrillic -----STORED----- ... "current settings"</p> <p>IBM/DOS Japanese -----STORED----- ... "current settings"</p> <p>ANSI/WIN Latin 1 -----STORED----- ... "current settings"</p> <p>ANSI/WIN Cyrillic -----STORED----- ... "current settings"</p> <p>ANSI/WIN Japanese -----STORED----- ... "current settings"</p> <p>Balance Feature: Enable -----STORED----- ... "current settings"</p> <p>Balance Feature: Disable -----STORED----- ... "current settings"</p>

Menu Option
Factory Setting

Printer P26 / P28

Examples if balance display is used:

Settings	
Reset	rESEt
List	LlSt
Header	
*No Header	no Hdr
Date/Time	d t Hdr
Date/Time/Balance	dtb Hdr
Date Day**	
18	dRY 18
Date Month**	
07	Month 07
Date Year**	
2008	YEAR2008
Time Minute**	
28	Min. 28
Time Hour**	
14	Hour 14
Date Format***	
DD.MMM.YYYY	ddMMY
MM DD.YYYY	MM dd Y
DD.MM.YYYY	ddMMY
MM/DD/YYYY	MM dd Y
YYYY-MM-DD	Y-MM-dd
Time Format***	
24:MM:SS	24:MM:SS
12:MM:SS	12:MM:SS
24.MM.SS	24MMSS
12.MM.SS	12MMSS
Baud rate	
*9600	bd 9600
19200	bd 19200
38400	bd 38400
600	bd 600
1200	bd 1200
2400	bd 2400
4800	bd 4800
Bit/Parity	
*8 Bit No	8b-no
7 Bit Even	7b-E
7 Bit No	7b-no
7 Bit Odd	7b-odd
Handshake	
*Xon/Xoff	HS OFF
HS Soft	HS Soft
HS Hard	HS HARD
Char Set	
*IBM/DOS	d05
IBM/DOS Cyrillic	d05 Cyr
IBM/DOS Japanese	d05 JAP
Ansi/Win Latin 1	11 in LAt
Ansi/Win Cyrillic	11 in Cyr
Ansi/Win Japanese	11 in JAP

English






Menu Option
Factory Setting

* factory settings

** instruments with built-in clock: not applicable, built-in clock has priority

*** instruments with built-in clock: format of printer has priority

Example for setting the day (P26 / P28 only):

 Press briefly	 Press and hold down	
P26 / P28	Balance Display	
	 	Open the menu: → Hold down the «  » key until "reset" appears.
		→ Select the menu option "Day" with the «  » key.
		→ Select the "Day" – scroll up with the «  » key. – scroll down with the «  » key.
		Save the settings: → Hold down the «  » key until "Stored" appears.
		

Applications (P28 only):

 Press briefly		 Press and hold down	
Statistics Application ¹⁾		Totaling Application ²⁾	
	→ Start "Statistics application" (LED on).		→ Start "Totaling application" (LED on).
	→ Add a new sample and press «M+».		→ Add a new sample and press «M+».
	→ If needed print subtotal press the «Print» key.		→ If needed print subtotal press the «Print» key.
	→ Print and end "Statistics application" (LED off).		→ Print and end "Totaling application" (LED off).

Important notes:

- ¹⁾ Sample weight must be in the range of 70% - 130% of the current average value.
- ¹⁾ Minimum dynamic weight changed must be 100 d.
- ^{1), 2)} Unit change during the application is not permitted.
- ^{1), 2)} Maximal 999 sample.
- ²⁾ Minimum weights + 1 d.

5 Printer Test

The Printer is equipped with a self-test that automatically prints out the character set of the printer. (standard IBM/DOS).

Starting the test

- Press the key «↑» and power on the Printer.
The printer prints out the current settings.

Terminating the test

- Power off the Printer.

6 Maintenance


Under normal conditions the Printer requires practically no maintenance. The METTLER TOLEDO service department is available if servicing is required. Please ask your responsible METTLER TOLEDO dealer for details.

6.1 Changing the Paper Roll

Insertion of a new roll of paper is described in section 2.3. See section 9 (Consumables) for the paper roll's order number.

Note: You may obtain the paper rotary axis from the old paper roll.

6.2 Replacing the Ribbon

- Remove the paper cover (1).
- Pull paper out of printing unit (2).

- To remove the ribbon cartridge, press in direction of arrow (3).
- Insert the new ribbon, and if necessary, take up slack with tensioning wheel (4). See section 9 (Consumables) for the ribbon cartridge order number.
- Insertion of a new roll of paper is described in section 2.3.

6.3 Cleaning

Since the printer housing is made of top grade, resistant materials, all commercially available mild cleaning agents may be used (see safety precautions in section 2.1).

6.4 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

If this device is passed on to other parties (for private or professional use), compliance with this regulation will be required as well.

Thank you for your contribution to environmental protection.

7 Troubleshooting

Problem / Message	Possible cause	Correction
Printout unreadable	<ul style="list-style-type: none"> Ribbon worn or entangled. Service life of printing unit at an end. 	<ul style="list-style-type: none"> Change ribbon cartridge or tense the ribbon properly. Please contact METTLER TOLEDO service.
Pilot lamp does not light	<ul style="list-style-type: none"> No power. Power supply fault. 	<ul style="list-style-type: none"> Switch printer on. Please contact METTLER TOLEDO service.
Pilot lamp flashes ---Bauderate not found--- *Please connect device and switch the printer off/on. *Alternatively you can change the settings of the printer or the connected device. Printer settings: Baudrate 9600 Bit/Parity 8 Bit None Handshake Off	<ul style="list-style-type: none"> No connection printer device. Printer did not find the settings of the connected instrument. 	<ul style="list-style-type: none"> Make sure that the printer and the instrument is connected. Switch the printer off/on. Alternatively change the settings of the printer or the connected instrument. (see chapter 8)
-----Timeout -----	<ul style="list-style-type: none"> Data transfer has not completed in the last 45 s. 	<ul style="list-style-type: none"> Move balance to more stable location. Check the balance.
----- Over Load -----	<ul style="list-style-type: none"> Balance in overload range. 	<ul style="list-style-type: none"> Check the weight of sample.
----- Under Load -----	<ul style="list-style-type: none"> Balance in underload range. Weighing pan is not in position. 	<ul style="list-style-type: none"> Check the balance.
-- Transmission Error --	<ul style="list-style-type: none"> Instrument does not understand printer instruction. 	<ul style="list-style-type: none"> Check interface parameter settings on printer and instrument.
---- Not Executable ----	<ul style="list-style-type: none"> Instrument can not execute printer instruction. 	<ul style="list-style-type: none"> Check the balance.
--- Date/Time Not Set --- ----- Memory Error -----	<ul style="list-style-type: none"> Printer error messages. 	<ul style="list-style-type: none"> Check the instrument. Reset the printer. Set date/time. Back-up battery depleted Please contact METTLER TOLEDO service.
----- Memory full -----	<ul style="list-style-type: none"> Statistics or Totaling Memory is full. 	<ul style="list-style-type: none"> Start a new Statistics or Totaling application.

Problem / Message	Possible cause	Correction
-- Load/Unload Weight --	<ul style="list-style-type: none"> No sample or the same sample is on the weighing pan. 	<ul style="list-style-type: none"> Put a new sample or remove it and put it again on the pan.
Unit change not allowed !	<ul style="list-style-type: none"> Changed unit during Statistics or Totaling application. 	<ul style="list-style-type: none"> Reset the balance. Change the unit back to the first one
-- Out of range --	<ul style="list-style-type: none"> Sample weight is out of range 	<ul style="list-style-type: none"> Sample weight must be in the range of 70%–130% of the current average value.


8 Recommended Settings and Printer Models

Connected Device	Device Settings:	Printer Settings:	Compatible Printer Models
AB-S/FACT, AB-S, AB-L, PB-S/FACT, PB-S, PB-L, AL, PL, PL-S, PL-L, JB-G/FACT, JB-C/FACT, JB-L-G, JB-L-C, JL-G, JL-C, JL-G/L, EL, AB-S/PH, PB-S/PH	Printer Mode	Factory settings	RS-P25, RS-P26, RS-P28
PG-S, PG-S/PH, SB	Printer Mode	Factory settings	RS-P25
XP, XS	Printer Mode	Factory settings	RS-P25
Katakana	Printer Mode and language ANSI/WIN Katakana	Factory settings	
Russian	Printer Mode and language ANSI/WIN Russian	CharSet IBM/DOS Cyrillic	
Chinese, Japanese	Printer Mode	Factory settings	
MJ33, HB-S43	Printer Mode	Factory settings	RS-P25
Nihongo	Printer Mode and language Nihongo	Factory settings	
Russian	Printer Mode and language Russian	CharSet IBM/DOS Cyrillic	

Connected Device	Device Settings:	Printer Settings:	Compatible Printer Models
HG63, HR83	Select on LocalCan Switchbox: 1 Printer, Baudrate:9600; Bit/Parity:8 Bit No, Handshake: Xon/Xoff	*Reset Auto baudrate off, Balancefeature off	RS-P25
Russian	Select on LocalCan Switchbox: 1 Printer, Baudrate:9600; Bit/Parity:8 Bit No, Handshake: Xon/Xoff and language Russian	*Reset Auto baudrate off, Balancefeature off, CharSet IBM/DOS Cyrillic	
QB1	Factory settings	Factory settings	RS-P25
SQC XP	Printer Mode	Factory settings	RS-P25
SQC16	Printer Mode	*Reset, Auto baudrate off	RS-P25
T50, T70, T90	USB compact printer	Factory settings	USB-P25
V20, V30, C20, C30	USB compact printer	Factory settings	USB-P25
DL15, DL22, DL28	RS-P42	*Reset Auto baudrate off, Balancefeature off	RS-P25 (Connect with RS9 F/F Cable Order number: 51 190362)
DL31, DL38, DL32, DL39	GA42	Auto baudrate off, Balancefeature off, Baudrate:1200; Bit/Parity:8 Bit No, Handshake: None	RS-P26 (Connect with RS9 F/F Cable Order number: 51 190362)
DE40, DE45, DE51, RE40D, RE50, DR40, DR45	GA42	Auto baudrate off, Balancefeature off, Baudrate:1200; Bit/Parity:8 Bit No, Handshake: Hardware	RS-P26 Connect with RS9 F/F Cable (Order number: 51 190362)
Refracto 30PX/GS, Densito 30PX	PrinterMode and Baudrate 9600	Auto baudrate off, Balancefeature off, Baudrate:9600; Bit/Parity:8 Bit, Handshake: Hardware	RS-P26
SevenEasy S20, S30, v 1.30 or higher	Factory settings	*Reset Auto baudrate off, Balancefeature off	RS-P26, RS-P25
SevenMulti S40, S47, S50, S70, S80, v1.1 or higher	Baudrate:9600; Bit/Parity:8 Bit No, Handshake: Xon/Xoff	*Reset, Auto baudrate off, Balancefeature off	RS-P25, RS-P26
SevenGo pro SG6, SG7, SG8	Factory settings of the RS-232 Infrared device	Auto baudrate off Baudrate: 1200; Bit/Parity:8 Bit No Handshake: off	RS-P26 Data transfer of SevenGo pro to the printer occurs via the Infrared-RS-232 adapter (Order number: 51302333)

*The Printer detects automatically the settings of the connected device. Therefore disconnect the Printer of the connected device and reset the Printer to set the factory settings again: Baudrate:9600; Bit/Parity:8 Bit No, Handshake: Xon/Xoff


9 Consumables and Spare Parts

Pos.	Item	Part number	Delivery quantity
1	Paper cover	12120734	1 pcs
2	Paper roll 57.5 x ø 50 to 60 mm	12120799	2 pcs
3	Paper rotary axis	12120735	1 pcs
4	Ribbon cartridge, black	12120798	1 pcs
-	RS9 M/F Cable, 1.5 m	11101051	1 pcs
-	RS9 F/F Cable, 1.5 m	51190362	1 pcs
-	USB Cable	12130716	1 pcs

10 RS232C Interface

RS-P25, RS-P26 and RS-P28 printers are equipped with an RS232C interface to connect METTLER TOLEDO instruments.

- 9-pin male connector
- Matching to other device (transmission parameters) see section 4.

11 Technical Data

Power line voltage, frequency	100 – 240 VAC, 50/60 Hz
Power consumption	150 mA (100 V), 60 mA (240 V)
Altitude	up to 4000 m
Ambient temperature	0 °C – 40 °C
Relative humidity	max. 80 % at 31 °C, decreasing linearly to 50 % at 40 °C, noncondensing
Overvoltage category	Class II
Pollution degree	2
Printing unit	Dot matrix printer 5x7, 24-character line length
Printing speed	2.3 lines per second
Ribbon cartridge	Exchangeable, black
Paper roll	Standard paper 57.5 mm x ø 50 to 60 mm, integrated in housing, commercial size
Print quality	Light resistant and thermally stable printing (GLP, GMP, ISO 9001)
Interface: RS-P25, RS-P26, RS-P28 USB-P25	RS232C USB
Printer dimensions	(WxDxH) 120 mm x 203 mm x 73 mm
Packaging dimensions	(WxDxH) 255 mm x 205 mm x 150 mm
Net weight	740 g (paper roll incl.)
Compatibility limits	<ul style="list-style-type: none"> • Function Date / Time: not applicable on balances with built-in clock (built-in clock has priority) • Function Zeroing / Taring: not applicable on all models of balances

12 Declaration of Conformity

The undersigned declare on behalf of

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that the Printer **METTLER TOLEDO P25, P26** and **P28** to which this declaration relates (serial number specified on the product) is in compliance with the below mentioned EEC Directives (including all amendments)

73/23/EEC Low Voltage Directive

89/336/EEC Electromagnetic compatibility

and that the following standards have been applied

IEC/EN61010-1:2001,

EN61326:1997+ A1:98+A2:01+A3:03 (class B)


EN61326:1997+ A1:98+A2:01+A3:03 (Minimal requirements)

for Canada, USA and Australia

**CAN/CSA-C22.2 No.61010-1-04, UL Std. No.61010A-1, FCC, Part 15, class A,
AS/NZS CISPR 22, AS/NZS 61000.4.3**

Greifensee, 15.08.2008

Mettler-Toledo AG
Laboratory & Weighing Technologies

René Lenggenhager
General Manager

Marcel Strotz
Manager SBU LAB Basic Weighing

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
4 Menü

Drucker P25

 Kurz drücken	 Lange drücken	Druckausgabe
<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p style="text-align: center;">Reset</p> <hr/> <div style="border: 1px dashed gray; padding: 5px; margin: 5px 0;"> <p style="text-align: center;">Nur RS232 Interface</p> <p style="text-align: center;">Auto Bauderate</p> <p style="text-align: center;">On</p> <p style="text-align: center;">Auto Bauderate</p> <p style="text-align: center;">Off</p> </div> <p style="text-align: center;">Char Set</p> <p style="text-align: center;">IBM / DOS</p> <p style="text-align: center;">IBM / DOS Cyrillic</p> <p style="text-align: center;">IBM / DOS Japanese</p> <p style="text-align: center;">ANSI / WIN Latin 1</p> <p style="text-align: center;">ANSI / WIN Cyrillic</p> <p style="text-align: center;">ANSI / WIN Japanese</p> <p style="text-align: center;">Balance Feature</p> <p style="text-align: center;">Enable</p> <p style="text-align: center;">Balance Feature</p> <p style="text-align: center;">Disable</p> </div>	<p style="text-align: center;">Drucker rückstellen (Reset)</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p>	<p>Settings Reset</p> <p>----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS</p> <p>Auto Bauderate: On -----STORED----- ...<i>"current settings"</i></p> <p>Auto Bauderate: Off -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS Cyrillic -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS Japanese -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Latin 1 -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Cyrillic -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Japanese -----STORED----- ...<i>"current settings"</i></p> <p>Balance Feature: Enable -----STORED----- ...<i>"current settings"</i></p> <p>Balance Feature: Disable -----STORED----- ...<i>"current settings"</i></p>

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Bit/Parität

*8 Bit No 8b-no
7 Bit Even 7b-E
7 Bit No 7b-no
7 Bit Odd 7b-odd

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**To protect your METTLER TOLEDO product's future:
METTLER TOLEDO Service assures the quality, measuring
accuracy and preservation of value of all METTLER TOLEDO
products for years to come.**

**Please send for full details about our attractive terms of
service.**

Thank you.

* 1 1 7 8 1 2 0 4 *

Subject to technical changes and to the availability
of the accessories supplied with the instruments.

© Mettler-Toledo AG 2008 11781204B Printed in China 0809/2.45

Mettler-Toledo AG, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland
Phone +41-44-944 22 11, Fax +41-44-944 30 60, Internet: <http://www.mt.com>