

M420 2-wire Transmitter, for pH, O₂

Transmitters for reliable measurements and for harsh conditions with HART® interface.

Technical Data



Short description

The M420 transmitter series is METTLER TOLEDO's solution for your most demanding process conditions in hazardous area applications. Thanks to the mixed-mode input functionality, it accepts any conventional (analog) or ISM® sensor of your choice. The easy-to-use interface with large backlit display allows for intuitive and straightforward application.

Features

- ATEX/FM* approved for X versions
- Advanced ISM® technology for easy and reliable maintenance
- HART®-communication available as a standard
- Sensocheck® real-time sensor diagnostics and Sensoface® display information
- Internal log-book (100 entries) expandable to 200 (with AuditTrail®)
- Multi-level password protection
- 2 analog outputs available as a standard
- Selection of TAN software options available

*pending




Contents

Drawings	2
Specifications M420 pH HART®	4
Specifications M420 O ₂ HART®	6
Terminal assignment M420 transmitters	8
General specifications M420 transmitters	10
Ordering information	14

METTLER TOLEDO


Drawings

Assembly



- 1 Sealing plugs (3 pieces)
- 2 Washer (1 piece),
for pipe mounting:
insert between case and nut
- 3 Pg cable glands (3 pieces)
- 4 Hinge pin (1 piece)
- 5 Enclosure screws (4 pieces)
- 6 Pg plug (1 piece)
- 7 Rubber reducer (1 piece)
- 8 Cable ties (3 pieces)
- 9 Filler plugs (3 pieces)
- 10 Hexagon nuts (5 pieces)
- 11 Jumper (2 pieces)

Mounting




- 1 Cable gland (3 pieces)
- 2 Breakthroughs for cable
gland or conduit 1/2",
Ø 21.5 mm
(2 breakthroughs)
Conduits not included!
- 3 Holes for post mounting
(4 breakthroughs)
- 4 Holes for wall mounting
(2 breakthroughs)

All dimensions in mm


Drawings

Pipe mounting with ZU 0274 bracket kit




- 1 Hose clamps with worm gear drive to DIN 3017 (2 pieces)
- 2 Pipe mount plate (1 piece)
- 3 For vertical or horizontal post/pipe mounting
- 4 Self-tapping screws (4 pieces)

Protective hood ZU 0737 for wall and pipe mounting



Panel-mount kit ZU 0738



- 1 Seal (1 piece)
- 2 Screws (4 pieces)
- 3 Position of the panel
- 4 Span pieces (4 pieces)
- 5 Threaded sleeves (4 pieces)

Panel cutout 138 x 138 mm (DIN 43700)

All dimensions in mm

Specifications

M420 transmitters for pH measurement

Transmitter	2-wire HART®	M420 pH
pH/mV input	Input pH or ORP electrodes or ISFET	
	Input	Glass electrode or ISFET
	Input	Reference electrode
	Input	ORP electrode (e. g. platinum) or working electrode for impedance measurement
Measurement range	-1500 ... +1500 mV	
Display range	pH value	-2,00 ... +16,00
	ORP	-1999 ... +1999 mV
Glass electrode input ⁴⁾	Input resistance	> 1 x 10 ¹² Ohms
	Input current	< 1 x 10 ⁻¹² A ²⁾
	Impedance measurement range	0,5 ... 1000 MOhms (± 20%)
Reference electrode input ⁴⁾	Input resistance	> 1 x 10 ¹⁰ Ohms
	Input current	< 1 x 10 ⁻¹⁰ A ²⁾
	Impedance measurement range	0,5 ... 200 kOhms (± 20%)
Measurement error ^{1,2,3)}	pH value	< 0,02 TC (temperature coefficient): 0,002 pH/K
	mV value	< 1 mV TC: 0,1 mV/K
Sensor Standardization pH^{*)}	pH calibration	
Operating modes	BUF	Calibration with Calimatic automatic buffer recognition
	MAN	Manual calibration with input of individual buffer values
	DAT	Data entry of premeasured electrodes
Calimatic Buffer Sets ^{*)}	-01- METTLER TOLEDO	2,00/4,01/7,00/9,21
	-02- Merck/Riedel de Haen	2,00/4,00/7,00/9,00/12,00
	-03- Ciba (94)	2,06/4,00/7,00/10,00
	-04- NIST Technisch	1,68/4,00/7,00/10,01/12,46
	-05- NIST Standard	1,679/4,006/6,865/9,180
	-06- HACH	4,00/7,00/10,01
	-07- WTW techn. Puffer	2,00/4,01/7,00/10,00
	-08- Hamilton	4,01/7,00/10,01
	-09- Reagecon	2,00/4,00/7,00/9,00/12,00
Zero Offset	± 200 mV (only ISFET)	
Max. calibration range	Asymmetry potential:	± 60 mV
	Slope:	80 ... 103% (47,5 ... 61 mV/pH)
Sensor Standardization ORP^{*)}	ORP calibration (Zero Offset)	
Max. calibration range	-700 ... +700 mV	

Adaptive Calibration Timer*)	Preset interval	0000 ... 9999 h
Sensocheck	Automatic monitoring of glass and reference electrode, can be disabled	
Delay time	Approx. 30 s	
Sensoface	Provides information on the electrode status Evaluation of zero/slope, response, calibration interval, wear monitor, Sensocheck, can be disabled	
Temperature Input*)	Pt 100/Pt 1000/NTC 30 kOhms*) 2-wire connection, adjustable	
Measurement range	Pt 100/Pt 1000 NTC 30 kOhms	-20,0 ... +200,0 °C (-4,0... +392,0 °F) -20,0 ... +150,0 °C (-4,0... +302,0 °F)
Adjustment range	10 K	
Resolution	0,1 °C / 1 °F	
Measurement error ^{1,2,3)}	<0,5 K (<1 K bei Pt100; <1K bei NTC >100 °C)	
Temp. compensation of process medium	Linear -19,99 ... +19,99 %/K Reference temperature 25 °C	
Power output	For operating an ISFET adapter +3 V/0,5 mA -3 V/0,5 mA	

*) User-defined
 1) According to EN 60746-1, at nominal operating conditions
 2) ± 1 count
 3) Plus sensor error
 4) at environment temperature




Transmitter	2-wire HART®	M420 O ₂
Standard Device	Sensors: InPro 6800	
Input Range	Measuring current 0 ... 600 nA resolution 10 pA	
Measurement error	<0,5 % v. M. + 0,05 nA + 0,005 nA/K	
Operation Modes	GAS DO	Measurement in gases Measurement in liquids
Display Ranges	Saturation (-10 ... 80 °C) Concentration (-10 ... 80 °C) (Dissolved oxygen) Volume concentration in gas	0,0 ... 600 % 0,00 ... 99,99 mg/l 0,00 ... 99,99 ppm 0,00 ... 99,99 Vol-%
Polarization voltage	-400 ... -1000 mV Presetting -675 mV (resolution < 5 mV)	
Accepted Guard Current	≤ 20 µA	
Traces Device	Sensors: InPro 6800/6900/6950	
Input Range I ¹⁾	Measuring current 0 ... 600 nA	Resolution 10 pA
Measurement error	<0,5 % v. M. + 0,05 nA + 0,005 nA/K	
Input Range II ¹⁾	Measuring current 0 ... 10000 nA	Resolution 166 pA
Measurement error	<0,5 % v. M. + 0,8 nA + 0,08 nA/K	
Operation Modes	GAS DO	Measurement in gases Measurement in liquids
Measurement Ranges StandardSensors «10»	Saturation (-10 ... 80 °C) Concentration (-10 ... 80 °C) (Dissolved oxygen) Volume concentration in gas	0,0 ... 600,0 % 0,00 ... 99,99 mg/l 0,00 ... 99,99 ppm 0,00 ... 99,99 Vol %
Measurement Ranges StandardSensors «01»	Saturation (-10 ... 80 °C) Concentration (-10 ... 80 °C) (Dissolved oxygen) Volume concentration in gas	0,000 ... 150,0 % 0,000 ... 9999 µg/l / 10,00 ... 20,00 mg/l 0,000 ... 9999 ppb / 10,00 ... 20,00 ppm 0,000 ... 9999 ppm / 1,000 ... 50,00 Vol %
Measurement Ranges StandardSensors «001»	Saturation (-10 ... 80 °C) Concentration (-10 ... 80 °C) (Dissolved oxygen) Volume concentration in gas	0,0 ... 150,0 % 0,000 ... 9999 µg/l / 10,00 ... 20,00 mg/l 0,000 ... 9999 ppb / 10,00 ... 20,00 ppm 0,000 ... 9999 ppm / 1,000 ... 50,00 Vol %
Polarization voltage	0 ... -1000 mV Presetting -675 mV (resolution < 5 mV)	
Accepted Guard Current	≤ 20 µA	

Sensor Adjustment ^{*)}	
Operation Modes ^{*)}	AIR Automatic calibration on air WTR Automatic calibration in air saturated water Product calibration Zero point calibration
Calibration Range Standard Sensor «10»	Zero point ±2 nA Slope 25 ... 130 nA (at 25 °C, 1013 mbar)
Calibration Range Standard Sensor «01»	Zero point ±2 nA Slope 200 ... 550 nA (at 25 °C, 1013 mbar)
Calibration Range Standard Sensor «001»	Zero point ±3 nA Slope 2000 ... 9000 nA (at 25 °C, 1013 mbar)
Calibration Timer ^{*)} Pressure Correction ^{*)}	Preset interval 0000 ... 9999 h manual 0,000 ... 9,999 bar / 999,9 kPa / 145,0 psi)
Sensocheck	Monitoring of membrane, electrolyte and sensor feed cable for short circuits/open circuits (can be disabled)
Delay time	Approx. 30 s
Sensoface	Provides information on the condition of the sensor, evaluation of zero point/slope, response time, calibration interval, Sensocheck (also wear monitor with digital sensors), can be disabled
Temperature Input	NTC 22 kOhm / NTC 30 kOhm ^{*)} 2-wire connection, adjustable
Measurement Range	-20,0 ... +150,0 °C (-4,0 ... +302,0 °F)
Adjustment Range	10 K
Resolution	0.1 °C / 1 °F
Measurement error ^{2,3,4)}	< 0,5 K (< 1 K at >100 °C)




^{*)} User-defined
 1) Automatic range switch
 2) According to EN 60746-1, at nominal operating conditions
 3) ±1 count
 4) Plus sensor error

Transmitter M420 pH


METTLER TOLEDO pH
M420 pH H
 Art. No. 52 121 405
 No. 76473 / 0000000 / 0835
 -20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany  <http://www.mt.com/pro>  



METTLER TOLEDO pH
M420 pH H OUT2
 Art. No. 52 121 406
 No. 00000 / 0000000 / 0835
 -20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany  <http://www.mt.com/pro>  

METTLER TOLEDO pH
M420 pH XH
 Art. No. 52 121 407
 No. 76474 / 0000000 / 0835
 -20 ≤ T_a ≤ +65 °C


KEMA 08 ATEX xxxxx 
 II 2(1) G Ex ib [ia] IIC T4/
 II 1 G Ex ia IIC T4
 II 1 D Ex iaD 20 IP6x T85°C/
 II 2 D Ex iaD 21 IP6x T85°C

IECEX KEM 08.xxxxx
 Ex ib [ia] IIC T4/
 Zone 0 Ex ia IIC T4
 Ex iaD 20 IP6x T85°C



Electrical data see
 Control drawing 212.002-230   0499

CH - 8902 Urdorf Made in Germany

METTLER TOLEDO pH
M420 pH XH OUT2
 Art. No. 52 121 408
 No. 00000 / 0000000 / 0835
 -20 ≤ T_a ≤ +65 °C


KEMA 08 ATEX xxxxx 
 II 2(1) G Ex ib [ia] IIC T4/
 II 1 G Ex ia IIC T4
 II 1 D Ex iaD 20 IP6x T85°C/
 II 2 D Ex iaD 21 IP6x T85°C

IECEX KEM 08.xxxxx
 Ex ib [ia] IIC T4/
 Zone 0 Ex ia IIC T4
 Ex iaD 20 IP6x T85°C

Electrical data see
 Control drawing 212.002-230   0499

CH - 8902 Urdorf Made in Germany

METTLER TOLEDO Module
pH M420i pH
 Art. No. 52 121 451
 No. 0000000 / 0835


75310 / 0000000 / 0835 Made in Germany 

pH Sensor

[Temp]


SHIELD RTD RTD (GND) ISM (DATA) ISM (GND) - 3V + 3V SG REF EL. MEAS EL.

K I H G F E D C B A

pH M420i 

METTLER TOLEDO Module
pH M420i X pH
 Art. No. 52 121 452
 No. 0000000 / 0835

M420
 KEMA 08 ATEX xxxxx
 IECEX KEM08.xxxxx
 see Control drawing 212.002-240


CH - 8902 Urdorf 75311 / 0000000 / 0835  0499


pH Sensor

[Temp]

SHIELD RTD RTD (GND) ISM (DATA) ISM (GND) - 3V + 3V SG REF EL. MEAS EL.

K I H G F E D C B A

pH M420i X 

Electrical data
 see control dwg.
 212.002-240 

Transmitter M420 O₂

METTLER TOLEDO O₂

M420 O₂ H

Art. No. 52 121 415

No. 76475 / 0000000 / 0835

-20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany

Internet <http://www.mt.com/pro>

⚠ CE

METTLER TOLEDO O₂

M420 O₂ H OUT2

Art. No. 52 121 416

No. 00000 / 0000000 / 0835

-20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany

Internet <http://www.mt.com/pro>

⚠ CE

METTLER TOLEDO O₂

M420 O₂ XH

Art. No. 52 121 417

No. 76476 / 0000000 / 0835

-20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany

KEMA 08 ATEX xxxxx

II 2(1) G Ex ib [ia] IIC T4/
II 1 G Ex ia IIC T4

II 1 D Ex iaD 20 IP6x T85°C/
II 2 D Ex iaD 21 IP6x T85°C

IECEX KEM 08.xxxxx

Ex ib [ia] IIC T4/
Zone 0 Ex ia IIC T4

Ex iaD 20 IP6x T85°C

Electrical data see Control drawing 212.002-230

⚠ CE 0499

METTLER TOLEDO O₂

M420 O₂ XH OUT2

Art. No. 52 121 418

No. 00000 / 0000000 / 0835

-20 ≤ T_a ≤ +65 °C

CH - 8902 Urdorf Made in Germany

KEMA 08 ATEX xxxxx

II 2(1) G Ex ib [ia] IIC T4/
II 1 G Ex ia IIC T4

II 1 D Ex iaD 20 IP6x T85°C/
II 2 D Ex iaD 21 IP6x T85°C

IECEX KEM 08.xxxxx

Ex ib [ia] IIC T4/
Zone 0 Ex ia IIC T4

Ex iaD 20 IP6x T85°C

Electrical data see Control drawing 212.002-230

⚠ CE 0499

METTLER TOLEDO Module

O₂ M420i

Art. No. 52 121 453

No. 0000000 / 0835

75378 / 0000000 / 0835 Made in Germany

O₂ Sensor

[Temp]

SHIELD RTD (GND) ISM (DATA) ISM (GND) GUARD ANODE REFERENCE CATHODE

K I H G F E D C B A

⚠ CE

O₂ M420i

⚠

METTLER TOLEDO Module

O₂ M420i X

Art. No. 52 121 454

No. 0000000 / 0835

75379 / 0000000 / 0835 Made in Germany

O₂ / °C

M420

KEMA 08 ATEX xxxxx

IECEX KEM08.xxxxx

see Control drawing 212.002-250

CH - 8902 Urdorf 75379 / 0000000 / 0835 0499

O₂ Sensor

[Temp]

SHIELD RTD (GND) ISM (DATA) ISM (GND) GUARD ANODE REFERENCE CATHODE

K I H G F E D C B A

⚠ CE

O₂ M420i X

Electrical data see control dwg. 212.002-250

⚠

ISM Input	«One wire»-interface for operation with ISM (digital Sensors) (6 V/Ri = approx. 1,2 kOhm)	
I Input	Supply voltage 0/4 ... 20 mA / 50 Ohm for external pressure compensation	
Measurement range	Configurable 0 ... 9,999 bar	
Characteristic	Linear	
Measurement error ^{2,4)}	< 1 % of current value + 0,1 mA	
HOLD Input	Galvanically isolated (OPTO coupler)	
Function	Switches the device into HOLD	
Switching Voltage	0 ... 2 V (AC/DC) 10 ... 30 V (AC/DC)	Inactive HOLD active
CONTROL Input	Galvanically isolated (OPTO coupler)	
Function	Switch parameter set A/B	
Switching voltage	0 ... 2 V (AC/DC) 10 ... 30 V (AC/DC)	Parameter set A Parameter set B
Output 1	loop current 4 ... 20 mA, floating, protected against wrong polarity, HART communication, feeding voltage 14 ... 30 V	
Measurement Value ^{*)}	M420 pH M420 O ₂	pH, ORP, temperature DO saturation / DO concentration / temperature
Characteristic	Linear	
Alarm ^{*)}	22 mA on error messages	
Output Filter ^{*)}	PT ₁ -filter, time constant: 0 ... 120 s	
Measurement error ³⁾	< 0,25 % of current value + 0,025 mA	
Measurement range ^{*)}	Configurable within the chosen measurement range	
Permissible Measurement Span M420 pH	pH 2,00 ... 18,00 / 200 ... 3000 mV / 20 ... 320 K / 36 ... 576 °F	
Min. Measurement Range M420 O ₂	Standard: 5 % / 0,5 mg/l (ppm) / 2 Vol % Traces: 2 % / 0,1 mg/l (ppm) / 100 ppm	
Output 2	loop current 4 ... 20 mA, floating, protected against wrong polarity, feeding voltage 14 ... 30 V	
Measurement Value ^{*)}	M420 pH M420 O ₂	pH, ORP, temperature DO saturation / DO concentration / temperature

Characteristic	Linear
Alarm*)	22 mA on error messages
Output Filter*)	PT ₁ -filter, time constant: 0 ... 120 s
Measurement error ^{2, 3, 4)}	<0,25 % of current value + 0,025 mA
Measurement range*)	Configurable within the chosen measurement range
Permissible Measurement Span M420 pH	pH 2,00 ... 18,00 / 200 ... 3000 mV / 20 ... 320 K / 36 ... 576 °F
Min. Measurement Range M420 O ₂	Standard: 5 % / 0,5 mg/l (ppm) / 2 Vol % Traces: 2 % / 0,1 mg/l (ppm) / 100 ppm
Real-time Clock	several time and date formats selectable
Power Reserve	> 5 days
Display	LC display, 7-segment with icons, backlit (white)
Main display	Character height ca. 22 mm, unit symbols 14 mm
Secondary display	Character height ca. 10 mm
Text line	14 characters, 14 segments
Sensoface	3 status indicators (friendly, neutral and sad smiley)
Mode indicators	Meas, cal, config, diag Further icons for configuration and messages
Alarm indication	Alarm icon on display, blinking
Keypad	Keys: meas, info, 4 cursor keys, enter
HART-Communication	Digital communication via FSK modulation of the output current I, device identification, measured values, status and messages, parameters, calibration, protocols.
IrDA-Interface	Infrared interface for data transmission as protocols and log book, parameterization, calibration, firmware update.
FDA 21 CFR Part 11	Access control through configurable pass codes. When the configuration is changed, a log book entry is generated a flag is set on the HART protocol. Message and log book entry when the transmitter is opened.

Diagnose Functions

Calibration Data	Calibration date, zero, slope and response time
Device self-test	Display test, automatic memory test (RAM, FLASH, EEPROM), module test
Log Book	100 events with date and time
Ext. Log Book (TAN)	AuditTrail: 200 events with date and time

Service Functions

Sensor Monitor	Display of direct, uncorrected sensor signal
Current Source	Current can be defined for output 1 and 2 (00,00 ... 22,00 mA)
IrDA	Unlocking of the IrDA functionality
Password Protection	Password assignment for menu access
Factory settings	Reset of all settings to the factory values Exception: Calibration data
TAN	Release of optional additional functions

Data retention Parameters and calibration data >10 years (EEPROM)

EMC DIN EN 61326-1 (general requirements)

Emitted interference	Class B
Immunity to interference	Industrial sector

Explosion Protection	M420	USA FM CI I Div 2 **) Canada CSA CI I Div 2 **)
	M420X	ATEX / IECEx / FM / CSA Zone 1 / CI 1 Div 1 **)

Nominal Operation Conditions

Ambient Temperature	-20 ... +65 °C (-4,0... + 149,0°F)
Transport/Storage Temperature	-20 ... +70 °C (-4,0... + 158,0°F)
Relative Humidity	10... 95 % non condensating
Feeding Voltage	14 ... 30 V

Enclosure	Molded enclosure made of PBT (polybutylene terephthalat), fiber-glass reinforced
Assembly	Wall, pipe and panel mounting
Color	Gray RAL 7001
Protection	IP67
Combustibility	UL 94 V-0
Dimensions	148 mm x 148 mm
Panel Cutout	138 mm x 138 mm nach DIN 43 700
Weight	1,2 kg (1,6 kg incl. accessories and packaging)
Cable Glands	3 breakthroughs for cable glands M20 x 1.5 2 breakthroughs for NPT 1/2" or Rigid Metallic Conduit
Wiring	Terminal block, wire cross section max. 2,5 mm ²

*) User defined

**) Pending

1) Automatic range switch

2) According to EN 60746-1, at nominal operating conditions

3) ± 1 count

4) Plus sensor error

Description	Designation	Order no.
2-wire instruments		
Transmitter M420 pH H	M420 pH H	52 121 405
Transmitter M420 pH H OUT2	M420 pH H OUT2	52 121 406
Transmitter M420 pH XH	M420 pH XH	52 121 407
Transmitter M420 pH XH OUT2	M420 pH XH OUT2	52 121 408
Transmitter M420 O ₂ H	M420 O ₂ H	52 121 415
Transmitter M420 O ₂ H OUT2	M420 O ₂ H OUT2	52 121 416
Transmitter M420 O ₂ XH	M420 O ₂ XH	52 121 417
Transmitter M420 O ₂ XH OUT2	M420 O ₂ XH OUT2	52 121 418
Software options		
Log Book	SW 420-002	52 121 466
Extended Log Book (AuditTrail)	SW 420-003	52 121 467
Oxygen trace measurement	SW 420-004	52 121 468
Current input & 2 digital inputs	SW 420-005	52 121 469
Mounting accessories		
Bracket kit	ZU 0274	52 120 741
Panel-mount kit	ZU 0738	52 121 471
Protective hood	ZU 0737	52 121 470
Replacement modules		
pH M420i module	pH M420i	52 121 451
pH M420i X module	pH M420i X	52 121 452
O ₂ M420i module	O ₂ M420i	52 121 453
O ₂ M420i X module	O ₂ M420i X	52 121 454

Sales and Service:

Australia

Mettler-Toledo Ltd.
220 Turner Street
Port Melbourne
AUS-3207 Melbourne/VIC
Phone +61 1300 659 761
Fax +61 3 9645 3935
e-mail info.mtaus@mt.com

Austria

Mettler-Toledo Ges.m.b.H.
Südrandstraße 17
A-1230 Wien
Phone +43 1 607 43 56
Fax +43 1 604 28 80
e-mail infoprocess.mtat@mt.com

Brazil

Mettler-Toledo Ind. e Com. Ltda.
Alameda Araguaia, 451
Alphaville
BR-06455-000 Barueri/SP
Phone +55 11 4166 7444
Fax +55 11 4166 7401
e-mail mettler@mettler.com.br
service@mettler.com.br

China

Mettler-Toledo Instruments
(Shanghai) Co. Ltd.
589 Gui Ping Road
Cao He Jing
CN-200233 Shanghai
Phone +86 21 64 85 04 35
Fax +86 21 64 85 33 51
e-mail mtcs@public.sta.net.cn

Croatia

Mettler-Toledo d.o.o.
Mandlova 3
HR-10000 Zagreb
Phone +385 1 292 06 33
Fax +385 1 295 81 40
e-mail mt.zagreb@mt.com

Czech Republic

Mettler-Toledo s.r.o.
Trebohosticka 2283/2
CZ-100 00 Praha 10
Phone +420 2 72 123 150
Fax +420 2 72 123 170
e-mail sales.mtcz@mt.com

Denmark

Mettler-Toledo A/S
Naverland 8
DK-2600 Glostrup
Phone +45 43 27 08 00
Fax +45 43 27 08 28
e-mail info.mtdk@mt.com

France

Mettler-Toledo
Analyse Industrielle S.A.S.
30, Boulevard de Douaumont
BP 949
F-75829 Paris Cedex 17
Phone +33 1 47 37 06 00
Fax +33 1 47 37 46 26
e-mail mtpro-f@mt.com

Germany

Mettler-Toledo GmbH
Prozeßanalytik
Ockerweg 3
D-35396 Gießen
Phone +49 641 507 333
Fax +49 641 507 397
e-mail prozess@mt.com

Great Britain

Mettler-Toledo LTD
64 Boston Road, Beaumont Leys
GB-Leicester LE4 1AW
Phone +44 116 235 7070
Fax +44 116 236 5500
e-mail enquire.mtuk@mt.com

Hungary

Mettler-Toledo Kereskedelmi KFT
Teve u. 41
HU-1139 Budapest
Phone +36 1 288 40 40
Fax +36 1 288 40 50
e-mail mthu@t-online.hu

India

Mettler-Toledo India Private Limited
Amar Hill, Saki Vihar Road
Powai
IN-400 072 Mumbai
Phone +91 22 2803 1111
Fax +91 22 2857 5071
e-mail sales.mtin@mt.com

Italy

Mettler-Toledo S.p.A.
Via Vialba 42
I-20026 Novate Milanese
Phone +39 02 333 321
Fax +39 02 356 2973
e-mail customercare.italia@mt.com

Japan

Mettler-Toledo K.K.
Process Division
4F Izumikan Sanbancho Bldg.
3-8 Sanbancho
Chiyoda-ku
JP-102-0075 Tokyo
Tel. +81 3 3222 7103
Fax +81 3 3222 7118
e-mail helpdesk.ing.jp@mt.com

Malaysia

Mettler-Toledo (M) Sdn Bhd
Bangunan Electroscon Holding, U 1-01
Lot 8 Jalan Astaka U8/84
Seksyen U8, Bukit Jelutong
MY-40150 Shah Alam Selangor
Phone +60 3 78 44 58 88
Fax +60 3 78 45 87 73
e-mail
MT-MY.CustomerSupport@mt.com

Mexico

Mettler-Toledo S.A. de C.V.
Pino No. 350, Col. Sta.
MA. Insurgentes, Col Atlampa
MX-06430 México D.F.
Phone +52 55 55 47 57 00
Fax +52 55 55 41 22 28
e-mail mt.mexico@mt.com

Poland

Mettler-Toledo (Poland) Sp.z.o.o.
ul. Poleczki 21
PL-02-822 Warszawa
Phone +48 22 545 06 80
Fax +48 22 545 06 88
e-mail polska@mt.com

Russia

Mettler-Toledo Vostok ZAO
Sretenskij Bulvar 6/1
Office 6
RU-101000 Moscow
Phone +7 495 651 98 86
Fax +7 495 621 63 53
+7 495 621 78 68
e-mail inforus@mt.com

Singapore

Mettler-Toledo (S) Pte. Ltd.
Block 28
Ayer Rajah Crescent #05-01
SG-139959 Singapore
Phone +65 6890 00 11
Fax +65 6890 00 12
+65 6890 00 13
e-mail precision@mt.com

Slovakia

Mettler-Toledo s.r.o.
Hattalova 12/A
SK-831 03 Bratislava
Phone +421 2 4444 12 20-2
Fax +421 2 4444 12 23
e-mail predaj@mt.com

Slovenia

Mettler-Toledo d.o.o.
Peske 12
SI-1236 Trzin
Phone +386 1 530 80 50
Fax +386 1 562 17 89
e-mail keith.racman@mt.com

South Korea

Mettler-Toledo (Korea) Ltd.
Yeil Building 1 & 2 F
124-5, YangJe-Dong
SeCho-Ku
KR-137-130 Seoul
Phone +82 2 3498 3500
Fax +82 2 3498 3555
e-mail Sales_MTKR@mt.com

Spain

Mettler-Toledo S.A.E.
C/Miguel Hernández, 69-71
ES-08908 L'Hospitalet de Llobregat
(Barcelona)
Phone +34 93 223 76 00
Fax +34 93 223 76 01
e-mail bcn.centralita@mt.com

Sweden

Mettler-Toledo AB
Virkesvägen 10
Box 92161
SE-12008 Stockholm
Phone +46 8 702 50 00
Fax +46 8 642 45 62
e-mail info.se@mt.com

Switzerland

Mettler-Toledo (Schweiz) GmbH
Im Langacher
Postfach
CH-8606 Greifensee
Phone +41 44 944 45 45
Fax +41 44 944 45 10
e-mail salesola.ch@mt.com

Thailand

Mettler-Toledo (Thailand) Ltd.
272 Soi Soonvijai 4
Rama 9 Rd., Bangkok
Huay Kwang
TH-10320 Bangkok
Phone +66 2 723 03 00
Fax +66 2 719 64 79
e-mail
MT-TH.CustomerSupport@mt.com

USA/Canada

Mettler-Toledo Ingold, Inc.
36 Middlesex Turnpike
Bedford, MA 01730, USA
Phone +1 781 301 8800
Freephone +1 800 352 8763
Fax +1 781 271 0681
e-mail mtprous@mt.com
ingold@mt.com



Management System
certified according to
ISO 9001 / ISO 14001

Subject to technical changes.
09/2008 © Mettler-Toledo AG
Printed in Switzerland. 52 121 473

Mettler-Toledo AG, Process Analytics
Im Hackacker 15, CH-8902 Urdorf, Switzerland
Phone + 41 44 729 62 11, Fax +41 44 729 66 36

www.mt.com/pro