

# Easy, Ergonomic Control

## Low-force Aspiration and Dispensing



### Variable Speed, Low-Forces

Easily adjustable aspiration/dispense speeds offer greater flexibility with various solutions and media, and protection for sensitive cells.



### Telescoping Head

The telescoping head on the SP+ offers up to an additional 45 mm of height to reduce arm fatigue – especially useful in fume hoods where the sliding door can limit range of motion.



### Easy to Use

The ergonomic design makes it well balanced and easy to hold. Included support frame allows the instrument to rest on the bench between uses.



### Powerful Battery

Quick-charge technology delivers 15 hours of operation on a single charge and up to one hour of operation after just 10 minutes. An LED display indicates charge level.



## Rainin SP+ Pipette Controller

### Exceptionally Easy To Use

The Rainin SP+ elegantly blends superior ergonomics and state-of-the-art technology with a host of convenient features. The unique, extendable head offers greater flexibility and improves ergonomics, particularly in cramped spaces, like under a flow hood. The SP+ can be used while attached the charging cable or set in the handy wall/shelf Hang-Up™, which includes a charging port. A built-in check valve and replaceable 0.22 micron hydrophobic membrane filter protect internal components. The adapter housing and pipette adapter are autoclavable.

# Rainin SP+ Pipette Controller

The Rainin SP+ helps reduce costs by transferring liquid volumes ranging from 0.1 mL to 100 mL, and it works with any type of serological pipette (1 mL, 2 mL, 5 mL, 10 mL, 15 mL, 25 mL, 50 mL and 100 mL).

## Specific Applications

- Cell culture
- Bacteria and yeast culture
- Sample prep (DNA/RNA/ protein extraction)

## Liquid Transfer

- Buffer
- Chemicals
- Media

## Environment

- On the bench
- In the biosafety- and chemical-hood



## Ordering Information

Product	Description	Material Number
<b>Rainin SP+ pipette controller</b>	Rainin SP+ pipette controller	30671415
<b>Magnetic Hang-up™</b>	Magnetic Hang-up	30671310
<b>Support frame</b>	Support frame	30621248
<b>Power adapter</b>	Power adapter	30671311
<b>Aspiration cone</b>	Aspiration cone	30621241
<b>Filter (non-sterilized)</b>	Filter (non-sterilized) (1 pc)	30621244
<b>Filter (sterilized)</b>	Filter (sterilized) (1 pc)	30695285
<b>Pipette adapter</b>	Pipette adapter	30621247

RAININ, Pipetting 360 and Hang-up™ are trademarks of Mettler-Toledo Rainin, LLC.

**METTLER TOLEDO Group**  
 Laboratory Division  
 Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

[www.mt.com/Rainin-SPplus](http://www.mt.com/Rainin-SPplus)

For more information

Subject to technical changes  
 © 4/2022 METTLER TOLEDO. All rights reserved.  
 30770821B (EN)  
 Rainin MarCom