

# Certificate of Quality

Rainin pipette tips are certified to be free of **RNase, DNase, DNA, Pyrogens and ATP** based on the testing specifications listed below.

Contaminants tested	Testing detection levels
RNase	$\leq 5 \times 10^{-7}$ Kunitz units/ $\mu$ L
DNase	$\leq 5 \times 10^{-5}$ Kunitz units/ $\mu$ L
DNA	< 20 pg human DNA
Pyrogens	$\leq 0.01$ EU/mL
ATP	$< 3 \times 10^{-12}$ mg/ $\mu$ L