

METTLER TOLEDO Use Case

Pharma - Analytical Laboratory



Customer

- Salutas Pharma GmbH in Magdeburg, Germany
- No.2 in production of generic drugs
- Analytical Laboratory



Application

- Quantitative analysis, according to USP



Need


- Avoiding transfer errors
- Increase yield of precious substances

Solution

- Weighing 10mg of substance directly into 50ml volumetric flasks
- No more sample transfer, when using XP56 (compared to AX26MC)
- Minimal weight qualified by a Service Technician

Conclusion

- Reducing substance cost due to less transfer error and a low minimal weight
- Increasing safety when weighing directly into volumetric flasks



Mettler-Toledo GmbH
Im Langacher
PO Box LabTec
8606 Greifensee
Switzerland

Subject to technical changes
©05/2006 Mettler-Toledo

www.mt.com/micro

For more information