InLab® 742-5m


Cond °C ISM







For portable applications in the lab and outdoor. Ideal for low conductivity range samples like pure water and ultra pure water.



Sensor type	conductivity cell
Order number	51344128
Specifications	
Specifications measuring range	0.001500 μS/cm
zero point (pH)	-
zero pomi (pm)	
temperature range	0100 °C
temperature probe	NTC 30 kΩ
ala auth ara alba a' al	alada a alad
shaft material	stainless steel
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	
cable length	5 m
connector	LTW (IP67)
cell constant (cm ⁻¹)	0.105 (certified)
cell type	2 steel poles
indicating metal	-
Storage	
dry	