


InLab® Semi-Micro

Micro sensor for small sample volumes like test tubes, vials

pH Ref °C ISM



Scale 1:1



Sensor type combined pH electrode

Order number 51343165

Specifications

pH range 0...12

zero point (pH) 7.00 ± 0.25

temperature range 0...100 °C

temperature probe -

shaft material glass

membrane glass A41

membrane resistance < 600 MΩ

reference system ARGENTHAL™

reference diaphragm open junction

reference electrolyte XEROLYT®EXTRA polymer

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)
or FRISCOLYT-B® (51340053)