


**InLab®605**

**dissolved oxygen**



**Order number** 51340291

former order number -

**Specifications**

measurement range 0 ... 200 %; 0 ... 20 mg/L  
compatible with SevenExcellence™ with DO/BOD module  
Seven2Go™ S4,  
SevenGo Duo pro™  
SG68

temperature range 0 ... 60 °C

temperature probe NTC 22 kΩ

shaft material PPS

shaft length 120 mm

shaft diameter 12 mm

connector

electrode head -

cable length 1.8 m

connector BNC/RCA (cinch)

**Storage**

dry

**Special features**

IP67

fixed cable

membrane is protected by a steel mesh

shaft made of glass fiber reinforced PPS

**Choice of accessories**

Membrane kit: 51340293 (containing 3 membrane bodies and 25 mL of electrolyte)

Electrolyte: 51340294