
LE427

pH Ref °C ISM

For semi-solid and solid samples as meat, fruits, cheese and soil





Sensor type	combined pH electrode
Order number	51340333
Specifications	
pH range	111
zero point (pH)	7.00 ± 0.25
temperature range	050 °C
temperature probe	-
shaft material	glass
membrane glass	LoT
membrane resistance	< 250MΩ
reference system	Ag / AgCl
reference diaphragm	ceramic
reference electrolyte	gel
bridge electrolyte	-
electrode head	=
cable length	1.0 m
connector	BNC
cell constant (cm ⁻¹)	-
cell type	-
indicating metal	-
Storage	
3M KCI (51350072)	