HE53

Simple Moisture Analysis

Reliable and Affordable Performance

Bright, Clear Display

Thanks to large numbers and the brilliant LCD display, information and results are easy to read in all lighting conditions.



METTLER TOLEDO

Simple 3-Step Operation

So straightforward everyone can perform moisture analysis:

- 1. Tare the empty weighing pan
- 2. Evenly distribute your sample
- 3. Press start and see the result in just a few minutes

Affordable Quality

The HE53 provides excellent value for money with halogen drying technology and a high performance weighing unit.



HE53 Halogen Moisture Analyzer

The HE53 halogen moisture analyzer is a robust and affordable instrument for determining the moisture content of common everyday samples. With on-screen assistance and easy 3-step moisture measurement process, it is simple to use, even for untrained operators. It is perfect for when you have a tight budget but require reliable performance.

With the HE53, you also get:

- Reliable results
- Intuitive operation
- Easy printer connection



Rugged Construction

With a robust metal housing and easy-to-clean stainless steel surface, the HE53 has been designed to provide you with years of reliable service.



HE53 Halogen Moisture Analyzer



Technical Data

Repeatability	0.15% with 2 g sample	
	0.05% with 10 g sample	
Drying Programs	Standard, Rapid	
Readability	1 mg, 0.01 %MC	
Method Memory	1	
Results Storage	1	
Recommended %MC	1.0-100 %MC	
Results Handling	Progress indicator	
Switch-Off Criteria	1 automatic, timed (1–120 min)	
Drying Temperature	50–160 °C (1 °C increments)	
Test and Adjust	Temperature adjustment 100 °C/160 °C	
	External test weight 50 g	
Display Modes	%MC, %DC, %AM, %AD, g	
Interfaces	RS232	
Capacity	54 g	
Display	Backlit LCD	
Dimensions ($W \times D \times H$)	183 × 334 × 161 mm (341 mm cover open)	

Accessories for Moisture Analyzers

		Material No.
Sample	Aluminium sample pans	00013865
Handling	Glass fiber filters	00214464
Quality	Certified temperature kit – HE	30134141
Management	Certified adjustment weight 50 g	30402574
Other	Protective cover	30209145
Accessories	Dust filters – HC/HE	11113883
	Dust filters housing – HC/HE	30216118
	Printer, e.g. RS-P25/01	11124300

SmartCal[™] Test Substance

	Material No.
SmartCal StarterPac	30005917
cSmartCal StarterPac	30005918
SmartCal12	30005792
cSmartCal12	30005793
SmartCal24	30005790
cSmartCal24	30005791

www.mt.com/smartcal

METTLER TOLEDO Group

Laboratory Weighing Local contact: www.mt.com/contacts

www.mt.com/moisture .

For more information

Subject to technical changes © 08/2020 METTLER TOLEDO. All rights reserved 30609986A Group MarCom RITM656512 MB