

ASTORtec AG
 Zürichstrasse 59
 8840 Einsiedeln
 Schweiz

T +41 55 418 37 37
 F +41 55 418 37 38
 info@astortec.ch
 www.astortec.ch

Certificate of Analysis

Product description: ELECTRODE CLEANER
 Package batch No.: **AT004_A / AT004_P**
 Production item No.: K23812
 Production batch No.: P00043300001
 Production date: 08.09.2021
 Expiry date: 07.09.2024
 Date / Examiner: 08.09.2021 / RWÄ

Casein Test	Test method / measuring devices	Temp.	Unit	Min	Max	Results
Aspect	Visual	transparent and very low viscous liquid				Ok
Cleaning effect	Mortar and magnetic stirrer. Process: Dissolve protein powder Casein in the electrode cleaner. Dispense 0.01 g Casein powder into a reagent glass. Add 5 ml of the cleaner mixture. Close the test tube tightly and mix under the magnetic stirrer. The powder Casein must be completely dissolved in the cleaner after 24 h! Please confirm with pass!	25°C	h	- - -	< 24	pass