



Filter Weighing Solutions

Unique Accurate Measuring

Excellence Balances equipped with a Filter Weighing Kit provide a state-of-the-art solution for emissions determination. The METTLER TOLEDO measurement system enables fast and accurate weight determination down to 0.1 μg , meeting strict automotive and environmental requirements.

METTLER TOLEDO

Filter Weighing Solutions

XP Micro and Ultra Micro Balances



Part	Filter Weighing Kit	XP-W5 Draftshield (optional)	Large Filter Weighing Kit
Filter size	Ø 30-70 mm		Ø max. 110 mm
Balances	XP Micro and Ultra Micro Balances	XP Micro and Ultra Micro Balances	XP/XS Micro and Ultra Micro Balances
Order No.	11122136	11116043	00211227

XP/XS Analytical Balances



Part	Filter Weighing Kit for Analytical Balances	ErgoClip Filter Holder
Filter size	Ø max. 110 mm	Ø max. 150 mm
Balances	XP/XS Analytical Balances	XP/XS Analytical Balances
Order No.	11140000	11140185



AntiStatic Kit

Equip your Filter Weighing Solutions with an AntiStatic Kit. The U-Electrode discharges your sample automatically before weighing – absolutely safe to the touch.

Mettler-Toledo AG
Laboratory & Weighing Technologies
CH-8606 Greifensee, Switzerland
Phone +41 44 944 22 11
Fax +41 44 944 30 60

Subject to technical changes
© 09/2009 Mettler-Toledo AG
Printed in Switzerland 11796347
MCG MarCom Greifensee

www.mt.com/filter

For more information