

Easy, Ergonomic Control

Low-force Aspiration and Dispensing



Variable Speed, Low-Forces

Easily adjustable aspiration/dispense speeds offer greater flexibility with various solutions and media, and protection for sensitive cells.



Telescoping Head

The telescoping head on the SP+ offers up to an additional 45 mm of height to reduce arm fatigue – especially useful in fume hoods where the sliding door can limit range of motion.



Easy to Use

The ergonomic design makes it well balanced and easy to hold. Included support frame allows the instrument to rest on the bench between uses.



Powerful Battery

Quick-charge technology delivers 15 hours of operation on a single charge and up to one hour of operation after just 10 minutes. An LED display indicates charge level.



Rainin SP+ Pipette Controller Exceptionally Easy To Use

The Rainin SP+ elegantly blends superior ergonomics and state-of-the-art technology with a host of convenient features. The unique, extendable head offers greater flexibility and improves ergonomics, particularly in cramped spaces, like under a flow hood. The SP+ can be used while attached the charging cable or set in the handy wall/shelf Hang-Up™, which includes a charging port. A built-in check valve and replaceable 0.22 micron hydrophobic membrane filter protect internal components. The adapter housing and pipette adapter are autoclavable.

Rainin SP+ Pipette Controller

The Rainin SP+ helps reduce costs by transferring liquid volumes ranging from 0.1 mL to 100 mL, and it works with any type of serological pipette (1 mL, 2 mL, 5 mL, 10 mL, 15 mL, 25 mL, 50 mL and 100 mL).

Specific Applications

- Cell culture
- Bacteria and yeast culture
- Sample prep (DNA/RNA/ protein extraction)

Liquid Transfer

- Buffer
- Chemicals
- Media

Environment

- On the bench
- In the biosafety- and chemical-hood



Ordering Information

Product	Description	Material Number
Rainin SP+ pipette controller	Rainin SP+ pipette controller	30671415
Magnetic Hang-up™	Magnetic Hang-up	30671310
Support frame	Support frame	30621248
Power adapter	Power adapter	30671311
Aspiration cone	Aspiration cone	30621241
Sterilized filter	Sterilized filter (1 pc)	17011745
Pipette adapter	Pipette adapter	30621247

RAININ, Pipetting 360 and Hang-up™ are trademarks of Mettler-Toledo Rainin, LLC.

METTLER TOLEDO Group
 Laboratory Division
 Local contact: www.mt.com/contacts

www.mt.com/Rainin-SPplus

For more information

Subject to technical changes
 © 4/2022 METTLER TOLEDO. All rights reserved.
 30695217B (US)
 Rainin MarCom