


LE408 - DIN

2-in-1

pH



Order number 51340342

former order number -

Specifications

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C

temperature probe -

shaft material PSU

membrane glass U

membrane resistance $< 250 \text{ M}\Omega$

reference system Ag / AgCl

reference diaphragm Ceramic

reference electrolyte 3 mol/L KCl,
AgCl saturated

bridge electrolyte -

electrode head -

cable length 1.0 m

connector DIN

cell constant (cm^{-1}) -

cell type -

indicating metal -

Storage 3 mol/L KCl

Special features

Suitable for dirty solutions

Refillable reference electrolyte

Scale 1:1