

LabX™ software comes conveniently in one DVD, with all configurations and options handled via license keys. Choose your system based on a Starter Pack. In addition, LabX is licensed per connected instrument, therefore you only need a license per instrument and not for every computer. Further Option Licenses are added based on your needs.

Starter Pack LabX Express		Starter Pack LabX Server	
LabX Express		LabX Server	
Instrument License No. 1 *		Instrument License No. 1 *	
		User Management	
		Auto Import / Export	
Included			
Optional			
User Management		Product Database	
Auto Import / Export		Statistical Evaluation	
Product Database		Report Designer	
Statistical Evaluation		System Integration	
Instrument License No. 2 (optional)		Regulation	
Instrument License No. 3 (optional)		Instrument License No. 2 (optional)	
		⋮	
		Optional up to 30 Licenses	
		⋮	
		Instrument License No. 30 (optional)	


* Starter Pack Titration includes both, one instrument license for a titrator and one for a balance.

LabX™ Software Licensing Components



Starter Pack

Each installation of LabX begins with a Starter Pack, which consists of an edition (LabX Express or LabX Server), instrument license(s), and in some cases Option Licenses. The Starter Pack includes the functionalities of method and report editing as well as data- and instrument management.



Editions

There are two editions available, LabX Express and LabX Server. LabX Express edition is an installation of LabX software on one computer. Up to three instruments can be connected.

LabX Server edition is capable of a distributed installation on multiple computers and advanced options are available. Up to 30 instruments can be connected.



Instrument Licenses

LabX is licensed per connected instrument, i.e. each instrument in the installation requires an Instrument License. Independent of the Starter Pack, all types of instruments can be added to the system.



Option Licenses

Customize the LabX system to fit your application. There are several options available. For example, the "Product Database" option allows you to handle products in a professional manner and to run your methods based on specific product data. The option "Regulation" ensures full support of 21 CFR part 11 compliance and is available for LabX Server Editions only.

LabX™ Descriptions and Order Numbers

LabX™ Starter Packs

Express Edition

Full version of LabX for installation on a single computer. Up to 3 instruments can be connected to this edition.

Name	Description	Part. No.
Balance	1 Balance instrument license.	11153120
LiquiPhysics	1 LiquiPhysics instrument license.	11153130
Thermal Values	1 Thermal Values instrument license.	30005779
Quantos	1 Quantos instrument license.	30062403
Titration	1 Titration instrument license. 1 Balance instrument license.	30097754
UV/VIS	1 UV/VIS instrument license.	30247984
SevenExcellence	1 SevenExcellence instrument license.	30247987

Server Edition

Full version of LabX for installation on multiple computers and available in distributed network. Up to 30 instruments can be connected to this edition.

Name	Description	Part. No.
Balance	1 Balance instrument license. User Management Option. Auto Import/Export Option.	11153121
LiquiPhysics	1 LiquiPhysics instrument license. User Management Option. Auto Import/Export Option.	11153131
Thermal Values	1 Thermal Values instrument license. User Management Option. Auto Import/Export Option.	30005810
Quantos	1 Quantos instrument license. User Management Option. Auto Import/Export Option.	30062404
Titration	1 Titration instrument license. 1 Balance instrument license. User Management Option. Auto Import/Export Option.	30097755
UV/VIS	1 UV/VIS instrument license. User Management Option. Auto Import/Export Option.	30247985
SevenExcellence	1 SevenExcellence instrument license. User Management Option. Auto Import/Export Option.	30247988

Instrument Licenses

For each additional instrument to be connected to the LabX installation, an Instrument License is needed.

Name	Description	Part. No.
1 Excellence Balance	License for 1 Excellence balance	11153220
5 Excellence Balances	License for 5 Excellence balances	11153221
1 LiquiPhysics	License for 1 LiquiPhysics instrument	11153250
1 Thermal Values	License for 1 Thermal Values instrument	30005811
1 Quantos	License for 1 Quantos instrument	30062405
1 Titrator	License for 1 Titration instrument	30097756
3 Titrators	License for 3 Titration instruments	30097757
1 UV/VIS	License for 1 UV/VIS instrument	30247986
1 SevenExcellence	License for 1 SevenExcellence instrument	30247989

Options

The LabX installation is fully extendable. See the options & functionalities below that can be added to the installation. Some options require a LabX Server edition.

Name	Description	Part. No.
User Management	Centrally maintain and control users and roles within LabX.	11153102
Product Database	Manage product data within LabX centrally – run methods based on product data.	11153103
Auto Import / Export	LabX allows you to import and export data into and from CSV and XML formats.	11153105
Statistical Evaluation	Analyze measurements and monitor process changes, taking into account applicable regulations.	30262374
Regulation ¹⁾	LabX fully supports regulatory requirements, including 21 CFR Part 11.	11153104
Report Designer ¹⁾	Freely customize your reports with all details and possibilities including charts, graphics, and tables.	11153106
System Integration ¹⁾	Integrate LabX installation with other software packages including LIMS and ERPs using API webservice.	11153107

¹⁾ Option can be combined with LabX Server editions only

LabX™ Update and Upgrade Options

LabX systems can be updated to a newer version or upgraded from previous system.

Name	Description	Part. No.
LabX Express - Server	Sidegrade a LabX Express edition installation into a LabX Server edition.	11153117
LabX Update DVD Pack	Most recent software DVD and Quick Guide for updating an existing LabX system.	30005812
LabX Titration Upgrade Light to Express	Upgrade existing LabX Light version 3.1 to LabX Express edition.	30094214
LabX Titration Upgrade Pro to Server	Upgrade existing LabX Pro version 3.1 to LabX Server edition.	30094215

LabX™ Validation Manuals

LabX validation manuals are the tools for efficient LabX system validation.

Name	Description	Part. No.
LabX Validation Manual 1	Description of the organization, development procedures and quality management of METTLER TOLEDO.	30003640
Starter Pack Balance	Starter Pack Balance Facilitates the IQ/OQ process of a LabX installation with connected balance(s).	30003619
Starter Pack LiquiPhysics	Facilitates the IQ/OQ process of a LabX installation with connected LiquiPhysics instrument(s).	51710825
Starter Pack Thermal Values	Facilitates the IQ/OQ process of a LabX installation with connected Thermal Values instrument(s).	51710900
Starter Pack Quantos	Facilitates the IQ/OQ process of a LabX installation with connected Quantos instrument(s).	30068906
Starter Pack Titration	Facilitates the IQ/OQ process of a LabX installation with connected Titration instrument(s).	30097758
Starter Pack UV/VIS	Facilitates the IQ/OQ process of a LabX installation with connected UV/VIS instrument(s).	30253826
Starter Pack Seven-Excellence	Facilitates the IQ/OQ process of a LabX installation with connected SevenExcellence instrument(s).	30253828
Chapter Balance	Add the chapter balance to an existing non-balance Validation Manual 2 Starter Pack.	30003641
Chapter LiquiPhysics	Add the chapter LiquiPhysics to an existing non-LiquiPhysics Validation Manual 2 Starter Pack.	51710826
Chapter Thermal Values	Add the chapter Thermal Values to an existing non-Thermal Values Validation Manual 2 Starter Pack.	51710899
Chapter Quantos	Add the chapter Quantos to an existing non-Quantos Validation Manual 2 Starter Pack.	30068907
Chapter Titration	Add the chapter Titration to an existing non-Titration Validation Manual 2 Starter Pack.	30097759
Chapter UV/VIS	Add the chapter UV/VIS to an existing non-UV/VIS Validation Manual 2 Starter Pack.	30253827
Chapter SevenExcellence	Add the chapter SevenExcellence to an existing non-SevenExcellence Validation Manual 2 Starter Pack.	30253829
IPac LabX	Start working right away and be sure your system is properly configured with the Starter Pac solution.	51710898

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 06/2016
51725281H



51725281