


pH



Order number 30014096

former order number -

Specifications

pH range 0 - 14
zero point (pH) 7.00 ± 0.25

temperature range 0 - 100 °C
temperature probe NTC 30 kΩ

shaft material PEEK
membrane glass U
membrane resistance < 250 MΩ

reference system ARGENTHAL™
reference diaphragm double open junction
reference electrolyte XEROLYT® polymer

bridge electrolyte -

electrode head -
cable length 1.2 m
connector BNC / RCA (Cinch)

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage FRISCOLYT-B®
(51340053)

Special features

unbreakable, chemically resistant PEEK shaft
two open reference junctions
integrated temperature probe NTC 30 kΩ
compatible with appropriate METTLER TOLEDO
meters with Intelligent Sensor Management
(ISM) functionality

Scale 1:1