

Simple Data Management EasyDirect Balance Software



Collect data automatically

Collect weighing data from up to 10 balances via Ethernet, USB-C or RS232. Minimize manual transcription errors, save time and store data securely.



Efficient results handling

Easily review results and filter by date, instrument, user or sample. Visualize results in charts to assess target and tolerance ranges and perform statistics for efficient trend and production analysis.



Easy reporting and export

Generate simple and clear reports of weighing results. Export data to PC in various formats (XML, CSV, XLSX or PDF) or print on a network printer.



Simple instrument overview

By providing adjustment, test and service status information, EasyDirect Balance gives you a simple operational overview of all your connected instruments.



Weighing Data-Handling Made Easy

Connect your balances to EasyDirect Balance PC software and improve the management and handling of your data. Installation and setup are simple and fast, and operation is intuitive.

Easy data management will help improve the efficiency of your weighing tasks. With EasyDirect Balance, you collect weighing data automatically, making it easier to create reports, visualize results in control charts, perform statistics and have a clear overview about the status of your balances.

Download your free 45-day trial version at:

► www.mt.com/EasyDirectBalance

METTLER TOLEDO

Technical Description

	MX, MR ¹⁾	MA, LA, MS, ML, ME-T, ME, PL-E, JP, JS, JE, JET, JL-GE, XP, XPE, XS, XSE	MS-TS, ML-T ¹⁾
Supported Instruments			
Connectivity	Ethernet, RS232 ²⁾ , USB-C	RS232	RS232, Ethernet
Data Collection			
Balance Information			
Basic balance information	•	•	•
Level, adjustment, test and service status			•
Balance usage	•		•
Weighing Data			
Weighing, counting, percent weighing ³⁾ , factor weighing ³⁾	•	•	•
Check weighing, dynamic weighing, formulation, totaling, back weighing, differential weighing ⁴⁾ , density	•		•
Weight, unit, date/time, sample and task IDs	•	•	•
Target and tolerances	•		•
User name	•		•
Application-specific results and parameter	•		•
Maximum Number of Instruments	Up to 10 balances		
Report Formats	XML, CSV, XLSX, PDF		
Protection	Password protection of settings and results		
EasyDirect Balance Menu Languages ⁵⁾	English, German, French, Spanish, Italian, Russian, Polish, Czech, Hungarian, Dutch, Portuguese, Turkish, Japanese, Chinese, Korean		

1) Full functionality only via Ethernet.

2) RS232 is applicable only for MR.

3) Except for MX/MR balance percent weighing, factor weighing is applicable.

4) Except for MR balance, differential weighing is applicable.

5) Except for MS-TS and ML-T balances with software version 4.20 or higher, when transmitting data, the language of the balance menu must be English.



METTLER TOLEDO Group
Laboratory Weighing
Local contact: www.mt.com/contacts

www.mt.com/EasyDirectBalance

For more information

Subject to technical changes
© 10/2024 METTLER TOLEDO. All rights reserved
30553484C
Group MarCom RITM1258517 MB/JS