

Benchtop Pipettor Service Order Form

800 4 Rainin (800 472 4646) RaininService.com

Step 1 Service Plans and Specifications

Product Name	Standard Calibration	Advanced Calibration	Accredited* ISO 17025 Calibration
As-found Calibration (volumes x weighings)	n/a	2 x 4	3 x 4
As-returned Calibration (volumes x weighings)	2 x 4	2 x 4	3 x 10
Pipettor Type	QTY		
Liquidator 96	17800153 <input type="checkbox"/>	17800154 <input type="checkbox"/>	30415512 <input type="checkbox"/>
BenchSmart 96	30451098 <input type="checkbox"/>	30451099 <input type="checkbox"/>	30451100 <input type="checkbox"/>
MicroPro 300	30848615 <input type="checkbox"/>	30848616 <input type="checkbox"/>	30848617 <input type="checkbox"/>
MicroPro 20 (5 channels only)	30931986 <input type="checkbox"/>	30931987 <input type="checkbox"/>	30932018 <input type="checkbox"/>
MicroPro (Avidien)**	<input type="checkbox"/> 30850283		
Specifications	As-Found	As Returned	
ISO 8655	<input type="checkbox"/>	<input type="checkbox"/>	
Manufacturer	<input type="checkbox"/>	<input type="checkbox"/>	
Custom***	<input type="checkbox"/>	<input type="checkbox"/>	

*Includes Measurement Uncertainty, which is not Included in PASS/FAIL determination.
 **Legacy plan for Avidien customers with prepaid service.
 ***Must be supplied with each order and may not exceed the capability of the pipette.

Step 2 Service Interval

Annual 6 Months 3 Months Other _____

Step 3 Optional Services (Additional fees apply)

96-channel DNA and RNase-off decontamination (30932020)

Step 4 Decontamination Policy

If your organization maintains a current master service agreement with Mettler-Toledo Rainin, LLC and an acceptable decontamination process is specified within such agreement, no signature is required on this form for certification of decontamination.

I certify that my benchtop pipettor is free of biological, chemical and radiological contaminants.

Signed _____ Date _____

Credit card information will be destroyed after processing.

Step 5 Payment Method

Visa Amex MasterCard Purchase Order (please attach)

Credit Card or PO No. _____

Credit Card Expiration	Your Reference Number (optional)
Month: _____ Year: _____	
Rainin Use Only	Promo Code/Voucher No.

Step 6 Billing

Company	_____
Cardholder Name	_____
Address 1	_____
Address 2	_____
Address 3	_____
City, State, ZIP	_____
Quote No.	_____

For current pricing, visit www.mt.com/RaininServicePricing

I need an estimate if additional parts/repairs are needed.

Note: We will gladly design a custom plan to meet your unique needs.
 For complete Terms & Conditions, visit www.mt.com/CaITC

Orders without payment information will not be processed.

Step 7 Shipping

Contact Name	_____
Phone Number	_____
Email	_____

Shipping address same as Billing address

Organization	_____
Address 1	_____
Address 2	_____
Address 3	_____
City, State, ZIP	_____

By submitting this service request I agree to the METTLER TOLEDO standard Terms and Conditions of Sale, which can be accessed of www.mt.com/legal and are incorporated herein by reference. If your organization maintains a current master service agreement with Mettler-Toledo, LLC, the Terms and Conditions of that master agreement will supersede the standard terms and Conditions of Sale as described above and referenced within this pipette calibration order form.

Step 8 Service Centers

7500 Edgewater Drive, Oakland, CA 94621
 150 Wells Avenue, Newton, MA 02459
 5955 Mira Mesa Boulevard, Suite A, San Diego, CA 92121
 200 Rittenhouse Circle East, Unit 2, Bristol, PA 19007

Standard Turnaround
 3-4 days (in house)

Special Instructions:

RaininService

Convenient Mail-in Calibration Service




Rainin makes it easy to get fast reliable service and calibration, guaranteed!

- Quick turnaround, with four ISO 17025 accredited service labs across North America
- We service all brands of pipettes – most with manufacturers’ original parts (Rainin®, Eppendorf®, Gilson®, Thermo Fisher®, Sartorius®)
- We service all pipettes – single channels, multichannels, Liquidator® 96, BenchSmart™ 96 and MicroPro™
- Free shipping and return box for Liquidator and MicroPro – To order a shipping label and return box, contact our Sales team by email at IS@rainin.com or by phone at [800-472-4646](tel:800-472-4646) option 2

Service Plans

All plans include:

- Full preventive maintenance (PM) with seal replacement or restoration
- Liquid leak test on all nozzles
- Adjustment to manufacturer’s tolerances
- Calibration certificate with detailed results
- Calibration label

	Standard Plan	Advanced Plan	Accredited Plan
	Labs with basic service needs 	When prior performance matters 	For complete ISO 17025 documentation 
As Found	N/A	4 Test weighings @ 10, 100 %	4 Test weighings @ 10, 50, 100 %
As Returned	4 Test weighings @ 10, 100 %	4 Test weighings @ 10, 100 %	10 Test weighings @ 10, 50, 100 %
Preventive Maintenance	Yes	Yes	Yes
Calibration Certificate	Non-accredited	Non-accredited	Accredited
Compliance	ISO 8655-7 Annex A	ISO 8655-7 Annex A	ISO 8655-7:2022

Rainin, Pipetting 360+, Liquidator, BenchSmart and MicroPro are trademarks of Mettler-Toledo Rainin, LLC. Eppendorf, Gilson, Thermo Fisher and Sartorius are trademarks of their respective owners.