

Expert Service

For All Pipette Brands...

Gilson, Eppendorf, Thermo Fisher, Sartorius...




Your pipettes need regular check ups. Whatever your service need, you can trust Rainin for exceptional quality and attention to detail.

As one of the largest pipette service providers, Rainin makes it easy to get fast, reliable pipette service and calibration, guaranteed.

- Quick turnaround; four ISO 17025 accredited service centers in the U.S.
- We service all brands of pipettes – most with manufacturers' original parts (Gilson, Eppendorf, Thermo Fisher, Sartorius)
- Free inbound shipping
 - 4 or more pipettes, request our free return kit at mt.com/rainin-kit
 - 3 or less, pack them yourself and print a shipping label at mt.com/RaininServiceLabels



Pipette service is available for these three plans:

Standard Plan	Advanced Plan	Accredited Plan
Labs with basic service needs	When prior performance matters	For complete ISO 17025 documentation
 Perfect for Academic Research R&D	 Pharma Biotech Clinical research labs	 GLP/GMP Regulated Labs

See reverse page for detail information

► www.ShopRainin.com

RAININ

Feel the difference!

Ordering Information

STANDARD PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For labs with basic service needs	✓	n/a	4 test weighings @ 10% and 100%	Standard Certificate

• **Free inbound shipping**

- Full preventive maintenance (PM) with seal replacement and leak test
- Adjustment to manufacturer’s tolerances
- Calibration certificate with detailed results and historical trending

ORDER No.	DESCRIPTION	PRICE
17800031	Single-channel 4 weighings at 2 volume settings, calibration certificate & sticker	\$ 39.00
17800096	Multichannel All – 8 channels	108.40
17800131	Multichannel All – 12 channels	126.00

ADVANCED PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
When prior performance matters	✓	4 test weighings @ 10% and 100%	4 test weighings @ 10% and 100%	Standard Certificate

• **Free inbound shipping**

- Also:
- As found calibration results – including failure confirmation

ORDER No.	DESCRIPTION	PRICE
17800040	Single-channel 4 weighings at 2 volume settings, calibration certificate & sticker	\$ 47.20
17800105	Multichannel All – 8 channels	145.60
17800140	Multichannel All – 12 channels	177.60

ACCREDITED PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For complete ISO 17025 documentation	✓	4 test weighings @ 10%, 50% and 100%	10 test weighings @ 10%, 50% and 100%	ISO 17025 Certificate

• **Free inbound shipping**

- Also:
- As found NVLAP calibration results and failure confirmation
 - NVLAP calibration certificate incl. uncertainty of measurement (ISO 17025 Compliant), historical trending and Calibration label

ORDER No.	DESCRIPTION	PRICE
17800048	Single-channel As Found: 4 weighings at 3 volume settings at 10%, 50% & 100% As Returned: 10 weighings at 3 volume settings	\$ 91.00
17800113	Multichannel All – 8 channels	273.80
17800148	Multichannel All – 12 channels	306.40

For more information on service selection,

please visit

www.RaininService.com

Online orders:

www.ShopRainin.com

Phone orders:

800 4 RAININ (800 472 4646)

Contact your Rainin representative: