

Order status – 800 662 7027 Information – 800 543 4030

To expedite order, please provide the following information: Account No. _____ Contract No. _____

Company/ Institution: _____ **Principal User's Name:** _____

BILLING ADDRESS:

SHIPPING ADDRESS (not PO Box):

MAILING ADDRESS (for reminder notice):

INSTITUTION _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

() _____

PHONE _____

INSTITUTION _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

() _____

FAX _____

INSTITUTION _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

() _____

E-MAIL _____

P.O. Number: _____ **Card Holder Name:** _____

☐ Visa ☐ M/C ☐ AmEx Credit Card Number: _____ Expiration Date: _____

SELECT NEXT DUE DATE: ☐ 3 Months ☐ 6 Months ☐ Annual ☐ Other: _____

ISO 17025 COMPLIANT CERTIFICATES AVAILABLE : ☐ Add \$11 per pipette

CHOOSE YOUR PLAN (all plans include inspection, cleaning, scheduled preventive maintenance, pipette labels, and Calibration Certificates)

PIPETTE QTY	PLAN	CAT. No.	PROGRAM	PIPETTE CHANNELS	AS FOUND DATA			AS RETURNED DATA			LIST PRICE**
					VOLUME SETTINGS	WEIGHINGS	CHANNELS	VOLUME SETTINGS	WEIGHINGS	CHANNELS	
	PLAN A	PA	Single-channel	1	N/A			2	4	1	\$ 37.00
		PA-M	Multichannel ONE SM	8 or 12	N/A			2	4	Center*	59.00
		PA-M8	Multichannel ALL SM	8	N/A			2	4	8	85.50
		PA-M12	Multichannel ALL	12	N/A			2	4	12	102.00
	PLAN B	PB	Single-channel	1	2	4	1	2	4	1	43.50
		PB-M	Multichannel ONE	8 or 12	2	4	Center*	2	4	Center*	65.50
		PB-M8	Multichannel ALL	8	2	4	8	2	4	8	112.00
		PB-M12	Multichannel ALL	12	2	4	12	2	4	12	143.00
	PLAN B w/ EMO	PE	Single-channel	1	3	4	1	3	10	1	55.00
		PE-M	Multichannel ONE	8 or 12	3	4	Center*	3	10	Center*	77.50
		PE-M8	Multichannel ALL	8	2	4	8	2	10	8	165.00
		PE-M12	Multichannel ALL	12	2	4	12	2	10	12	195.50
	PLAN B w/ EMO, EP	PE-M8-P	Multichannel ALL	8	3	4	8	3	10	8	226.50
		PE-M12-P	Multichannel ALL	12	3	4	12	3	10	12	257.50

Please ensure complete decontamination of radioactive and biohazards of pipettes.

Ship in a secure package by FedEx, UPS, or other insurable means to your nearest RAININ Service Center:

7500 Edgewater Drive, Oakland, CA 94621 5955 Mira Mesa Blvd, Ste B, San Diego, CA 92121

36 Middlesex Turnpike, Bedford, MA 01730 200 Rittenhouse Circle East, Unit 2, Bristol, PA 19007

- ☐ Customer SOP on file
☐ Provide Estimate (may delay turnaround)
☐ Special Instructions

EMO = Extended Measurement Option **EP** = Plus Plan

*Full gravimetric evaluation is done on (1) center channel only. All other channels receive one weighing at nominal volume.

**Includes PM parts and labor. 16- and 24-channel pipettes: request a quote. Additional charges: If incoming pipettes fail functional inspection for reasons not covered under preventive maintenance, repairs will be made for additional charge of \$20 labor plus parts.

Terms: Net 30 days, FOB shipping point. Pipettes are returned by overnight service.

Turnaround: 48 hours in Service Center (Multichannel All Plan requires 72 hours).