

**English**

**Installation information**

**Weighing kit, 11106707**

**Deutsch**

**Installationsanleitung**

**Wägekit, 11106707**

**Français**

**Informations d'installation**

**Kit de pesage, 11106707**

**Español**

**Informaciones de instalación**

**Kit de pesaje, 11106707**

**Italiano**

**Informazioni di installazione**

**Kit di pesata, 11106707**

**METTLER TOLEDO**




Before installing "ErgoClips", remove the grid weighing pan from the guides and place it next to the balance.

Für die Montage der "ErgoClips" nehmen Sie die Gitterwaagschale aus der Führung und stellen sie neben die Waage.

Pour le montage des "ErgoClips", retirez le plateau à grille de sa fixation et posez-le à côté de la balance.

Para montar los "ErgoClips", saque el plato de pesada reticular de la guía y colóquelo junto a la balanza.

Per effettuare il montaggio degli "ErgoClips", togliere il piatto a griglia dalle guide e posizionarlo vicino alla bilancia.




Attachment for test tubes

Aufsatz für Reagenzgläser

Support pour tubes à essais

Elemento especial para pesar tubos de ensayo

Elemento per flaconi di reagenti



Attachment for small plastic tubes (side opening)


Aufsatz für kleine Plastikröhrchen (seitliche Öffnung)

Support pour petits tubes plastiques (ouverture latérale)


Elemento especial para pesar tubitos de plástico (abertura lateral)

Elemento per provette in plastica (apertura laterale)






Attachment for round-bottom flasks (various sizes)  
 Aufsatz für Rundkolben (diverse Größen)  
 Support pour ballons de chimie (diverses tailles)  
 Elemento especial para pesar matraces de fondo redondo (diferentes tamaños)  
 Elemento per palloni a fondo sferico (diverse grandezze)



Attachment for weighing boats  
 Aufsatz für Einmal-Wägeschalen  
 Support pour coupelles de pesée jetables  
 Elemento especial para pesar cápsulas de un solo uso  
 Elemento per piatti di pesata monouso



Weighing boat, 500 boats	Part no. 11106712
Einmal-Wägeschalen, 500 Stück	Bestell-Nr. 11106712
Coupelles de pesée jetables, 500 pièces	Référence 11106712
Cápsulas para pesar de un solo uso, 500 unidades	N.º de pedido 11106712
Piatti di pesata monouso, 500 pz.	Numero d'ordine 11106712



Single-use aluminum pan, 10 pans	Part no. 11106711
Fold the side flaps under the bars of the grid weighing pan.	
Einweg Aluminium Wägeschalen, 10 Stück	Bestell-Nr. 11106711
Falzen Sie die seitlichen Lappen unter die Stäbe der Gitterwaagschale.	
Plateaux aluminium jetables, 10 pièces	Référence 11106711
Repliez les languettes latérales sous les barreaux du plateau à grille.	
Platos de pesada de aluminio desechables, 10 unidades	N.º de pedido 11106711
Plique las solapas debajo de las varillas del plato de pesada reticular.	
Piatti di pesata in alluminio, monouso, 10 pz.	Numero d'ordine 11106711
Ripiegare le linguette laterali sotto la griglia del piatto.	