Accessories for Pharmacy Line



LC-P45 Printer

Order No. 00229119 Applications:

- Statistical functions
- Totaling results
- Multiplication and division




Density kit Order No. 00033360 For determining the density of solids.



Anti-theft device

Order No. 00590101 Steel cord with lock secures the balance.





RS/LC-BLD aux. display

Order No. 00224200 On bench stand with a backlit LC display.



RS232 aux. display

Order No. 12120057 The backlit auxiliary display can be placed elsewhere.



B-S AccuPac

Order No. 21254691 For portable weighing. Rechargeable external power supply up to 18 hours of operation.



RS-P26 Printer

Order No. 12120788

Applications:

• Ergonomically designed operation keys like Zero, Date/Time, Feed, Print



USB-232 Cable

Order No. 11103691

The USB-232 converter is an intelligent expansion module that connects to a PC USB port. Supporting baud rates from



LabX direct balance

Order No. 11120340 Simple data input by connecting the balance to a PC.



Weight for tests and adjustments. Guaranteed traceability. For more Information: www.mt.com/weights

Worldwide Service

Swiss Quality



Global Service

Our global service network ensures maximum equipment up-time and accurate results at all times.



Mettler-Toledo AG

Im Langacher

Switzerland

PO Box LabTec

8606 Greifensee

Subject to technical changes

© 03/2007 Mettler-Toledo AG

Printed in Switzerland 11795718

Swiss Quality

Rely on weighing technology from Switzerland.



Pharmacy Line

Accurate

Selection Guide

Easy to use

Dedicated

www.mt.com/pharmacy

For more information



Quality Certificate.

Development, production and tests according to ISO 9001. Environmental management system according to ISO 14001.



"Communauté Européenne".

This mark assures you that our products correspond to the latest guidelines.



For the Right Pharmacy Balance



