# Pharmacy Line

# Balances for pharmacies



# Dedicated

The METTLER TOLEDO Pharmacy Line has been specifically tailored to the needs of pharmacies – to make your everyday weighing as easy as possible.



### Accurate and reliable METTLER TOLEDO is a byword for

absolute accuracy for every weighing job. All our balances are robust and designed for a long service life.

Simple formulation All the balances in the Pharmacy

at the press of a button.


Line have the "Formula" function.

Formulate recipes quickly and easily



# Pharmacy – The new range of balances

Applications in pharmacies are manifold. Handling high-quality active ingredients, for example eye drops, requires precision work and reliable instruments. Dispensing of a wide variety of ointments, creams and medications must be simple and accurate, infusions are often mixed according to the customer's requirements. The proven quality of METTLER TOLEDO meets all the needs of your pharmacy – to your satisfaction and that of your customers. www.mt.com/pharmacy



### Automatic adjustment With FACT fully automatic adjustment,

you can be sure that your balances are always accurate – under all conditions.




## Always well advised


Even good balances must be checked regularly. Take advantage of the expertise of our service engineers. A regular service will extend the service life of your balance. www.mt.com/service




# Technical data – Pharmacy Line

1000 1000	Readability		Max. load	Adjustment		,	Balance dimensions (W/D/H)	Platform dimension mm
	0.0001g	AB104-S/PH	110 g	FACT	0.1 mg	0.2 mg	245/321/344 mm	≥ 80
		AB204-S/PH	210 g	FACT	0.1 mg	0.2 mg	245/321/344 mm	
1	0.001g	PB303-S/PH	310 g	FACT	0.001 g	0.002 g	245/321/236 mm	
		PB503-S/PH	510 g	FACT	0.001 g	0.002 g	245/321/236 mm	> 100
in the		PB503-SDR/PH	100 g/510 g	FACT	0.001 g/0.01 g	0.01 g	245/321/236 mm	≥100
	0.01g	PB3002-S/PH	3100 g	FACT	0.01 g	0.02 g	245/321/89 mm	. 100
		PB3002-SDR/PH	600 g/3100 g	FACT	0.01 g/0.8 g	0.01 g	245/321/89 mm	≥180
	0.001g	PG1003-S/PH	1010 g	FACT	0.0015 g	0.003 g	226/360/239 mm	128x128
	0.01g	PG4002-S/PH	4100 g	FACT	0.008 g	0.02 g	226/360/95 mm	
1 1	0.019	PG4002-SDR/PH	800 g/4100 g		0.008 g/0.03 g	5	226/360/95 mm	165x165
		PG6002-S/PH	6100 g	FACT	0.000 g/0.03 g 0.01 g	0.02 g/0.03 g 0.02 g	226/360/95 mm	
		PG6002-SDR/PH	1200 g/6100 g	FACT	0.01 g/0.03 g		226/360/95 mm	
		- 50002-5000111	1200 g, 0100 g	1101	5.01 g/ 0.03 y	5.02 gi 0.05 g	2201000170 11111	
	0.01g	PG6001-S/MPH	6100 g	FACT	0.03 g	0.05 g	245/321/89 mm	180x180







# Technical data - Special products



Model	PL1501-S	PL6001-S	BBA422-35LA	BBK422-35DLA
Readability	0.1 g	0.1 g	5 g/10 g	0.1 g/1 g
Maximum load	1510 g	6100 g	15 kg/35 kg	7 kg/35.1 kg



LabX direct balance Order no. 11120340 Simple data input by connecting the balance to a PC. If GMP and traceability are issues in your pharmacy, you can find more models on www.mt.com



Quality certificate ISO 9001 Environment certificate ISO 14001 Internet: http://www.mt.com Global service network

Subject to technical changes © 08/2005 Mettler-Toledo GmbH Printed in Switzerland 11795584 MCG MarCom Greifensee

More information www.mt.com/pharmacy

Backlit or fluorescent display for easy-to-read results.

LC-P45 printer.

RS232 interface for connection to a PC or

High-quality, chemical-resistant cast aluminium casing protects the weighing cell and electronics.



METAL

RS 232

Protective cover. Covers the balance casing and protects against dirt and scratches.