
8857


Strip Printer



**Powerful, Compact and Reliable
in Weighing Applications**

METTLER TOLEDO

A Powerful Printer in a Compact Package



Assured compatibility with METTLER TOLEDO products.

A Durable and Effective Printer for Many Applications

- **Fast print speed: 3.5 lines per second**
- **Compact: Occupies minimal space**
- **Dependable: Long-life, 9-pin printhead**

The 8857 Strip Printer is compact, fast, and reliable. Impact printing will not fade over time like thermal printing. So whatever the application, important results are safe for later reference. The 8857 printer is quiet and compact for easy integration into any work area. Drop-in auto-paper loading and few moving parts make this a durable and easy-to-use printer.

Standard Features:

Drop-in auto loading of paper

Universal power supply

Paper-out sensor

Ribbon and paper roll included

| 8857 Specifications | |
|------------------------|--|
| Type | Dot matrix impact strip printer |
| Dimensions (h x w x d) | |
| Printer | 5.2 x 6.3 x 10 (in) 133 x 160 x 250 (mm) |
| Power Supply | 1.7 x 3.2 x 6.5 (in) 45 x 78 x 165 (mm) |
| Weight | |
| Printer | 5 lb (2.2 kg) |
| Power Supply | 1.1 lb (0.5 kg) |
| Print Speed | Approximately 3.5 lines per second @ 16 characters per inch |
| Copies | Up to 1 original + 2 copies / Up to total thickness of 0.008 in |
| Printhead | 9-pin |
| Line Length | 2.5 in (63 mm) (35 columns) |
| Line Spacing | 6 lines per inch |
| Paper Width | 3 in (76 mm) |
| Roll Diameter | 3.27 in (83 mm) |
| Buffer | 1 k |
| Fonts | 9 x 9 dot matrix default (expanded fonts through control codes) |
| Character Size (W x H) | 1.6 x 3.1 (mm) |
| Paper Feed | Friction |
| Power | 100 – 240 VAC / 50-60Hz |
| Operating Environment | 41° to 122° F (5° to 50° C) 10% to 90% relative humidity, non-condensing |
| Ribbon | Snap-in endless loop cartridge |
| Ribbon Life | 1.5 million characters (typical) |
| Interface | RS232 serial, DB25 female connector |
| Approvals | UL, CUL, CE, TUV, FCC Class A |



| | TEL | FAX | TEL | FAX | | | |
|-----------------------|---------------------------|--------------------|--------------------|------------------------|---------------------|---------------------|---|
| METTLER TOLEDO | Australia | (61-3) 9644 5700 | (61-3) 9645 3935 | Korea | (82-2) 518-2004 | (82-2) 518-0813 | Internet |
| | Austria | (43-1) 604 1980 | (43-1) 604 2880 | Malaysia | (603) 7041773 | (603) 7031772 | www.mf.com |
| | Belgium | (32-2) 334 0211 | (32-2) 378 1665 | Mexico | (52-5) 547 5700 | (52-5) 541 2228 | Specifications subject to change without notice. |
| | Brazil | (55-11) 421 5737 | (55-11) 7295-1692 | Netherlands | (31-344) 638363 | (31-344) 638390 | © 2000 Mettler-Toledo, Inc. |
| | China (Industrial) | (86-519) 6642040 | (86-519) 6641991 | Norway | (47-22) 30 44 90 | (47-22) 32 70 02 | METTLER TOLEDO® is a registered trademark of Mettler-Toledo, Inc. |
| | China (Systems) | (86-21) 6485-0435 | (86-21) 6485-3351 | Poland | (48-22) 651 9232 | (48-22) 651-7172 | |
| | Croatia | (385-1) 23 36 317 | (385-1) 23 04 147* | Russia | (7-095) 921 9211 | (7-095) 921 6353 | |
| | Czech Republic | (420-2) 254962 | (420-2) 24247583 | Singapore | (65) 8900011 | (65) 8900012 | |
| | Denmark | (45-43) 27 08 00 | (45-43) 27 08 28 | Slovak Republic | (421-7) 43 42 74 96 | (421-7) 43 33 71 90 | |
| | France | (33-1) 30 97 17 17 | (33-1) 30 97 16 16 | Slovenia | (386-61) 162 1801 | (386-61) 162 1789* | |
| | Germany | (49-641) 50 70 | (49-641) 52951 | Spain | (34-93) 223 2222 | (34-93) 223 0271 | |
| | Hong Kong | (852) 27441221 | (852) 27446878 | Sweden | (46-8) 702 50 00 | (46-8) 642 45 62 | |
| | Hungary | (36-1) 257 9889 | (36-1) 256 2175 | Switzerland | (41-1) 944 45 45 | (41-1) 944 45 10 | |
| | India | (91-22) 857-0808 | (91-22) 857-5071 | Taiwan | (886-2) 2579-5955 | (886-2) 2579-5977 | |
| | Italy | (39-02) 333321 | (39-02) 356-2973 | Thailand | (66-2) 719 6480 | (66-2) 719 6479 | |
| | Japan | (81-6) 6949-5901 | (81-6) 6949-5944 | United Kingdom | (44-116) 235 7070 | (44-116) 236 6399 | |
| | Kazakhstan | (7-3272) 50 63 69 | (7-3272) 60 88 35 | | | | 10M1000 Printers PT03356.0E |

*Call for fax.