

# Compatibility List

Compatibility List for Printer		Category	P-2X Line				P-5X			Last Update
		Article number					30237290	30094673	30094674	
Status	Class	Product Model	USB-P25	RS-P25	RS-P26	RS-P28	P52	P56	P58	Date
Active	Jewellery	JE, JL-G, JP, JS		•	•	•		•	•	Mai 2019
Active	Standard Balances	PL-E		•	•	•		•	•	Mai 2019
Active	Standard Balances	ME, TLE		•	•	•	•	•	•	Mai 2019
Active	Standard Balances	ME-T	•	•	•	•	•	•	•	Mai 2019
Active	Advanced Balances	ML, MS, MS-S, MS-L		•	•	•	•	•	•	Mai 2019
Active	Advanced Balances	ML-T, MS-TS	•	•	•	•	•	•	•	Mai 2019
Active	Excellence Balances	XPR					•	•	•	Mai 2019
Active	Excellence Balances	XSR					•	•	•	Mai 2019
Active	Excellence Balances	XPE		•	•	•	•	•	•	Mai 2019
Active	Comparator	AX						•	•	Mai 2019
Active	Moisture Analyzer	HE53, HE73, HE83		•	•	•	•	•	•	Mai 2019
Active	Moisture Analyzer	HG63, HR83		•	•	•	•	•	•	Mai 2019
Active	Moisture Analyzer	MJ33, HB43-S		•	•	•	•	•	•	Mai 2019
Active	Moisture Analyzer	HS153, HX204		•	•	•	•	•	•	Mai 2019
Active	Moisture Analyzer	HC103		•	•	•	•	•	•	Mai 2019
Active	Density Meter / Refractometer	Excellence D4, D5, D6, R4, R5,	•					•	•	Mai 2019
Active	Density Meter / Refractometer	Easy D40, Easy D30, Easy R40, Easy Brix, Easy Bev	•							Mai 2019
Active	Density Meter / Refractometer	Densito, DensitoPro	•							Mai 2019
Active	Density Meter / Refractometer	Densito 30PX, Refracto 30PX, Refracto 30GS			•					Mai 2019
Active	Titrator	DL15, DL22		•	•	•				Mai 2019
Active	Titrator	DL31, DL38, DL32, DL39		•	•	•				Mai 2019
Active	Titrator	T5, T7, T9	•				•	•	•	Mai 2019
Active	Titrator	V20S, V30S, C20S, C30S, G30S	•				•	•	•	Mai 2019
Active	Melting Point System	MP50, MP70, MP80	•							Mai 2019
Active	Dropping Point System	DP70, DP90	•							Mai 2019
Active	Spectrometer	UV5, UV5Bio, UV7, UV5 Nano	•					•	•	Mai 2019
Active	Industrial Balances	ICS4_5		•	•	•	•			Mai 2019
Active	Industrial Balances	ICS685		•	•	•	•			Mai 2019
Active	Industrial Balances	ICS4_9		•	•	•	•			Mai 2019
Active	Industrial Balances	ICS689		•	•	•	•			Mai 2019
Active	Industrial Balances	IND570		•	•	•	•			Mai 2019
Active	PH	SevenCompact S210, S220, S230, V1.1 or higher		•	•	•	•	•	•	Mai 2019
Active	PH	SevenExcellence (all models)	•	•	•	•	•	•	•	Mai 2019
Active	PH	SevenGo Pro Duo SG68/78/98		•	•	•	•	•	•	Mai 2019
Active	PH	Five EasyPlus		•	•	•	•	•	•	Mai 2019
Active	PH	EL			•	•		•	•	Mai 2019

Compatibility List for Printer		Category	P-2X Line				P-5X			Last Update
		Article number					30237290	30094673	30094674	
Status	Class	Product Model	USB-P25	RS-P25	RS-P26	RS-P28	P52	P56	P58	Date
Phased-out	Jewellery	JB-G/FACT, JB-C/FACT		•	•	•		•	•	-
Phased-out	Jewellery	JB-L-G, JB-L-C		•	•	•		•	•	-
Phased-out	Pharmacy Balances	PG-S, PG-S/PH		•	•	•		•	•	-
Phased-out	Pharmacy Balances	PHL, PHS		•	•	•		•	•	-
Phased-out	Light Balances	PL-S, PL-L		•	•	•	•	•	•	-
Phased-out	Analytical Balances	AB-S/FACT, AB-S, AB-L		•	•	•		•	•	-
Phased-out	Analytical Balances	AB-S/PH, PB-S/PH		•	•	•		•	•	-
Phased-out	Analytical Balances	AB-S/PH, PB-S/PH		•	•	•		•	•	-
Phased-out	Advanced Balances	NewClassic MS, NewClassic ML, NewClassic ME		•	•	•	•	•	•	-
Phased-out	Advanced Balances	PB-S/FACT, PB-S, PB-L		•	•	•		•	•	-
Phased-out	Excellence Balances	XS		•	•	•	•	•	•	-
Phased-out	Excellence Balances	XSE		•	•	•	•	•	•	-
Phased-out	Excellence Balances	XA		•	•	•	•	•	•	-
Phased-out	Excellence Balances	XP		•	•	•	•	•	•	-
Phased-out	Ultra-microbalance	UMT, UMX		•	•	•	•	•	•	-
Phased-out	Comparator	XP Comparator					•	•	•	-
Phased-out	Weighing Platform	QA1, QA3		•						-
Phased-out	Weighing Platform	QB1, QB3, QB5, QB7		•						-
Phased-out	Weighing Platform	QD		•						-
Phased-out	Titrator	DL28		•	•	•				-
Prod-out	Titrator	V20, V30, C20, C30, G20	•				•	•	•	-
Prod-out	Titrator	T50, T70, T90	•					•	•	-
Phased-out	Density Meter	DE40, DE45, DE51, RE40D, RE50, DR40, DR45		•	•	•				-
Phased-out	Density Meter / Refractometer	Densito 30PX								-
Phased-out	Density Meter / Refractometer	LiquiPhysicsDM40, DM45, DM50, RM40, RM50	•				•		•	-
Phased-out	PH	SevenEasy S20, S30, V 1.30 or higher		•	•	•	•	•	•	-
Phased-out	PH	SevenMulti S40, S47, S50, S70, S80, V 1.1 or higher		•	•	•	•	•	•	-
Phased-out	PH	SevenGo pro SG6, SG7, SG8		•	•	•	•	•	•	-
Phased-out	Industrial Balances	SQC IND890						•	•	-
Phased-out	Industrial Balances	SQC XP		•	•	•	•	•	•	-
Phased-out	Industrial Balances	SQC16		•	•	•				-

**METTLER TOLEDO Group**

Laboratory Weighing  
Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

[www.mt.com](http://www.mt.com)

For more information

Subject to technical changes  
© 06/2019 METTLER TOLEDO. All rights reserved

Global MarCom 2761 JC/AG