


LE611

dissolved oxygen



**Order number** 51340401

former order number -

**Specifications**

pH range -

zero point (pH) -

temperature range 0...50 °C

temperature probe -

shaft material PVC

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.1 m

connector BNC

cell constant ( $\text{cm}^{-1}$ ) -

cell type -

indicating metal -

**Storage** 10% NaCl

**Special features**

Measurement range: 0...200 %, 0...20mg/L

Compatible with: FiveGo FG4

Scale 1:1