






NOTES:


1.Expansion+vibration kit may be rotated to any angle.



110kg-1.1t
[250lb-2500lb]



4.4t
[10000lb]



2.2t
[5000lb]

| |
|--|
| <input checked="" type="radio"/> Metrology |
| <input type="radio"/> Safety |
| <input type="radio"/> Ex |

| | | | | | | | |
|----------|--|--|--|-------------------------|--|--|--|
| To Scale | | | | METTLER TOLEDO | | | |
| | | | | MTCN /MTCT / MTCT | | | |
| | | | | EXPANSION/VIBR KIT EVKT | | | |
| | | | | 30136325A | | | |