

## IND360 Precision

To upgrade IND360 Precision firmware to version 2.05.0008, choose one of the two methods below based on current firmware version.

Current firmware version can be found in the web interface under Device --> Info recall --> Main SW.

If the current version is <b>2.02.0000 or higher</b> , use method 1	
	Method 1:  upgrade the firmware using the file IND360_WebUpgrade_Precision_Method1_2.05.0008.mtb
If the current version is <b>lower than 2.02.0000</b> , use method 2 (2 steps)	
	Method 2:  Step 1: Upgrade the firmware using the file IND360_WebUpgrade_Precision_Method2_Step1_2.05.0008.mtb  Step 2: Upgrade the firmware using the file IND360_WebUpgrade_Precision_Method2_Step2_Stack_1.02.0003.mtb